

2.0 INTRODUCTION

The objective of this study was to validate Bayer Method 01056, entitled "Analytical method 01056 for the determination of fenoxaprop-p-ethyl (AE F046360) and its metabolites AE F088406 and AE F054014 in soil using LC/MS/MS" [1]. This method was not successful.

An updated method, entitled "An Analytical Method for the Determination of Fenoxaprop-P-Ethyl, AE F088406 and AE F054014 in Soil Using LC/MS/MS" [2], was generated and successfully validated.

This study was designed to fulfill the requirements of the US EPA Test Guidelines OPPTS 850.7100 [3] and OPPTS 860.1340 [4]. In addition, this study was conducted in compliance with US EPA FIFRA (40 CFR Part 160) GLP standards [5].

3.0 MATERIALS AND METHODS

3.1 Test Substance

Standard name: Fenoxaprop-P-Ethyl (AE F046360)
Bayer Code: K-1822
CAS name: Ethyl (2R)-2-[4-[(6-Chloro-2-benzoxazolyl)oxy]phenoxy]propanoate
CAS number: 71283-80-2
Reference no.: 0218200312
GLP purity: 99.5%
Expiration date: 02 October 2012
Storage conditions: Frozen

Standard name: Fenoxaprop-P (AE F088406)
Bayer Code: K-1821
CAS name: (2R)-2-[4-[(6-Chloro-2-benzoxazolyl)oxy]phenoxy]propanoic acid
CAS number: 113158-40-0
Reference no.: 0626200801
GLP purity: 98.3%
Expiration date: 24 August 2014
Storage conditions: Frozen

Standard name: 6-Chlorobenzoxazolone (AE F054014)
Bayer Code: K-1957
CAS name: 6-Chloro-2(3H)-benzoxazolone
CAS number: 19932-84-4
Reference no.: 1017200312
GLP purity: 98.5%
Expiration date: 09 November 2017
Storage conditions: Frozen

3.2 Test Matrices

The matrix (soil) used for the validation was obtained from Bayer CropScience, 17745 South Metcalf, Stilwell, Kansas, 66085. The matrix was stored at room temperature in a fume hood until needed for analysis.

3.3 Equipment and Reagents

The equipment and reagents used for the method validation were as outlined in Bayer Method 01056 for Trial 1 and 2 (Appendix 2: Apparatus and Reagents Apparatus) and Bayer Method FP-002-S11-02 for Trial 3 (Section 2: APPARATUS and Section 3: REAGENTS). Identical or equivalent apparatus and materials were used.

3.3.1 Equipment and Apparatus

Agilent 1200[®] HPLC System (Agilent Technologies)
Analytical Balance (Mettler Toledo)
API 4000[™] Tandem Mass Spectrometer, MS/MS (Applied Biosystems[™])
EDP3-Plus[®] Electronic Micro Pipettor 200 μ L (Rainin Instrument)
EDP3-Plus[®] Electronic Micro Pipettor 1000 μ L (Rainin Instrument)
Electronic Micro Pipettor 100 μ L (Eppendorf Research pro)
Electronic Micro Pipettor 1000 μ L (Eppendorf Research pro)
Luna C₁₈ (2)-HST HPLC Column (Phenomenex[®])
Manual Micro Pipettor 1000 μ L (VWR International)
Microcentrifuge 16 (Beckman Coulter[®])
Microwave Extractor (CEM-MARSXpress)
N-EVAP (Organomation Associates)
Refrigerator/Freezer (Nor-lake[®] Scientific)
Security Guard Cartridge, C₁₈ (Phenomenex[®])
Top-loading Balance (Mettler Toledo)
Ultrasonic Cleaner 5510 (Branson)
Wrist Action Shaker Model 75 (Burrell)

3.3.2 Reagents

Acetonitrile (EMD)
Acetonitrile (Fisher Scientific)
Ammonium formate, 99% (ACROS)
Formic Acid (EMD[®])
Methanol (EMD)
Methanol (Pharmco-AAPER)
Millipore Water (Millipore[™])
Sodium Chloride (J.T. Baker)
Water (EMD)

3.4 Experimental Design

3.4.1 Establishment of the Method

Prior to performing the ILV, the analyte retention times, instrument detection limits, and linearity of instrument responses to a range of analyte concentrations were determined.

3.4.2 Sample Validation Sets, Fortification, and Extraction Procedure

Sample Validation Sets

Each analytical set consisted of 13 samples: one reagent blank, two untreated controls, five samples fortified with 0.010 ppm of Fenoxaprop-P-Ethyl (AE F046360) and metabolites AE F088406 and AE F054014, and five samples fortified with 0.10 ppm of Fenoxaprop-P-Ethyl (AE F046360) and metabolites AE F088406 and AE F054014.

Data are summarized in Table 1 to Table 3 for Trial 1, Table 4 to Table 6 for Trial 2, and Table 7 to Table 9 for Trial 3. Residue data sheets are included in Appendix 1A. Confirmatory residue data sheets are include in Appendix 1B.

Standard solutions (0.0005 to 0.100 µg/mL) and reagent blanks were also included in each sample set respectively.

Fortification

The control LOQ and 10× LOQ samples were fortified with 0.100 mL and 1.00 mL of the appropriate fortification standard solution respectively.

Extraction and Workup

The following extraction steps were followed for each sample.

1. Weighed 20 g of sample, except for the reagent blank, into a tared 125-mL glass jar and recorded weight.
2. Added the appropriate amount of fortification solution to the sample (0 mL for the blank and the untreated controls, 0.100 mL for the LOQ samples, and 1.00 mL for the 10× LOQ samples).
3. Added 30 mL acetonitrile:water (80:20, v/v) to each sample. The glass jars were capped and shaken.
4. The jars were placed on a mechanical shaker for 20 minutes with the setting on high.
5. The jars were removed and centrifuged at approximately 2000 rpm (not exceeding 2000 rpm) for approximately 3 minutes.
6. After centrifugation, each sample was decanted into glass jars labelled with the sample ID numbers.
7. Added 30 mL acetonitrile:water (80:20, v/v) to each sample. The soil cake at the bottom of each jar was broken up by vigorously shaking.
8. The jars were placed on a mechanical shaker for 20 minutes with the setting on high.

9. The jars were removed and centrifuged at approximately 2000 rpm (not exceeding 2000 rpm) for approximately 3 minutes.
10. The liquid was decanted into the same jars as used in step 6.
11. Added 20 mL of 10-g/L sodium chloride aqueous solution. The soil cake at the bottom of each jar was broken up by vigorously shaking.
12. The jars were placed on a mechanical shaker for 20 minutes with the setting on high.
13. The jars were removed and centrifuged at approximately 2000 rpm (not exceeding 2000 rpm) for approximately 3 minutes.
14. The liquid was decanted into the same jars as used in step 6, and the contents were mixed by swirling.
15. A 0.75-mL aliquot of each sample solution was pipetted into an LC vial containing 0.75 mL methanol:water (70:30, v/v). The vials were capped and vortexed.

3.4.3 Sample Processing and Analysis

The samples were processed and analyzed as described by Bayer Method 01056 [1] for Trial 1 and Trial 2 and by Bayer Method FP-002-S11-02 [2] for Trial 3.

3.4.4 Fortification and Standard Solution Preparation

Trial 1 and Trial 2

Primary stock solutions were prepared as follows. Primary stock solutions for the three standards were prepared by weighing approximately 10.0 mg of each standard into a 100-mL volumetric flask and diluting to volume with acetonitrile. One fortification solution was prepared at a concentration of 20.0 µg/mL by adding an appropriate amount of each primary stock solution to a 10-mL volumetric flask and diluting to volume with acetonitrile. A second fortification solution was prepared at a concentration of 2.00 µg/mL by measuring 1 mL of the first fortification solution into a 10-mL volumetric flask and diluting to volume with acetonitrile.

Secondary standards solutions were prepared at eight concentrations ranging from 0.001 to 0.100 µg/mL by adding an appropriate amount of fortification solution to a 2.0-mL glass autosampler vial and added an appropriate amount of (10:90) methanol/water containing 10 mM ammonium formate and 120 µL/L formic acid.

All solutions were stored in a refrigerator when not in use.

Trial 3

Primary stock solutions were prepared as follows. Primary stock solutions for the three standards were prepared by weighing approximately 5.0 mg of each standard into a 50-mL volumetric flask and diluting to volume with acetonitrile. Fortification solutions were prepared at concentrations of 2.00 µg/mL and 20.0 µg/mL by adding an appropriate amount of primary stock solution to a 10-mL or 100-mL volumetric flask and diluting to volume with acetonitrile.

Secondary standards solutions were prepared at six concentrations ranging from 0.0005 to 0.100 µg/mL by adding an appropriate amount of fortification solution to a 2.0-mL glass autosampler vial and added an appropriate amount of (70:30) methanol/water.

All solutions were stored in a refrigerator when not in use.

3.5 LC-MS/MS Instrumentation

Instrumentation

HPLC System (Agilent 1200[®])

Tandem Mass Spectrometry, MS/MS (Applied Biosystems API 4000[™])

Software: Applied Biosystems, Analyst[®] 1.5

Column: Phenomenex Luna C₁₈ (2)-HST HPLC Column 2.00 × 50 mm, 2.5 µm

3.6 Data Acquisition and Reporting

Peak integration was performed by Analyst[®] software version 1.5. The MS detector responses (peak area) for various injected standard concentrations were used to generate an external calibration curve for the analytes of interest. The overall purpose for the external calibration curve was to display acceptable linearity ($r \geq 0.99$) of the assigned calibration range. The recoveries of the analyte from the fortified samples were calculated by single point calibration.

Recovery results were computed for each sample. The equation used for quantification is presented in Appendix 2. A statistical treatment of the data includes the calculation of averages, standard deviations, RSDs, and the 95% confidence intervals. Mean percent recoveries, standard deviations, RSDs, and 95% confidence intervals were calculated using Microsoft[®] Office Excel 2003. Results were rounded off for reporting purposes but not during calculations.

The Fenoxaprop-P-Ethyl transitions from m/z 362.10 to 288.10 and from m/z 362.10 to 77.10 were used to quantitate and confirm the analyte, respectively. Fenoxaprop-P transitions from m/z 332.00 to 259.90 and from 332.00 to 151.90 were used to quantitate and confirm the analyte, respectively. 6-Chlorobenzoxazolone transitions from m/z 167.70 to 131.90 and from 167.70 to 76.00 were used to quantitate and confirm the analyte, respectively.

Prior to performing the ILV, the analyte retention times, instrument detection limits, and linearity of instrument responses to a range of analyte concentrations were determined.

**TABLE 10 HPLC SYSTEM OPERATING PARAMETERS FOR BAYER
METHOD 01056 (TRIAL 1 AND 2)**

HPLC System: Agilent Model 1200[®]
Software: Applied Biosystems, Analyst[®] 1.5
Analytical Column: Phenomenex Luna C₁₈ (2)-HST HPLC Column 2.00 × 50 mm,
2.5 μm
Column Temperature: 60°C
Injection Volume: 10.0 μL
Run Time: 11.0 minutes

Mobile Phase: (A—Aqueous): 10:90 (v/v) methanol:water (containing 10 mM
ammonium formate and 120 μL/L formic acid)
(B—Organic): 90:10 (v/v) methanol:water (containing 10 mM
ammonium formate and 120 μL/L formic acid)

Gradient:

Time (min)	A (%)	B (%)	Flow (μL/min)
0.00	95.0	5.0	500
4.50	5.0	95.0	500
6.00	5.0	95.0	500
6.10	95.0	5.0	500
11.00	95.0	5.0	500

**TABLE 11 HPLC SYSTEM OPERATING PARAMETERS FOR BAYER
METHOD FP-002-S11-02 (TRIAL 3)**

HPLC System: Agilent Model 1200[®]
Software: Applied Biosystems, Analyst[®] 1.5
Analytical Column: Phenomenex Luna C₁₈ (2)-HST HPLC Column 2.00 × 50 mm,
2.5 μm
Column Temperature: 60°C
Injection Volume: 10.0 μL
Run Time: 10.5 minutes

Mobile Phase: (A—Aqueous): 10:90 (v/v) methanol:water (containing 10 mM
ammonium formate and 120 μL/L formic acid)
(B—Organic): 90:10 (v/v) methanol:water (containing 10 mM
ammonium formate and 120 μL/L formic acid)

Gradient:

Time (min)	A (%)	B (%)	Flow (μL/min)
0.00	95.0	5.0	500
4.00	5.0	95.0	500
7.00	5.0	95.0	500
7.10	95.0	5.0	500
10.50	95.0	5.0	500

TABLE 12 MS/MS OPERATING PARAMETERS

Tandem Mass Spectrometry System, Applied Biosystems, API 4000™
 Software: Applied Biosystems, Analyst® 1.5

The following parameters were used for operation of the mass spectrometer:

Parameter	Setting
Ion Source:	TurboSpray
Scan Type:	MRM
Polarity:	Fenoxaprop-P-Ethyl: Positive Fenoxaprop-P: Negative 6-Chlorobenzoxazolone: Negative
Curtain Gas (CUR):	40.00
Temperature (TEM):	500.00
Ion Spray Voltage (IS):	Period 1(start 0 min): -4000.00 Periode 2 (start 6.10 min): 4000
Collision Gas (CAD):	7.00
Ion Source Gas 1 (GS1):	70.00
Ion Source Gas 2 (GS2):	50.00
Interface Heater (ihe):	ON
Entrance Potential (EP):	-10.00
Transitions Monitored:	(Q1) 362.10→(Q3) 288.10 m/z quantitative (Q1) 362.10→(Q3) 77.10 m/z confirmatory (Q1) 332.00→(Q3) 259.90 m/z quantitative (Q1) 332.00→(Q3) 151.90 m/z confirmatory (Q1) 167.70→(Q3) 131.90 m/z quantitative (Q1) 167.70→(Q3) 76.00 m/z confirmatory
(Q1) 362.10→(Q3) 288.10 m/z	
Collision Energy (CE):	38.00
Collision Cell Exit Potential(CXP):	8.00
Declustering Potential (DP):	78.00
(Q1) 362.10→(Q3) 77.10 m/z	
Collision Energy (CE):	67.00
Collision Cell Exit Potential(CXP):	8.00
Declustering Potential (DP):	78.00
(Q1) 332.00→(Q3) 259.90 m/z	
Collision Energy (CE):	-20.00
Collision Cell Exit Potential(CXP):	-4.00
Declustering Potential (DP):	-63.00
(Q1) 332.00→(Q3) 151.90 m/z	
Collision Energy (CE):	-20.00
Collision Cell Exit Potential(CXP):	-4.00
Declustering Potential (DP):	-63.00
(Q1) 167.70→(Q3) 131.90 m/z	
Collision Energy (CE):	-29.00
Collision Cell Exit Potential(CXP):	-6.00
Declustering Potential (DP):	-63.00
(Q1) 167.70→(Q3) 76.00 m/z	
Collision Energy (CE):	-29.00
Collision Cell Exit Potential(CXP):	-6.00
Declustering Potential (DP):	-63.00

APPENDIX 2 CALCULATIONS

For calculation of the concentrations, calibration curves were used. These curves were calculated automatically after each sequence run with the Applied Biosystems, Analyst[®] software version 1.5 using linear regression. Further calculations were performed using the software EXCEL 2003 (Office 2003[®]).

The linear equation is expressed as:

$$y = MX + B$$

where y = Native peak area
 M = Calibration line slope
 X = Concentration of the reference standard in ng/mL
 B = Calibration line intercept

Trial 1 and Trial 2

By means of the linear equation, the content of Fenoxaprop-P-Ethyl and its metabolites, AE F088406 and AE F054014, in dry soil or recoveries can be calculated as follows:

Area Analyte	y	
Standard Concentration	x	(µg/L) / (µg/L)
Sample Weight	G	kg
Fortified Amount	A	mg
Final Volume	VEND	L
Residue in Dry Soil (only used for residue sample)	R	mg/kg
Recovery	Rec	%
Moisture	M	

$$R = \frac{y - \text{intercept}}{\text{slope}} \times \frac{\text{VEND}}{G} \times \frac{1}{1 - M}$$

$$\text{Rec} = \frac{y - \text{intercept}}{\text{slope}} \times \frac{\text{VEND}}{A} \times 100\%$$

Trial 3

With each sample set, analyze an untreated control sample and one or more fortified control samples. Calculate recoveries using the following equation:

$$\text{Recovery (\%)} = \frac{(R - S)}{T} \times 100$$

where R = ppb of target analyte found in fortified sample
S = ppb of target analyte found in control sample, real or apparent
T = theoretical ppb in fortified sample

The dilution factor was calculated using the following equation:

$$\text{Dilution factor} = \frac{V1}{W} \times \frac{V3}{V2}$$

where V1 = Initial volume (mL)
W = Initial sample weight (g)
V3 = Final volume (mL)
V2 = Aliquot (mL)

Bayer CropScience

FP-002-S11-02

1. PRINCIPLE

The residues of Fenoxaprop P-Ethyl and its metabolites AE F088406, and AE F054014 in soil are determined by shaking a 20 gram soil sample twice with acetonitrile/water and once with water containing sodium chloride. An aliquot of the sample is diluted in methanol/water and analyzed by LC/MS/MS. Quantification is based on a comparison of peak areas with those of known standards.

The method limit of quantitation (LOQ) for Fenoxaprop-P-Ethyl and its metabolites AE F088406 and AE F065025 is 10ng/g in soil.

2. APPARATUS

Use as a guide; equivalent apparatus may be substituted.

- VWR Pyrex® Brand volumetric pipets, glass class A (Assorted Volumes)
- Eppendorf Reference Series 2000 pipettes (Cat. No.: 05-402-48 and 05-402-50)
- VWR Pyrex® Brand volumetric flasks, glass class A (Assorted Volumes)
- VWR Pyrex® Brand disposable Pasteur pipets (Cat. No.: 53283-910 & 53283-914)
- National Scientific LC vials, Snap-Its (Cat. No.: C4011-5)
- National Scientific LC vial Snap-It Seals, (Cat. No.: C4011-55)
- Phenomenex® Luna C18(2)-HST 50 x 2.00 mm Column (Part No.: 00B-446-B0)
- ABSciex API 4000 chromatograph/mass spectrometer (LC-MS/MS) equipped with electrospray ionization (ESI) interface, Shimadzu HPLC pumps and a CTC PAL autosampler, and Analyst 1.4.1 data collection software or higher version, or equivalent
- Fisherbrand 125-mL glass jars (Cat. No. 02-911-455)
- Centrifuge
- Mechanical Shaker
- Various general laboratory glassware and utensils

3. REAGENTS

Use as a guide; equivalents or different manufactures (brands) may be substituted.

- Acetonitrile, (HPLC Grade)
- Water (HPLC Grade or Millipore)
- Methanol (HPLC Grade)
- Sodium Chloride crystals (Mallinckrodt)
- 10 grams/liter sodium chloride aqueous solution. Dissolve 10g of sodium chloride in 1L of water (HPLC Grade or Millipore)
- Methanol/water solution (70/30). Combine 700mL of methanol with 300mL of water and mix well
- Acetonitrile/water solution (80/20). Combine 800mL of acetonitrile with 200mL of water and mix well

Bayer CropScience

FP-002-S11-02

- Certified analytical reference standards of Fenoxaprop-P-Ethyl, AE F088406 and AE F054014

4. PREPARATION OF ANALYTICAL STANDARDS

NOTE: The following procedure is an example description of how standard solutions may be prepared. Standards may be prepared as mixed solutions by dilution from individual stock solutions or prepared individually. Alternate or additional standards of appropriate weight and volume may be prepared as needed.

Class "A" volumetric glassware or calibrated pipets should be used in the preparation of all analytical standards. All standard solutions should be stored in a refrigerator in amber glass bottles when not in use. Solutions should be allowed to warm to room temperature prior to use.

4.1 Primary Stock Standard Solutions

Prepare individual 100µg/mL stock solutions of Fenoxaprop-P-Ethyl, AE F088406, and AE F054014 by transferring 0.0100 grams of each analyte in separate 100mL volumetric flasks. Dilute to volume with acetonitrile and mix well. Store at <-10°C when not in use.

Prepare individual 20µg/mL solutions of Fenoxaprop-P-Ethyl, AE F088406, and AE F054014 by taking a 20.0mL aliquot of the 100µg/mL stock solutions and diluting to 10mL with acetonitrile. Store at <-10°C when not in use.

Prepare individual 2.0µg/mL solutions of Fenoxaprop-P-Ethyl, AE F088406, and AE F054014 by taking a 2.0mL aliquot of the 100µg/mL stock solutions and diluting to 100mL with acetonitrile. Store at <-10°C when not in use.

NOTE: Corrections for standard purities should be applied when expressing standard concentrations. AE F088406 and AE F054014 should be expressed as parent equivalent.

4.2 Fortification Standard Solutions

The individual 2.0µg/mL solutions of Fenoxaprop-P-Ethyl, AE F088406 and AE F054014 prepared in section 4.1 above may be used as the fortification solutions.

Further dilutions of this mixed fortification solution may be made as needed.

4.3 Calibration Standard Solutions

Prepare working calibration solutions consisting of 0.5, 1, 5, 10, 50, and 100ng/mL of Fenoxaprop-P-Ethyl, AE F088406, and AE F054014 as described in the following table and diluting all solutions to 100mL with 70:30v/v methanol:water.

Concentration of Standard Solution used for dilution ($\mu\text{g/mL}$)	Aliquot Native mix Taken (mL)	Dilution Volume (mL)	Concentration of Calibration Solution (ng/mL)
2	0.025	100	0.5
2	0.05	100	1
2	0.25	100	5
2	0.5	100	10
20	0.25	100	50
20	0.5	100	100

Additional calibration solutions may be prepared as required. Refrigerate when not in use. All solutions are stable for approximately one week.

4.4 Extraction

NOTE: The analytical targets in this method are subject to rapid degradation. It is important to complete the extraction in a timely manner.

1. Transfer 20 ± 0.2 grams of soil into a 125ml glass jar.
2. Fortify the recovery samples at the desired fortification level using the individual standard solutions prepared in acetonitrile (see Section 4.2 Fortification Stock Solutions).
3. Add 30mL of acetonitrile/water solution (80/20). Cap and shake glass jar.
4. Place jar on a mechanical shaker for 20 minutes with the setting on high.
5. Remove the jar and centrifuge at approximately 2000 rpm for about 3 minutes. Because glass jars are used, the centrifuge speed should not exceed 2000 rpm because at higher rpm the glass jars could break.
6. After centrifugation, decant each sample into a 125ml glass jar labeled with the sample ID number.
7. Add 30mL of acetonitrile/water solution (80/20). Shake vigorously by hand to break up the soil cake at the bottom of the jar before placing on the shaker.
8. Place jar on a mechanical shaker for 20 minutes with the setting on high.
9. Remove the jar and centrifuge at approximately 2000 rpm for about 3 minutes.

Bayer CropScience

FP-002-S11-02

10. Decant the liquid into the same bottle used in step 4.4.6.
11. Add 20mL of the 10g/L sodium chloride aqueous solution. Shake vigorously by hand to break up the soil cake at the bottom of the jar before placing on the shaker.
12. Place jar on a mechanical shaker for 20 minutes with the setting on high.
13. Remove the jar and centrifuge at approximately 2000 rpm for about 3 minutes.
14. Decant the liquid into the same bottle used in step 4.4.6 and swirl to mix the contents.
15. Pipet a 0.75mL aliquot of the sample solution into a LC vial containing 0.75mL of 70:30v/v methanol:deionized water. Cap the vial and mix to await analysis by LC/MS/MS.

5. ANALYSIS

5.1 Sample Analysis

- Step 1. Using the recommended procedures listed below; analyze an aliquot of the 0.5, 1, 5, 10, 50, and 100ng/mL standard solutions (these are calibration solution analyses).
- Step 2. Analyze an aliquot of the analytical samples. Note: Up to 25 sample analyses can be made after the analysis of the standard solutions.
- Step 3. Again analyze an aliquot of the 0.5, 1, 5, 10, 50, and 100ng/mL calibration standard solutions.
- Step 4. When necessary, analyze additional samples and standard solutions. Always finish the procedure with the analysis of a set of standard solutions.

5.2 LC/MS/MS Standard Calibration and Residue Calculations

Standardize the LC/MS/MS response under the conditions outlined in Appendix 1 by injecting an aliquot of each LC/MS/MS calibration solution interspersed with samples.

The residues of Fenoxaprop-P-Ethyl and its metabolites AE F088406 and AE F054014 are quantified using external standard linear regression analysis. A separate calibration curve was produced for each set of samples analyzed on the LC/MS/MS. A calibration curve was generated by linear regression of the standard peak area versus the standard concentrations in ng/mL using Applied Biosystems Analyst Software (Version 1.4.1), a computer-programmed data capturing system. The Analyst Software uses the MS/MS standard responses to calculate the regression coefficients M and B, respectively called slope and intercept, for each analytical set.

The standards were fit to the linear equation: $Y = MX + B$

Bayer CropScience

FP-002-S11-02

where: X is the concentration of the reference standard in ng/mL
M is the calibration line slope
B is the calibration line intercept
Y is the native peak area

The calibration points are weighted 1/x to provide better fit near the limit of detection.

After regression coefficients were calculated, the residue in parts per billion was determined. The parts per billion (ppb) of fenoxaprop-p-ethyl in the soil was calculated using the following equation,

$$\text{Fenoxaprop-P-Ethyl found (ppb)} = \frac{(Y-B) \times D}{M}$$

$$\text{Where Dilution Factor (D)} = \frac{\text{Initial volume (V1)}}{\text{Initial sample wt. (W)}} \times \frac{\text{Final volume (V3)}}{\text{Aliquot (V2)}} = 8$$

Where: W = 20g
V1 = 80mL
V2 = 0.75mL
V3 = 1.5mL

Analyst software was used to calculate the residues of Fenoxaprop-P-Ethyl in ppb for each sample and the percent recovery for the fortified samples. AE F088406 and AE F054014 residue concentrations are determined in a similar fashion.

5.3 Fortification Experiments

Note: Fortification experiments may be performed as needed to monitor method efficiency and reproducibility, but are not required when analysis of samples is performed for tolerance enforcement. Fortification experiments are intended to be used for data collection methods or establishing & validating method efficiency.

With each sample set, analyze an untreated control sample and one or more fortified control samples. Calculate recoveries using the following equation:

$$\text{Recovery (\%)} = \frac{(R-S)}{T} \times 100$$

Where: R = ppb of target analyte found in fortified sample
S = ppb of target analyte found in control sample, real or apparent
T = theoretical ppb in fortified sample

Recoveries are determined by analyzing fortified control samples alone or in conjunction with a sample set. Samples may be fortified prior to extraction at the LOQ of 0.01ug/g or

Bayer CropScience

FP-002-S11-02

other appropriate level with fortification solutions. Calculate the final residue R for the control (S) and fortified control (R) samples.

Bayer CropScience

FP-002-S11-02

Appendix 1 Instrument Conditions For Fenoxaprop-P-Ethyl and its metabolites

Equipment with equivalent or better sensitivity and performance may be substituted.

LC/MS/MS Parameters

NOTE: Variations in equipment or sample characteristics may require slight modifications in the chromatographic or detector conditions listed in order to obtain adequate chromatographic peak shapes or sensitivity. Therefore, the given LC/MS/MS parameters listed below are guidelines and may be modified. These parameters should be optimized for the instrument and column actually used. Also, instrument parameters and mobile phase may be adjusted to improve separation from any observed interfering peaks.

The following conditions were used on an ABSciex API 4000 LC/MS/MS system.

HPLC Parameters

Pumps Used: Two Shimadzu LC-10ADVP (High Pressure Mixer) pumps with a Shimadzu SCL-10 controller
Autosampler CTC PAL
Column Temperature: 50°
Injection Volume: 10µL
Column: Manufacturer: Phenomenex®
Type: Luna C18 (2)-HST
Particle Size: 2.5 µm
Diameter: 2.0 mm
Length: 50 mm

Mobile Phase A: Water/Methanol 90:10 (v/v) with 10mmol/L ammonium formate and 120µL formic acid/L
Mobile Phase B: Methanol/Water 90:10 (v/v) with 10mmol/L ammonium formate and 120µL formic acid/L

HPLC gradient program:

Time (min.)	Module	Flow Rate (mL/min)	A(%)	B(%)
0.0	Pumps	0.50	95	5
0.5	Pumps	0.50	95	5
2.5	Pumps	0.50	5	95
6.0	Pumps	0.50	5	95
7.0	Pumps	0.50	5	95
7.1	Pumps	0.50	95	5
8.0	System Controller	Stop		

Bayer CropScience

FP-002-S11-02

Mass Spectrometer Instrument Conditions


Component:	Fenoxaprop-P-ethyl	Fenoxaprop-P-ethyl	AE F088406	AE F088406	AE F054014	AE F054014
Retention Time	5.2 minutes	5.2 minutes	5.1 minutes	5.1 minutes	3.4 minutes	3.4 minutes
Transition	MRM1	MRM2	MRM1	MRM2	MRM1	MRM2
Parent Ion	363	363	332	332	168	168
Product Ion	288	77	260	152	132	76
Ionization Mode	ESI	ESI	ESI	ESI	ESI	ESI
Polarity	+	+	-	-	-	-
Dwell Time (ms)	600	600	100	100	100	100
Declustering Potential (DP)	76	76	-70	-65	-65	-70
Entrance Potential (EP)	10	10	-10	-10	-10	-10
Collision Energy (CE)	27	83	-18	-30	-26	-34
Collision Cell Exit Potential (CXP)	10	4	-7	-11	-7	-5
Curtain Gas (CUR)	30	30	30	30	30	30
Collision Gas (CAD)	6	6	6	6	6	6
Ion Source Gas 1 (GS1)	30	30	30	30	30	30
Ion Source Gas 2 (GS2)	50	50	50	50	50	50
Source Temp (TEM)	500	500	500	500	500	500
Interface Heater (IHE)	On	On	On	On	On	On
Ion Transfer Voltage (IS)	4500	4500	-4500	-4500	-4500	-4500

Bayer CropScience

FP-002-S11-02


Appendix 2 Structures

Chemical Name: Fenoxaprop-P-Ethyl
(Parent Molecule)



CAS Name: Ethyl (2R)-2-[4-[(6-Chloro-2-benzoxazolyl)oxy]phenoxy]propanoate
CAS Number: 71283-80-2
Molecular Formula: C₁₈H₁₈ClNO₅
Molecular Weight: 361.7763

Chemical Name: AE F054014
(Metabolite)




CAS Name: 6-Chloro-2(3H)-benzoxazolone
CAS Number: 19932-84-4
Molecular Formula: C₇H₄ClNO₂
Molecular Weight: 169.5652

Bayer CropScience

FP-002-S11-02

Appendix 2 (continued)

Chemical Name: AE F088406
(Metabolite)



CAS Name: (2R)-2-[4-[(6-Chloro-2-benzoxazolyl)oxy]phenoxy]propanoic acid
CAS Number: 113158-40-0
Molecular Formula: C₁₆H₁₂ClNO₅
Molecular Weight: 333.7232