


PAMS Continuous VOC Monitoring

AUTOGE

Overview of Chromatographic Applications

Carol J. Meyer
carol@orsat.com



Orsat, LLC
www.orsat.com

Not-so-continuous VOC Monitoring: AutoGC ~1 Hour Composite Sample

- ▶ Sample Collection on Sorbent Trap (composite sample to correlate with hourly criteria measurements)
- ▶ Chromatographic separation of target compounds (generally separation occurs during collection of next sample)
- ▶ Detection of separated target compounds by one or more detection methods (FID, PID, MS)
- ▶ Data analysis to identify and quantitate separated targets compounds (Chromatographic data system supplied with each instrument)

Sample Collection Methods

- ▶ Collection on a single trap
 - ▶ Cryogenically cooled trap to capture C2-C12 VOCs
 - ▶ Dual adsorbent trap for full range collection
 - ▶ Single sample injected to multidimensional gas chromatograph
- ▶ Separate C2-C6 and C6+ traps
 - ▶ C2-C6 single adsorbent trap cryogenically cooled to trap C2 VOCs
 - ▶ C6+ single adsorbent trap for C6+ VOCs
 - ▶ Each trap injected into separate gas chromatograph

NMHC Monitoring Strategies: Analytical Choices

GC-FID

- ▶ Low cost
- ▶ Stable < 2% drift over 1 month
- ▶ Response relative to carbon content
- ▶ Linear detector response
- ▶ Not sensitive to O₂, N₂ or H₂O
- ▶ Possible interferences

GC-PID


- ▶ Low cost
- ▶ Uses lamp which can age and require cleaning
- ▶ Non-uniform detector response
- ▶ Reduced Linear range
- ▶ Sensitive to H₂O interference
- ▶ Not sensitive to interferences

GC-MS

- ▶ More complex data
- ▶ More complex operation
- ▶ Requires more frequent calibration
- ▶ Non-uniform detector response
- ▶ Sensitive to O₂, N₂ or H₂O interference
- ▶ Not sensitive to interferences


AutoGC Systems: Single Trap Systems

- ▶ Dual adsorbent cryogenically cooled trap
- ▶ Detection by Dual FID, MS or FID-MS
- ▶ Single Data System



AutoGC Systems: Dual Trapping Systems

- ▶ Dual Chromatographic Systems - separate samples
- ▶ Detection by Dual FID, PID
- ▶ Single Data System



Ancillary Equipment Requirements

- ▶ Carrier gas - Helium ~10-15 mL/min
- ▶ Air supply - Compressor
 - ▶ Purification Systems to supply
 - ▶ Dry air or N₂ < 1ppm moisture for Peltier cooler and nafion drier ~ 350 mL/min
 - ▶ Hydrocarbon-free air for FID support gas ~ 800 mL/min
 - ▶ Dilution gas for dilution system ~150 - 300 mL/min
- ▶ Hydrogen if FID
 - ▶ Cylinder gas or Hydrogen generator ~80 mL/min
- ▶ Sample manifold and sample pump
- ▶ Canisters and dilution system
 - ▶ Calibration curves
 - ▶ Check standards

Sampling System

- ▶ Standard glass manifold systems
 - ▶ Glass with blower to bring sample from outside
 - ▶ Heated to eliminate condensation
 - ▶ Cleaned regularly with only de-ionized water
- ▶ Sample lines from manifold
 - ▶ Flow rates < 50 sccm should be 1/8" or smaller
 - ▶ Heated stainless steel - chromatographic grade only
 - ▶ Silco treated stainless steel
- ▶ Problems associated with sample lines
 - ▶ Losses of heavy compounds - inadequate humidification
 - ▶ Contaminations - inadequately cleaned tubing
 - ▶ Carry-over of heavy compounds - inadequate humidification or condensation (may require heating)

Chromatography Data Systems

Requirements:

- ▶ Data portability
- ▶ Ability to reconstruct the original processing method from result
- ▶ Use of Retention time references to accommodate diurnal shifts
- ▶ Use of response factors and calibration by reference for unidentified HCs
- ▶ Ability to name files for easy identification of site, date, time, hour and sample type
- ▶ Ability to schedule and control introduction of routine quality control samples
- ▶ Ability to recover from simple power failures and continue hourly sampling
- ▶ Remote operation

Method Development TotalChrom Chromatography Data System

The screenshot displays the TotalChrom Method Editor interface, showing two windows. The top window is titled 'Method Editor - C:\Users\manager\Data Validation\TD300 Office Data\loc-meth\XX-PLOT.mth'. The bottom window is titled 'Method Editor - F:\Orsat\PE TD300 System\OZONE_A_031015.mth - [Component List]'. Both windows show a 'Component List' table, a 'Method Summary' panel, and a 'Data Processing and Reporting' section.

Component List (Top Window):

Number	Time	Component Name	Component Type	METHOD
1	6.752	METHANE	Single Peak Component	
2	7.437	ETHANE		
3	8.397	ETHYLENE		
4	10.860	PROPANE		
5	18.659	PROPYLENE		
6	21.812	ISOBUTANE		
7	23.085	N-BUTANE		
8	24.835	ACETYLENE		
9	28.574	TRANS-2-BUTENE		
10	29.138	1-BUTENE		
11	30.568	CIS-2-BUTENE		
12	31.793	CYCLOPENTANE		
13	32.193	ISOPENTANE		
14	33.018	N-PENTANE		
15	34.266	1,3-BUTADIENE		
16	35.918	TRANS-2-PENTENE		
17	36.529	2-METHYL-2-BUTENE		
18	36.904	1-PENTENE		
19	37.456	CIS-2-PENTENE		
20	38.969	2,2-DIMETHYLBUTANE		
21	39.711	2-METHYLPENTANE		
22	40.692	ISOPRENE		

Method Summary (Top Window):

Data Acquisition and Instrument Control

Instrument Name : Site-XX Injection : MANUAL
 Experiment Time : 20.00 min Injection Volume :
 Delay Time : 0.00 min Sampling Rate : 3.12500 pts/s
 Run Time : 20.00 min Channel : DUAL

Oven Temperature Program :
 Initial Temperature : 50 deg for 20.00 min

Data Processing and Reporting

Replot Pages : 3 BF : 3 User Prog

Component List (Bottom Window):

Number	Time	Component Name	Component Type	METHOD
1	11.593	N-Hexane 110-54-3	Single Peak Component	
2	12.298	Methylcyclopentane 96-37-7		
3	12.671	2,4-Dimethylpentane 108-08-7		
4	14.252	Carbon Tetrachloride 56-23-5		
5	15.186	Benzene 71-43-2		
6	15.551	Cyclohexane 110-82-7		
7	18.573	2-Methylhexane 591-76-4		
8	18.953	2,3-Dimethylpentane 565-59-3		
9	19.387	3-Methylhexane 599-34-4		
10	19.976	2,2,4-Trimethylpentane 540-84-1		
11	20.545	N-Heptane 142-82-5		
12	21.395	Methylcyclohexane 108-87-2		
13	23.747	2,3,4-Trimethylpentane 565-75-3		
14	25.447	Toluene 108-88-3		
15	26.296	2-Methylheptane 592-27-8		
16	26.854	3-Methylheptane 599-81-1		
17	27.966	Tetrachloroethylene 127-18-4		
18	28.556	N-Octane 111-65-9		
19	33.643	Ethylbenzene 100-41-4		
20	34.143	M/P-Xylene 108-38-3 & 106-42-3		
21	34.264	O-Xylene 95-47-6		
22	35.514	Styrene 100-42-5		
23	35.587	N-Nonane 111-84-2		
24	36.837	Isopropylbenzene 98-82-8		
25	38.328	N-Propylbenzene 103-65-1		
26	38.621	M-Ethyltoluene 620-14-4		
27	38.763	P-Ethyltoluene 622-96-8		
28	39.976	1,3,5-Trimethylbenzene 108-67-8		
29	39.020	O-Ethyltoluene 611-14-3		
30	39.634	1,2,4-Trimethylbenzene 95-63-6		
31	40.265	N-Decane 124-18-5		
32	41.541	1,2,3-Trimethylbenzene 526-73-8		
33	42.103	M-Diethylbenzene 141-93-5		
34	42.427	P-Diethylbenzene 105-05-5		
35	42.635	N-Undecane 1120-21-4		

Method Summary (Bottom Window):

Component Type : Benzene 71-43-2
 Retention : 15.186min
 Window : 2.00s 2.00%

Reference :
 ISTD component :
 Response :
 Curve : 1st Order - Include Origin
 Purity % : 100.0000
 Scaling : 1.0
 Weighting : 1.0

Data Processing and Reporting (Bottom Window):

Level	Amount	Area/Height
11	3.0000	2.720764e+05
12	6.0000	6.234720e+05
13	9.0000	1.127497e+06
14	12.0000	1.632762e+06

Calibration Curve Graph (Bottom Window):

R-squared = 0.999116
 $Y = (8828.853446)X + (134731.811035)X$

The graph plots Response (Y-axis, ranging from 311042 to 1811042) against Vol Adj Amt (X-axis, ranging from 2.5 to 10.6). Four data points are shown as blue squares, and a red linear regression line is fitted through them.

Calibration

Carbon Response Factor

- ▶ Single Response factor for each column
- ▶ Response factor based on response (peak area) per ppbC
- ▶ Requires certified standard for only Propane and Benzene
- ▶ Used only on FID systems
- ▶ Easy for validators to check for errors

Target Specific Regression

- ▶ Linear regression for each target component
- ▶ Response based on ppbv
- ▶ Requires standard containing all targets
- ▶ Required on PID and MS systems which do not have uniform response across targets
- ▶ Difficult to achieve and validate

Calibration


Mass Spectrometry or Photoionization Detector:

- ▶ May require calibration across larger dynamic range due to non-linear response (6 points)
- ▶ Generally requires target specific calibration
- ▶ Will require standards for each target to be calibrated


FID Detector

- ▶ Has mostly uniform response based on Carbon-Hydrogen bonds
- ▶ Extremely linear and may be calibrated with fewer points (1 - 3)
- ▶ Only requires standard with propane and benzene for calibration


C2-C6 Alumina PLOT Separation




C6+ Dimethylsiloxane (Boiling Point) Separation




Analytical Blank



Ambient Air - PLOT Column




Ambient Air - Boiling Point Column




Recoveries in Daily Check Standard at 0.5 ppbv

Check Standard % Recovery - C3/C6 Average Response Factor



Ambient Data - Carbon Response Factor Calibration



Sources of Data Losses and Poor QC Results

Data Losses

- ▶ Power failures
- ▶ Sample collection time errors
- ▶ Shelter AC failures
- ▶ Compressor failures
- ▶ “Wet” air

Poor QC Results

- ▶ Canisters
- ▶ Sample lines
- ▶ Sample pumps/flow
 - ▶ Trap icing
- ▶ Losses of light gases
 - ▶ Failure of trap cooling
 - ▶ Sample or dry purge flow too high
- ▶ Contaminants in Blank

Sources of Carry-over and Interferences


Carry-over

- ▶ Extremely high ambient
- ▶ Nafion Driers
- ▶ Sample lines
 - ▶ Heated SS
 - ▶ Treated Steel
- ▶ Trap failure resulting in trap material in system
- ▶ Ferrules


Interferences

- ▶ Canisters
- ▶ Failing Air system
- ▶ Site specific interferences
- ▶ Leak testing solutions
- ▶ Baseline issues
 - ▶ Detector ferrules
 - ▶ Air supply
 - ▶ Electrical
 - ▶ Vibration


Sources of Carry-over and Interferences




Sources of Carry-over and Interferences



Sources of Carry-over and Interferences



Sources of Carry-over and Interferences



Data Quality Objectives

- ▶ Data completeness
- ▶ Data representative of ambient concentration
- ▶ Minimize contributions from the system

Operations

- ▶ Well defined Operating Procedures
- ▶ Well documented instrumental parameters
- ▶ Fully automated system to reduce errors in operator activities
- ▶ Easily identifiable and transportable data files
- ▶ Fully automated Quality Checks

Data Validation

- ▶ Well defined Data Quality Objectives
- ▶ Real-time data transfer and review
- ▶ Well defined validation operating procedures
- ▶ Good annual audits to review instrument performance across network.

Network Quality Control

Quality Control Check	Composition	Purpose	Frequency	Acceptance Criteria
Retention Time Standard (RTS)	Mixture containing all target compounds ideally between 1-5 ppbC	To help assess retention time shifts and optimize processing methods	Twice a month or weekly	100% of the compounds are identified correctly in the multicomponent RTS
Calibration Verification Standard (CVS)	Mixture of 15 reference compounds including Propane and Benzene used for calibration	To assess the instrument drift and ensure continued instrument calibration	Daily	1) Propane and Benzene % recoveries within 75% - 125% and all other calibrants within 55 - 145% 2) Data must be bracketed by valid CVS
Method (Analytical) Blank	Humidified, clean air	To assess system contribution to the measurement	Daily	1) All target compounds < 2.0 ppbC 2) TNMHC < 20 ppbC 3) Data must be bracketed by valid blanks
Precision Check	Mixture used for CVS	To assess analytical precision	Weekly	Propane and Benzene %RPD < 20% in two consecutive CVS runs
Laboratory Calibration Standard (LCS)	Mixture of 15 reference compounds including Propane and Benzene used for calibration	Second source standard, statically blended 5 ppbv	Twice a month or weekly	Propane and Benzene % recoveries within 70-130%

Data Review - Data Validation

Data Review - Daily

- ▶ Site Operations
 - ▶ On-time collection
 - ▶ Correct identification
 - ▶ Equipment parameters
- ▶ Quality Controls
 - ▶ Passing Blanks
 - ▶ Passing check standard recovery

Data Validation - Monthly

- ▶ Review of Quality Controls
 - ▶ Passing check standard recoveries - flagging
 - ▶ Passing blanks - flagging failed targets
 - ▶ Retention time checks
- ▶ Chromatography review
 - ▶ Review of high hours or other issues

Electronic Logbooks

AutoGC Operator Log Entry

Site:

Operator:
Use Groupwise, Rhone, or AutoGC Validation ID, usually formatted as first initial followed by the first seven letters of last name.
Ex.) Melanie Hotchkiss = mhotchki

Pressures, Flows, and other Checks
Update the information as needed

Helium:	psig	Comment:	Nafion Flow:	slp	C
Hydrogen 1:	psig	Comment:	TOC Flowmeter:	slp	C
Hydrogen 2:	psig	Comment:	TD Sample Pump:	ml/min	C
Zero Air:	psig	Comment:	TD Outlet Split Flow:	ml/min	C
CVS:	psig	Comment:	TD Desorb Vent Flow:	ml/min	C
LCS:	psig	Comment:	TD Column Pressure:	psig	C
RTS:	psig	Comment:	GC Midpoint Pressure:	psig	C
Audit:	psig	Comment:	Other:		C
MDL:	psig	Comment:			
GC Cut Time:	min	Comment:	<input type="checkbox"/> Power Failure		C
Station Temp:	°F	Comment:	<input type="checkbox"/> Verified sample collection within the hour (sample timing)		C

Canister Installation
Update information below when canisters are installed / replaced; please include "CC", "ALM", etc. on cylinder information if applicable

<input type="checkbox"/> CVS	Can #:	Cyl #:	Tst #:	Comment:
<input type="checkbox"/> LCS	Can #:	Cyl #:	Tst #:	Comment:
<input type="checkbox"/> RTS	Can #:	Cyl #:	Tst #:	Comment:
<input type="checkbox"/> Audit	Can #:	Cyl #:	Tst #:	Comment:
<input type="checkbox"/> MDL	Can #:	Cyl #:	Tst #:	Comment:
<input type="checkbox"/> Other	Can #:	Cyl #:	Tst #:	Comment:

QC Sample Check
Update information below as needed. Select the checkboxes as appropriate when CVS, LCS, Blank, or Audit Failures occur. **Please note that Propane and Benzene % Recovery boxes may be used even if the QC did not fail.

<input type="checkbox"/> CVS Failed	Propane: %	Benzene: %	Comment:
<input type="checkbox"/> LCS Failed	Propane: %	Benzene: %	Comment:
<input type="checkbox"/> Blank Failed			Comment:
<input type="checkbox"/> Audit Failed			Comment:
<input type="checkbox"/> Other			Comment:

Hardware Checks
Update information as needed

<input type="checkbox"/> Nafion	<input type="checkbox"/> Plot Column
<input type="checkbox"/> BP Column	<input type="checkbox"/> Trap
<input type="checkbox"/> TOC	<input type="checkbox"/> CPU date t
<input type="checkbox"/> Other	

Preventive Maintenance (PM) Please document start date/time and end date/time of PM

<input type="checkbox"/> Annual Routine PM	Start date: mm/dd/yyyy 00-23	End date: mm/dd/yyyy 00-23	Comment:
--	------------------------------	----------------------------	----------

Calibration: Update the information below as needed; RF=Response Factor; RSD=Relative Standard Deviation

<input type="checkbox"/> Singlepoint	PLOT RF:	BP RF:	Updated Plot Method with RF	Comment:
<input type="checkbox"/> Multipoint Curve	PLOT RF:	BP RF:	Updated Plot Method with RF	Comment:
	% RSD:	% RSD:	Updated BP Method with RF	Comment:

Blender Settings:
Update the information below as needed. Select the Reset checkboxes as appropriate when SPAN (or Zero) is reset. **Please note that the "To" input box may be used to record current pressure even if the SPAN or ZERO was not reset.

<input type="checkbox"/> Reset SPAN	From: psig	To (current): psig	Comment:
<input type="checkbox"/> Reset ZERO	From: psig	To (current): psig	Comment:
<input type="checkbox"/> New Blender Serial#		Comment:	
<input type="checkbox"/> New Blend Ratio		Comment:	
<input type="checkbox"/> Other		Comment:	

Method & Sequence Optimization
Update the information below when a method/sequence is created or modified

<input type="checkbox"/> New Sequence	Start dt: mm/dd/yyyy 00-23	End dt: mm/dd/yyyy 00-23
<input type="checkbox"/> Plot Method	Comment:	
<input type="checkbox"/> BP-1 Method	Comment:	
<input type="checkbox"/> Modified Sequence	Comment:	
<input type="checkbox"/> Other:	Comment:	

Batch Reprocessing
Update the information below when data is batch reprocessed

<input type="checkbox"/> Plot Data	Start date: mm/dd/yyyy 00-23	End date: mm/dd/yyyy 00-23	Comment:
<input type="checkbox"/> BP Data	Start date: mm/dd/yyyy 00-23	End date: mm/dd/yyyy 00-23	Comment:


Other Comments

Electronic Logbooks

AutoGC Operator Log Report

Site	Date	Username	Source	Log Entry	Comment
2D	07/29/2016 07:22	lrice	Remote Visit	Pressures, Flows, and other Checks <i>Verified Sample Timing:</i> [system is on time with files ending @(:32) and correct row (175)] QC Sample Check <i>CVS:</i> Propane=100, Benzene=99 [2D*CG29C]	Methods and Blanks look good
2D	07/28/2016 10:40	lrice	Site Visit	Pressures, Flows, and other Checks <i>Helium:</i> 1100 psig [@60] <i>Hydrogen 1:</i> 2300 psig [@65] <i>Zero Air:</i> 2300 psig [@85] <i>CVS:</i> 25 psig [37806] <i>LCS:</i> 16 psig [35919] <i>RTS:</i> 30 psig [14440] <i>GC Cut Time:</i> 10.00 min <i>Station Temp:</i> 73 °F [High 79 Low 70 Humidity: 48% High 48% Low 25%] <i>Nafton Flow:</i> .25 slp <i>TOC Flowmeter:</i> .49 slp [@80 psi] <i>TD Sample Pump:</i> 15 ml/min <i>TD Column Pressure:</i> 42.4 psig <i>GC Midpoint Pressure:</i> 18.5 psig <i>Other:</i> 85 psi [zero air output] <i>Verified Sample Timing:</i> [system is on time with files ending @(:32) and correct row (155)] Canister Installation <i>RTS:</i> 12240, CC-328584 [Removed Can #12247] Blender Settings <i>SPAN Not Reset:</i> to 5.0 psig <i>ZERO Not Reset:</i> to 40.0 psig Method & Sequence Optimization <i>Modified Sequence:</i> [setup RTS (spike) in sequence]	Methods look good Router signal: 4 lights All green lights on router Magnihelic pressure: .55 RMS noise BP-1: 7.96 Plot: 8.84 A-Fid: .22 mV B-Fid: .22 mV GC inlet pressure: 44.5 psi TOC: green light Variac: 35 volts Manifold tube heater is working Compressor water trap drained Output: 95 psi Tank cutoff: 135 psi Auto drain is working All catalyst feel warm Setup RTS (spike) in sequence
2D	07/27/2016 08:32	lrice	Remote Visit	Pressures, Flows, and other Checks <i>Verified Sample Timing:</i> [system is on time with files ending @(:32) and correct row (128)] QC Sample Check <i>CVS:</i> Propane=102, Benzene=102 [2D*CG27C]	Methods and Blanks look good

Electronic Logbooks



The screenshot shows a web browser window with the title "Markes Logbook". The address bar shows "airqualitytexas@gmail.com". The browser's menu bar includes "File", "Edit", "View", "Insert", "Format", "Tools", "Table", "Add-ons", "Help", and "Insert DTS". The toolbar shows "100%", "Normal text", "Arial", "11", and "More". The main content area displays three log entries:

06/09/2016 14:15 /[gcoldwell](#)
Rerouted output of zero air to feed blender and gc detector flows
Measure relative humidity off the blender output after the zero air source was rerouted:
Measure as follows :
20 psi = 37.3
40 psi 31.0
60 psi 27.0
Checked for moisture in dropout and also ballast tank all looks good so far .

06/10/2016 15:00 [gcoldwel](#)
Performed 3 pt curve using 16 component CC-344433
results :Plot 6140 3 % rsd BP 5979 4 % rsd
Enter into methods I HOPE
When modifying sequence saved ezchrom sequence all appeared to be correct but when files
were written ezchrom did not change the file extension to reflect the file changes that were made
.



06/14/2016 12:15/[gcoldwell](#)
HE = 550 psi
System has lost power td was stuck trap purge and software was not responding so [cyc:ed](#) power
on TD and all looks good now
Corrected rows on both sequences and set up another round of QC

Tools to Monitor Performance

- ▶ Automation of routine quality control samples
- ▶ Remote Access via broadband connections
 - ▶ Applications for remote control (remote desktop)
 - ▶ FTP - polling and file transfer automation
- ▶ Email alerts
 - ▶ Power failures
 - ▶ No data (polling)

Tools to Monitor Performance: Control Charts


LEADS[®] Leading Environmental Analysis & Display System




Tools to Monitor Performance: Control Charts



Daily CVS



Weekly LCS



Tools to Monitor Performance QuickLook Email Report

Galena Park																			
CONCENTRATION PPB-C										CURRENT SEQUENCE/SEQ NAME					SEQUENCE (DAYS IN USE)				
150515-53.seq																			
file	date	time	methane	ethane	ethylene	propene	propylene	ic4	nc6	acetylene	1,2-butadiene	1-butene	1,3-butadiene	1,2-pentene	2-methylbutene				
Ce30a	30-May	152	0.6	16.3				25.4	86				1.3	26.9	23.3				
Be30b	30-May	232																	
Se30c	30-May	232	21.9	28.5	3.2	8.3	0.8	5.3	8.5	0.9	0.4	1.2	16.9	13.1					
Se30d	30-May	432	20.3	23.7	3.6	8.8	2.3	6.7	9.7	1.0		1.0	15.9	12.3					
Se30e	30-May	532	21.1	20.5	2.3	7.0	2.4	5.3	7.7	0.5	0.0	1.3	0.8	15.0	10.5				
Se30f	30-May	632	26.1	25.6	3.9	12.7	7.4	18.4	11.4	1.6	0.0	2.5	0.9	13.5	6.5				
Se30g	30-May	732	22.0	14.9	2.9	8.6	2.9	5.6	8.7	1.0	0.5	1.6	0.9	19.4	12.6				
Se30h	30-May	832	19.2	3.9	10.2	4.3	10.2	4.3	10.2	1.8	2.7	1.1	0.8	15.5	5.6				
Se30i	30-May	932	19.5	6.5	1.4	4.6	1.8	2.0	2.0	0.6			5.9	2.2					
Se30j	30-May	1032	20.2	6.5	1.4	2.0	1.9	1.0	2.8	0.5			3.9	1.4					
Se30k	30-May	1132	19.4	6.7	1.8	8.5				1.8	0.5		2.6	1.2					
Se30l	30-May	1232	19.3	6.4	1.3	5.4	0.6	2.9	8.1	0.6			0.6	20.9	13.9				
Se30m	30-May	1332	19.4	8.7	1.8	8.5		10.4	26.9	0.9	1.0	1.4	1.0	4.6	188.7	105.1			
Se30n	30-May	1432	19.1	7.0	1.0	2.9	0.9	3.6	17.9	0.6	0.4		4.3	188.7	129.4				
Se30o	30-May	1532	19.9	12.7	2.1	7.7			49.8	0.5	1.2		10.4	462.4	312.1				
Se30p	30-May	1632	18.8	20.6	3.5	12.3	4.7	15.2	103.8	0.8	2.1	1.4	2.0	26.4	1220.6				
Se30q	30-May	1732	18.7	15.2	3.9	2.7	1.7	1.3	3.2	0.4			11.1	8.7					
Se30r	30-May	1832	19.3	3.4	18.3	4.5	11.7	37.1	1.0	0.8	1.4	0.6	3.0	94.7	67.1				
Se30s	30-May	1932	20.3	18.6	2.5	17.2	2.2	17.3	58.4	1.5	1.2	1.4	1.0	4.7	158.9				
Se30t	30-May	2032	19.5	10.0	10.0	10.2	4.5	44.2	122.4	1.9	1.6	2.0	1.3	5.8	198.4				
Se30u	30-May	2132	19.6	14.2	2.6	9.4	2.9	3.2	6.9	0.9			0.5	9.5	4.7				
Se30v	30-May	2232																	
Se30w	30-May	2332	1.9	1.5	2.6	1.4	2.3	4.1			2.4	2.9	3.2	2.0	4.1				
Se30x	30-May	32	25.3	11.4	1.5	6.7	1.2	2.9	4.6	0.7			0.4	0.5	6.7				

file	nc6	mcyc5	2,4dmc5	benzene	cyc6	2-mc6
Ce30a	34.0			32.5	1.0	
Be30b						
Se30c	5.3	2.0	0.5	2.1	1.2	1.0
Se30d	4.0	2.0	0.8	2.0	1.3	1.0
Se30e	3.9	2.0	0.6	1.5	1.2	1.2
Se30f	4.9	3.0	0.9	1.4	1.8	1.7
Se30g	3.7	2.3	1.1	1.4	1.9	1.2
Se30h	3.4	2.1	0.6	1.3	1.6	1.2
Se30i	1.2	0.8		0.6	0.6	0.6
Se30j	0.8	0.6		0.5		0.4
Se30k	0.7					
Se30l	2.4	0.9	0.8	0.7	0.5	0.5
Se30m	13.8	7.0	8.1	3.7	3.3	2.7
Se30n	11.7	5.5	10.0	3.0	2.8	1.3
Se30o	28.2	12.9	22.7	10.3	6.5	2.7
Se30p	76.5	33.7	58.8	21.9	24.8	7.2
Se30q	1.4	0.6		0.7		0.5
Se30r	15.2	5.7	5.1	3.2	3.4	2.7
Se30s	19.6	7.8	9.1	5.0	4.0	3.8
Se30t	30.6	12.0	8.7	5.5	7.2	7.7
Se30u	2.5	1.4		1.1	0.7	1.0
Se30v	34.1			32.0		
Se30w	3.2	2.4	3.9	2.9	4.3	2.0
Se30x	2.4	0.9	0.7	0.8	0.5	0.7

DAILY															
5/30/2015					5/30/2015					5/30/2015					
DATE	5/30/2015	DATE	5/30/2015	DATE	5/30/2015	DATE	5/30/2015	DATE	5/30/2015	DATE	5/30/2015	DATE	5/30/2015	DATE	5/30/2015
Carbon	Conc	%Recov	Carbon	Conc	%Recov	Carbon	Conc	%Recov	Carbon	Conc	%Recov	Carbon	Conc	%Recov	Carbon
Ethane	2.00	1.06	11.99	92.8		Ethane	2.00	1.03	11.99	92.8		Ethane	2.00	1.03	11.99
Propane	3.00	1.01	16.28	16.27	99.9	Propane	3.00	1.02	16.77	17.49	104.3	Propane	3.00	1.02	16.77
m-Butane	4.00	1.02	21.92	25.98	119.8	m-Butane	4.00	1.04	22.80	25.71	104.0	m-Butane	4.00	1.04	22.80
Acetylene	2.00	1.06	11.99	9.62	84.5	Acetylene	2.00	1.07	11.73	11.01	93.9	Acetylene	2.00	1.07	11.73
1-butene	5.00	1.09	26.99	30.95	116.4	1-butene	5.00	1.08	28.99	30.02	101.5	1-butene	5.00	1.08	28.99
1,3-butadiene	4.00	1.03	22.13	23.34	105.5	1,3-butadiene	4.00	0.99	21.70	21.10	97.2	1,3-butadiene	4.00	0.99	21.70
2-methylbutane	6.00	1.00	32.23	34.95	108.4	2-methylbutane	6.00	1.01	33.21	35.76	107.7	2-methylbutane	6.00	1.01	33.21
PLUT BLANK	BLANK	B			2.43	PLUT BLANK	BLANK	B			2.43	PLUT BLANK	BLANK	B	
Hexane	6.00	1.03	33.20	34.05	102.6	Hexane	6.00	1.00	32.88	34.09	103.7	Hexane	6.00	1.01	33.21
Benzene	6.00	1.03	33.23	32.47	97.7	Benzene	6.00	1.01	33.21	32.03	96.5	Benzene	6.00	1.01	33.21
Toluene	7.00	1.05	35.84	30.8		Toluene	7.00	1.01	38.74	27.28	96.2	Toluene	7.00	1.01	38.74
m-Xylene	8.00	1.06	45.55	40.95	89.9	m-Xylene	8.00	0.97	42.52	43.28	101.8	m-Xylene	8.00	0.97	42.52
p-Propylbenzene	9.00	1.02	48.31	43.9	88.1	p-Propylbenzene	9.00	1.00	46.32	47.31	98.2	p-Propylbenzene	9.00	1.00	46.32
1,2,4-Trimethylbenzene	9.00	1.04	50.38	41.42	82.4	1,2,4-Trimethylbenzene	9.00	1.02	50.31	42.62	86.5	1,2,4-Trimethylbenzene	9.00	1.02	50.31
BP1 BLANK	BLANK	B			3.64	BP1 BLANK	BLANK	B			3.64	BP1 BLANK	BLANK	B	
BLEND RATIO	0.0037					BLEND RATIO	0.0037					BLEND RATIO	0.0037		
RF PLUT	5144					RF PLUT	5144					RF PLUT	5144		
RF BP1	4723					RF BP1	4723					RF BP1	4723		
DATE	5/27/2015					DATE	5/27/2015					DATE	5/27/2015		
BP1	6.00					BP1	6.00					BP1	6.00		

** NOTE: nc6c levels < 0.4 are not listed



Carol J. Meyer
carol@orsat.com



Orsat, LLC
www.orsat.com

FIFTY SIX
NON-METHANE
HYDROCARBONS
TWENTY FOUR
HOURS A DAY
SEVEN
DAYS A WEEK