

Independent Laboratory Validation for the Determination of Residues of Nitrapyrin in Ground,
Surface, and Drinking Water by Gas Chromatography with Mass Spectrometry and 6-
Chloropicolinic Acid (6-CPA) in Ground, Surface, and Drinking Water by Liquid
Chromatography with Tandem Mass Spectrometry

INTRODUCTION

Scope

Separate methods are applicable for the quantitative determination of residues of nitrapyrin and its metabolite 6-chloropicolinic acid (6-CPA) in ground, surface, and drinking water. The nitrapyrin method was validated over the concentration range of 0.05-1 µg/L with a validated limit of quantitation (LOQ) of 0.05 µg/L. The 6-CPA method was validated over the concentration range of 0.05-1 µg/L with a validated LOQ of 0.05 µg/L. Common names, chemical names, and molecular formulas for the analytes are given in Table 1.

This study was conducted to fulfill data requirements outlined in the EPA Residue Chemistry Test Guidelines, OCSPP 850.6100 (1). The independent validation also complies with the requirements of EU Council Regulation (EC) No. 1107/2009 with particular regard to SANCO/825/00 rev.8.1 (2).

Method Principle

The residue of nitrapyrin is extracted from water samples with two portions of 1:1 (v/v) hexane:toluene. The combined hexane:toluene extracts are evaporated to less than 1 mL and cleaned up on a deactivated silica solid-phase extraction (SPE) column. The SPE eluate containing nitrapyrin is analyzed by gas chromatography with electron-impact mass spectrometry (GC/EI-MS).

The residue of 6-CPA is extracted from water samples using C18 solid phase extraction (SPE) cartridges. The residue is eluted from the SPE cartridges with acetonitrile, the sample is brought to dryness, then reconstituted with 25/75 (v/v) methanol/DI water. The final sample is analyzed by liquid chromatography coupled with negative-ion electrospray tandem mass spectrometry (LC-MS/MS).

Test Substances/Reference Compounds/Analytical Standards

Test Substance	TSN	Percent Purity	Recertification Date	Reference
Nitrapyrin	AGR023966	99.9	11 July 2015	FAPC11-285444
6-CPA	TSN301485	99	07 July 2015	FAPC13-000236

The Certificates of Analysis for the test substances can be found in Figures 1 and 2. The above standards were obtained from Dow AgroSciences LLC.

EXPERIMENTAL

Sample Origin, Numbering, Preparation, Storage, and Characterization

Untreated control samples were obtained by Battelle from local sources. All samples were tracked in the Battelle-Norwell Laboratory Information Management System (LIMS) database. Unique reagent, sample and batch numbers were assigned to track reagents and samples during receipt, storage, and analysis following Battelle SOPs 6-007 "Chemistry Laboratory Sample Identification" and 5-015 "Solvent/Reagent Inventory and Contaminant Residue Checks."

Samples were prepared for nitrapyrin analysis following the method described in Appendix 1 "Determination of Nitrapyrin in Water by Gas Chromatography-Mass Spectrometry", which represents the method provided by the sponsor (3) with slight modifications to account for differences in laboratory equipment and changes glassware substitutions and concentration of fortification solutions. Samples were prepared for 6-CPA analysis following the method described in Appendix 2 "Determination of 6-Chloropicolinic Acid (6-CPA) in Water by Liquid Chromatography with Tandem Mass Spectrometry Detection," which represents the method provided by the sponsor (4) with slight modifications to account for differences in differences in laboratory equipment or instrumentation. Changes to both methods provided by the sponsor are described in the "Changes to Method" section.

Samples were stored refrigerated at approximately 4 °C after the time of sampling and during the course of the independent laboratory validation study, except when they were removed for taking aliquots for sample preparation and analysis.

Samples were characterized in compliance with 40 CFR Part 160. The following table lists critical water quality parameters. Complete source documentation and water quality reports are included in the study file.

Water Type	pH	Hardness (mg equiv. CaCO ₃ /L)	Total Suspended Solids (ppm)	Total Organic Carbon (ppm)	Dissolved Organic Carbon (ppm)
Drinking Water (Tap)	8.0	70	10	0.8	0.1
Ground Water (Monitoring Well)	6.9	32	12	0.5	0.5
Surface Water (Pond)	6.6	37	6	2.9	2.5

Instrumentation

For the detection of nitrapyrin, an Agilent Technologies 7890A GC System, equipped with a 7693 Autosampler and 5975C inert MSD was used. The following GC/EI-MS parameters were used to determine the concentration of nitrapyrin in sample extract solutions.

GC System: Agilent 7890A

Injector: Agilent Model 7693 with 150 sample capacity autosampler tray

Detector: Agilent Model 5975C mass spectrometer

Column: Agilent DB-5MS (30 m x 0.25 mm ID x 1.00 µm film)

Oven Method:

Column

130 °C for 1.0 minute
 130 °C to 210 °C at 10 °C/minute
 210 °C to 280 °C at 15 °C/minute
 280 °C for 2.0 minute

Carrier Gas:

Constant Flow

Helium

1 mL/min

Injection:

Injector Temperature

Splitless

180 °C

Injection Volume

2 µL

Purge flow

40.0 mL/min

Detection Method:

Source Temperature

Positive-ion Electron-Impact ionization mode

250 °C

Quad Temperature

180 °C

Ions Monitored:

Quantitation

Ion

Dwell

Confirmation (Primary)

m/z 194

50 msec

Confirmation (Secondary)

m/z 196

50 msec

m/z 198

50 msec

The following LC-MS/MS parameters are used to determine 6-CPA concentrations.

LC System: Agilent 1260 SL HPLC system (including a vacuum solvent degasser, binary HPLC pump, and column oven) and CTC Analytics HTC-Pal Autosampler

MS Detector: Applied Biosystems MDS Sciex API 5500 linear ion trap MS/MS system with TurboIonspray (ESI) source

HPLC Conditions:

Column

Waters XSelect® HSS T3 3.5 µm, 3.0 x 100 mm

Column Temperature

40°C

Injection Volume

10 µL for batches 14-0606, 15-0003, 14-0618,
 15-0002, 14-0625 and 15-0001

Mobile Phase

5 µL for batches 15-0016, 15-0027 and 15-0028

A: 0.01% (v/v) Acetic Acid

B: Methanol

HPLC Gradient Program:

Time (minutes)	%A	%B	Curve	Flow (mL/minute)
0.00	95	5	1	0.40
0.10	95	5	1	0.40
4.00	30	70	6	0.40
4.50	30	70	1	0.40
4.51	5	95	1	0.40
6.50	5	95	1	0.40
6.51	95	5	1	0.40
8.50	95	5	1	0.40

Divert flow to waste from 0 to 0.5 minutes. Direct flow to the MS detector from 0.5 to 8.5 minutes. Divert flow to waste again at end of run.

Ion Source Conditions:

Source temperature	350 °C
Gas supply (GS 1):	40
Gas supply (GS 2)	60
Curtain gas:	30
CAD gas:	Medium
ESI Polarity:	Negative
Entrance potential:	-10 V
Ionspray Voltage	-4000 V

MS/MS Conditions for 6-CPA:

Quantification MS/MS transition	155.8/111.9
Collision energy (CE):	-13 V
Cell exit potential (CXP):	-13 V
Dwell time:	100 ms
Declustering potential (DP):	-50 V
Confirmation MS/MS transition	157.9/114.0
Collision energy (CE):	-13 V
Cell exit potential (CXP):	-13 V
Dwell time:	100 ms
Declustering potential (DP):	-50 V

Note that samples analyzed before 13 January 2015, utilized slightly different HPLC conditions and samples analyzed before 08 January 2015 has slightly different MS conditions. These conditions are described below:

HPLC Conditions:

Column	Water HSS T3 1.8 µm, 2.1 x 100 mm
Column Temperature	Same as above
Injection Volume	Same as above
Mobile Phase	A: Same as above

B: Same as above

HPLC Gradient Program:

Time (minutes)	%A	%B	Curve	Flow (mL/minute)
0.00	95	5	1	0.20
0.10	95	5	1	0.20
4.00	30	70	6	0.20
4.50	30	70	1	0.20
4.51	5	95	1	0.20
6.50	5	95	1	0.20
6.51	95	5	1	0.20
8.50	95	5	1	0.20

Flow diversion same as above.

Ion Source Conditions:

Source temperature	Same as above
Gas supply (GS 1):	Same as above
Gas supply (GS 2)	Same as above
Curtain gas:	Same as above
CAD gas:	Same as above
ESI Polarity:	Same as above
Entrance potential:	Same as above
Ionspray Voltage	Same as above

MS/MS Conditions for 6-CPA:

Quantification MS/MS transition	Same as above
Collision energy (CE):	Same as above
Cell exit potential (CXP):	Same as above
Dwell time:	Same as above
Declustering potential (DP):	Same as above
Confirmation MS/MS transition	Same as above
Collision energy (CE):	-25 V
Cell exit potential (CXP):	Same as above
Dwell time:	Same as above
Declustering potential (DP):	Same as above

Calculation of Standard Calibration Curve

Calculation of a standard curve for nitrpyrin by GC-MS begins with the injection of a series of matrix-matched calibration standards described in Tables 2-4 and acquisition of peak areas for the following ions monitored in selective ion monitoring (SIM) mode.

Nitrapyrin	<i>m/z</i> 194 (quantitative)
	<i>m/z</i> 196 (confirmatory 1)
	<i>m/z</i> 198 (confirmatory 2)

Calculation of a standard curve for 6-CPA by LC-MS/MS begins with the injection of a series of matrix-matched calibration standards described in Tables 5-10 and acquisition of peak areas for the following transitions monitored in selected reaction monitoring (SRM) mode.

6-CPA	<i>m/z</i> Q1/Q3 155.8/111.9 (quantitative)
	<i>m/z</i> Q1/Q3 157.9/114.0 (confirmatory)

In order to generate a standard curve, the analyte concentration was plotted on the abscissa (x-axis) and the respective peak area on the ordinate (y-axis) in SigmaPlot for nitrapyrin data, or using MultiQuant for 6-CPA data. Using regression analysis, the equation for the curve fit was determined with respect to the abscissa using 1/x fit weight. Refer to Figures 3-11 for calibration plots for nitrapyrin and to Figure 18 for example calculations. Refer to Figures 12-17 for calibration plots for 6-CPA and to Figure 19 for example calculations.

Confirmation of Residue Identity

These methods are specific for the determination of nitrapyrin and 6-CPA by virtue of the chromatographic separations and selective detection systems used. To demonstrate further confirmation, two additional SIM ions were monitored for nitrapyrin and an additional MS/MS ion transition was monitored for 6-CPA.

Inject the series of calibration standards and acquire peak area for the following analytes.

Nitrapyrin	<i>m/z</i> 194 (quantitative)
	<i>m/z</i> 196 (confirmatory 1)
	<i>m/z</i> 198 (confirmatory 2)
6-CPA	<i>m/z</i> Q1/Q3 155.8/111.9 (quantitative)
	<i>m/z</i> Q1/Q3 157.9/114.0 (confirmatory)

Statistical Treatment of Data

Statistical treatment of data included but was not limited to the calculation of regression equations, correlation coefficients (*r*) for describing the linearity of calibration curves, and means, standard deviations, and relative standard deviations of the results for the fortified recovery samples.

Matrix Effects

Matrix effects were evaluated by comparing the response of the analyte fortified in a control extract after processing (for each matrix type) to the response of the analyte fortified in solvent. For the nitrapyrin method, the solvent consisted of toluene and for the 6-CPA method the solvent consisted of 25/75 (v/v) methanol/DI water. The calculation for the matrix effect is as follows:

$$\text{matrix effect} = \left(\frac{\text{peak area}_{\text{matrix matched calibration standard}}}{\text{peak area}_{\text{solvent calibration standard}}} - 1 \right) \times 100\%$$

A negative value for the matrix effect indicates matrix suppression, and a positive value for matrix effect indicates matrix enhancement.

The experimental details regarding determination of the matrix effects were recorded in the raw data file.

For the nitrapyrin method, the solvent calibration standards represent the fortifications of nitrapyrin in solvent (toluene). The matrix matched calibration samples were prepared by fortifying the eluent from the silica gel column elution for control samples with an appropriate amount of fortification standard (between 18 and 150 μL of solutions IF54 or IF55 as described in the study file) prior to the eluent being brought up to 3 mL with toluene. The matrix effects were determined at each calibration level by comparing the peak area of the matrix matched calibration standard to that of the peak area in solvent calibration standard. Matrix effects for the quantitative and two confirmatory SIM ions for nitrapyrin were calculated, and the results are summarized in Table 26.

For 6-CPA, the solvent calibration standards represent the fortifications of 6-CPA in solvent (25/75 (v/v) methanol/DI water). The matrix matched calibration samples were prepared by reconstituting the residue of an extracted control with 0.5 mL of solvent calibration standard. The matrix effects were determined at each calibration level by comparing the peak area of the matrix matched calibrations standard (i.e. the Peak Area_{matrix matched calibrations standard} value in the equation above) to that of the peak area in solvent calibration standard (i.e. the Peak Area_{solvent calibrations standard} value in the equation above).

Matrix effects for the quantitative and confirmatory transitions for 6-CPA were calculated, and the results are summarized in Tables 27-32.

Communication

There was no communication between the ILV laboratory and the method developers or others familiar with the method for the nitrapyrin component of this ILV.

There was communication between the ILV laboratory and the method developers and those at Dow AgroSciences familiar with the method during the 6-CPA component of this ILV. This communication revolved around understanding the reasons for the failed batches of samples.

Communication included sharing HPLC and MS/MS conditions, chromatograms and specifics of the blow-down step. Additionally, the validating laboratory was asked about alternative possible mechanisms for loss, such as absorption into plastic. Communication was initiated on January 6, 2015 after unsuccessful results were obtained. As a result of this communication, a different analytical column was purchased, more suitable to the Agilent 1260 HPLC system in use resulting in improved chromatography, and a series of tests was conducted to understand where loss was occurring. In preparing for this communication, it was noted that the collision energy for the confirmation transition was not similar to that of the quantitation transition, and adjusted the energy accordingly. A summary of this communication is provided in Appendix IV, and all communication is included in the study file. The results of the accuracy checks of the fortification solutions are included in Appendix III.



Changes to the Method

For the nitrpyrin component of this ILV, the only changes to the method between the ILV laboratory and the method developers were minor glassware substitutions (250 mL separatory funnels rather than 125 mL, collecting column eluate in 4 mL glass vials rather than 15 mL, measuring the final volume of the extract by comparing to a reference vial rather than using a vial with graduations marked on it) and changes to the concentrations of fortification solutions (note: the fortification amount was not changed). The concentrations of all standards are provided in the study file.

For the 6-CPA component of this ILV, minor deviations were made to account for differences in laboratory equipment or instrumentation:

- Polypropylene centrifuge tubes were used to collect eluate and blow-sample down for the drinking water batches 14-0606, 14-0616, 14-0627; groundwater batches 14-0608, 14-0618, 14-0626, and surface water 14-0607, 14-0617, 14-0625)
 - For the remaining sets (drinking water batches 15-0003 and 15-0016, groundwater batches 15-0002 and 15-0027, and surface water batches 15-0001 and 15-0028), glass test tubes were used
- The volume of SPE eluate was reduced using a nitrogen evaporator with a water bath, rather than a heating block, set at 30-35°C.
- After poor chromatography using the prescribed UPLC column (Waters HSS T3 1.8 µm, 2.1 x 100 mm) likely due to the fact that the HPLC system (Agilent 1260) was not a true UHPLC system, the project team switched to a similar column more amenable to HPLC flow rates and pressures (Waters HSS T3 3.5 µm, 3.0 x 100 mm).
 - Additionally, a gradient more amenable to this HPLC column was employed (see Experimental Section)

Table 1. Identities and Structures of Nitrapyrin and 6-CPA

Common Name	Structural Formula and Chemical Name
Nitrapyrin Molecular Formula: C ₆ H ₃ Cl ₄ N CAS Number: 1929-82-4	 2-chloro-6-(trichloromethyl)pyridine
6-CPA Molecular Formula: C ₆ H ₄ ClNO ₂ CAS Number: 4684-94-0	 6-chloro-2-pyridinecarboxylic acid

Calculation of Nitrapyrin in Sample Extract

Generate a standard curve by plotting the analyte concentration on the abscissa (x-axis) and the respective peak area on the ordinate (y-axis) in SigmaPlot. Using regression analysis, determine the equation for the curve with respect to the abscissa.

$$f = y_0 + a \cdot x$$

$$C_{\text{end}} = (f - y_0) / a$$

where

x = concentration of nitrapyrin in extract

f = analyte peak area

C_{end} = analyte concentration in extract (ng/mL)

y_0 = y-intercept of regression (ng/mL)

a = slope

Determine the concentration in the extract of each recovery sample by substituting the corresponding analyte peak area, y-intercept of regression and slope into the above equation and solving for the concentration. Determination of extract concentration in each sample proceeds similarly; therefore, only a single example calculation will be presented here. For example, using linear regression with 1/x weighting, for nitrapyrin (m/z 194) in drinking water (batch 14-0599, file name: Q0112; sample ID: CE042LOQ-AG(0)) the calculation is as follows (note: the chromatogram is shown in Figure 28):

$$x = (188 - 9.3315) / 128.4935$$

$$x = 1.39 \text{ ng/mL}$$

Calculation of Nitrapyrin in Original Sample

Convert the extract concentration (ng/mL) of nitrapyrin to the concentration in the original sample ($\mu\text{g/L}$).

The following equation was used to calculate the concentration of nitrapyrin (R, in $\mu\text{g/L}$) in the original fortified samples:

$$R = x \cdot \text{extract volume} / \text{sample volume} \cdot \text{mass unit conversion} / \text{volume unit conversion}$$

Determination of the nitrapyrin concentrations in all original sample proceeds similarly; therefore, only a single example calculation will be presented here. For example, for nitrapyrin

(*m/z* 194) in drinking water (batch 14-0599, file name: Q0112; sample ID: CE042LOQ-AG(0)) the required information is as follows:

Extract volume = 3.00 mL

Sample volume = 100 mL

Mass unit conversion (1000 ng = 1 µg)

Volume unit conversion (1000 mL = 1 L)

Therefore

$$R = 1.39 \text{ (ng/mL)} * 3.00 \text{ (mL)} / 100 \text{ (mL)} / 1000 \text{ (ng/}\mu\text{g)} * 1000 \text{ (mL/L)}$$

$$R = 0.0417 \text{ }\mu\text{g/L}$$

Calculation of Percent Recovery

Determine the percent recovery by dividing the concentration in the original sample (µg/L) for each recovery sample by the fortified concentration added as follows:

$$\text{Recovery} = R / \text{fortified concentration} * 100\%$$

Determination of all nitrapyrin recoveries proceed similarly; therefore, only a single example calculation will be presented here. For example, for nitrapyrin (*m/z* 194) in drinking water (batch 14-0599, file name: Q0112; sample ID: CE042LOQ-AG(0)) the calculation is as follows:

$$R = 0.0417 \text{ }\mu\text{g/L}$$

$$\text{Fortified concentration} = 0.0500 \text{ }\mu\text{g/L}$$

Therefore

$$\text{Recovery} = 0.0486 \text{ (}\mu\text{g/L)} / 0.0500 \text{ (}\mu\text{g/L)} * 100\%$$

$$\text{Recovery} = 83.4\%$$

Figure 19. Example Calculation for the Quantitative Determination of 6-CPA in Drinking Water, Groundwater and Surface Water

Calculation of 6-CPA in Sample Extract

Generate a standard curve by plotting the analyte concentration on the abscissa (x-axis) and the respective peak area on the ordinate (y-axis) in MultiQuant. Using regression analysis, determine the equation for the curve with respect to the abscissa.

$$y = y_0 + a*x$$

$$C_{\text{end}} = (f - y_0)/a$$

where

x = concentration of 6-CPA in extract

y = analyte peak area

C_{end} = analyte concentration in extract (ng/mL)

y_0 = y-intercept of regression (ng/mL)

a = slope

Determine the concentration in the extract of each recovery sample by substituting the corresponding analyte peak area, y-intercept of regression and slope into the above equation and solving for the concentration. Determination of extract concentration in each sample proceeds similarly; therefore, only a single example calculation will be presented here. For example, using linear regression with 1/x weighting, for 6-CPA (*m/z* 155.8/111.9) in drinking water (batch 15-0016, file name: 100057852-15-0016.wiff; sample ID: CE571LOQ-AG(0)) the required information is as follows (note: the chromatogram is shown in Figure 28):

$$x = (70612 - 20401.97811) / 25833.33768$$

$$x = 1.944 \text{ ng/mL}$$

Calculation of 6-CPA in Original Sample

Convert the extract concentration (ng/mL) of 6-CPA to the concentration in the original sample ($\mu\text{g/L}$).

The following equation was used to calculate the concentration of 6-CPA (R, in $\mu\text{g/L}$) in the original fortified samples:

$$R = x * \text{extract volume} / \text{sample volume} * \text{mass unit conversion} / \text{volume unit conversion}$$

Determination of the 6-CPA concentrations in all original samples proceeds similarly; therefore,

only a single example calculation will be presented here. For example, for 6-CPA (*m/z* 155.8/111.9) in drinking water (batch 15-0016), file name: 100057852-15-0016.wiff; sample ID: CE571LOQ-AG(0)) the required information is as follows:

Extract volume = 0.500 mL

Sample volume = 20.0 mL

Mass unit conversion (1000 ng = 1 µg)

Volume unit conversion (1000 mL = 1 L)

Therefore

$$R = 1.944 \text{ (ng/mL)} * 0.500 \text{ (mL)} / 20.0 \text{ (mL)} / 1000 \text{ (ng/}\mu\text{g)} * 1000 \text{ (mL/L)}$$

$$R = 0.0486 \text{ }\mu\text{g/L}$$

Determine the percent recovery by dividing the concentration in the original sample (µg/L) for each recovery sample by the fortified concentration added as follows:

$$\text{Recovery} = R / \text{fortified concentration} \times 100\%$$

Determination of all nitrapyrin and 6-CPA recoveries proceed similarly; therefore, only a single example calculation will be presented here. For example, for 6-CPA (*m/z* 155.8/111.9) in drinking water (batch 15-0016), file name: 100057852-15-0016.wiff; sample ID: CE571LOQ-AG(0)) the required information is as follows:

$$R = 0.0486 \text{ }\mu\text{g/L}$$

$$\text{Fortified concentration} = 0.0505 \text{ }\mu\text{g/L}$$

Therefore

$$\text{Recovery} = 0.0486 \text{ (}\mu\text{g/L)} / 0.0505 \text{ (}\mu\text{g/L)} \times 100\%$$

$$\text{Recovery} = 96.2\%$$

Determination of Nitrpyrin in Water by Gas Chromatography-Mass Spectrometry

Note: the method described herein was modified from EPL-BAS Method No. 205G907A to account for differences in differences in laboratory equipment or instrumentation. Changes to both methods provided by the sponsor are described in the "Changes to Method" section.

Safety Precautions

Each analyst must be acquainted with the potential hazards of the equipment, reagents, products, solvents, and procedures used in this method before commencing laboratory work. Sources of information include operation manuals, material safety data sheets, literature, and other related data. Safety information should be obtained from the supplier. Disposal of waste materials, reagents, reactants, and solvents must be in compliance applicable governmental requirements.

Acetonitrile, hexane, and toluene are flammable and should be used in well-ventilated areas away from ignition sources. Sulfuric acid is corrosive and can cause severe burns. It is imperative that proper eye and personal protection equipment be worn when handling these reagents.

Laboratory Equipment

Balance, analytical, Mettler AT201

Nitrogen evaporator, Organomation N-evap 115

Orbital shaker table, New Brunswick Scientific G-10

Pipet, 2-20 uL capacity, Rainin RL-20R

Pipet, 20-200 uL capacity, Rainin RL-200R

Pipet, 100-1000 uL capacity, Rainin R-1000R

Pipet, 200-2000 uL capacity, Rainin R-2000R

Pipet, 1-10 mL capacity, Rainin RL-10ML

Vortex mixer, Benchmark Scientific BV1000

GC-MS and LC-MS/MS systems

GC-MS system, Agilent 7890A

GC-MS Autosampler, Agilent 7693

GC-MS Detector, Agilent 5975C

Capillary GC column, Agilent DB-5MS (30m, 0.25mm ID, 1.00 μ m film)

GC-MS quantitation software, Agilent MSD ChemStation E.02.02.1431

Glassware and Materials

Vials, glass, 25 mL, catalogue number 06451173, Wheaton

Vials, glass, 4 mL, catalogue number C4015-1, Thermo Scientific

Vials, autosampler, 2 mL, catalogue number 03-375-3CC, National Scientific

Vial caps, PTFE-lined polypropylene, catalogue number B7815-24, National Scientific

Vial caps, PTFE-lined polypropylene, catalogue number B7815-13, Thermo Scientific

Vial cap, for autosampler vial, catalogue number 11-0030A-B, MicroLiter Analytical

Centrifuge tubes, polypropylene, 50 mL, catalogue number 05-539-9, Fisher Scientific

Disposable transfer pipette, glass, catalogue number 13-678-30, Fisher Scientific

Glass wool, 8 μ m fiber diameter, catalogue number 3950, Corning Life Sciences Glass

Separatory funnels, 125 mL

Graduated cylinders, glass, class A, 10-100 mL

Volumetric flasks, class A, 10-1000 mL

Beakers, glass, various sizes

Reagents

Hexane, Doe & Ingalls (catalogue number 9262-03), lot 0000082837

Toluene, Fisher Scientific HPLC grade (catalogue number T290SK-4), lot 144492

Acetonitrile, Fisher Optima LC/MS grade (catalogue number A955-4), lot 145298

Sulfuric acid, Doe & Ingalls (catalogue number 9673-00), lot 095813

Silica gel, Sigma-Aldrich (catalogue number 214477), lot BCBN1456V

Prepared Solutions

4:1 (v/v) Hexane:Toluene Solution

For each liter of solution prepared, combine 800 mL of hexane and 200 mL of toluene and mix well.

1:1 (v/v) Hexane:Toluene Solution

For each 4 liters of solution prepared, combine 2000 mL of hexane and 2000 mL of toluene and mix well.

6N Sulfuric Acid

For each 1 liter of solution prepared, slowly add 168 mL of concentrated sulfuric acid to a 1 L class A volumetric flask containing approximately 500 mL of deionized water. Slowly add deionized water to the base of the neck and mix. The solution will be very hot. Allow the solution to cool to room temperature and dilute to volume with additional deionized water. Mix well.

Stock Standard, 1000 $\mu\text{g}/\text{mL}$

Weigh 0.1 g of the neat nitrapyrin reference standard into a 100 mL volumetric flask and dilute to volume with 4:1 (v/v) hexane:toluene.

Working Standard, 50 $\mu\text{g}/\text{mL}$

Pipette 5 mL of the 1000 $\mu\text{g}/\text{mL}$ nitrapyrin stock solution into a 100 mL volumetric flask and dilute to volume with 4:1 (v/v) hexane:toluene.

Working Standard, 5 $\mu\text{g}/\text{mL}$

Pipette 10 mL of the 50 $\mu\text{g}/\text{mL}$ working solution into a 100 mL volumetric flask and dilute to volume with 4:1 (v/v) hexane:toluene.

Working Standard, 1 $\mu\text{g}/\text{mL}$

Pipette 10 mL of the 5 $\mu\text{g}/\text{mL}$ working solution into a 50 mL volumetric flask and dilute to volume with 4:1 (v/v) hexane:toluene.

Working Standard, 0.1 $\mu\text{g}/\text{mL}$

Pipette 1 mL of the 1 $\mu\text{g}/\text{mL}$ working solution into a 10 mL volumetric flask and dilute to volume with 4:1 (v/v) hexane:toluene.

Preparation of Spiking Solutions

Stock Standard, 100 µg/mL

Weigh 1 mg of the neat nitrapyrin reference standard into a 10 mL volumetric flask and dilute to volume with acetonitrile.

Fortification Standard, 5 µg/mL

Pipette 1.25 mL of the 100 µg/mL nitrapyrin stock solution into a 25 mL volumetric flask and dilute to volume with acetonitrile.

Fortification Standard, 0.5 µg/mL (500 ng/mL)

Pipette 0.125 mL of the 100 µg/mL nitrapyrin stock solution in acetonitrile into a 25 mL volumetric flask and dilute to volume with acetonitrile.

Fortification Standard, 0.05 µg/mL (50 ng/mL)

Pipette 2.5 mL of the 0.5 µg/mL nitrapyrin stock solution in acetonitrile into a 25 mL volumetric flask and dilute to volume with acetonitrile.

Preparation of Calibration Standards

Calibration standards are prepared by diluting the 1.0 µg/mL (in hexane/toluene) working solution with toluene as follows:

Concentration of Stock Solution (µg/mL)	Aliquot of Stock Solution (µL)	Final Soln. Volume (mL)	Calibration Soln. Final Conc. (ng/mL)	Equivalent Sample Conc. (µg/L) ^a
1.00	15	25.0	0.600	0.0180
1.00	37	25.0	1.48	0.0444
1.00	62	25.0	2.48	0.0744
1.00	187	25.0	7.48	0.224
1.00	375	25.0	15.0	0.450
1.00	625	25.0	25.0	0.750
1.00	875	25.0	35.0	1.05
1.00	1250	25.0	50.0	1.50

^aThe equivalent sample concentration is based on spiking a 100 mL matrix sample.

Analytical Procedure*Silica Gel Deactivation*

1. Fill a 500-mL Erlenmeyer flask with silica gel and store in an oven at 130°C ($\pm 10^\circ\text{C}$) for at least 24 hours before use.
2. Weigh out approximately 7.9 g of 6N sulfuric acid into a 250 mL Teflon bottle.
3. Add 150 g of silica to the bottle.
4. Immediately cap the bottle and place on a mechanical shaker to shake vigorously for approximately 2 hours.
5. Store the deactivated silica, tightly capped, in the Teflon bottle inside a desiccator with fresh desiccant.

Sample Extraction

1. Transfer a 100 mL aliquot of water sample into a 250 mL separatory funnel.
2. Fortify laboratory QC samples with the appropriate spiking solution as described in the table below:

To spike nitrapyrin to 100 mL of water:			
Description	Spiking Volumes (μL)	Spiking Solutions ($\mu\text{g}/\text{mL}$)	Fortification Level ($\mu\text{g}/\text{L}$)
Control	---	---	---
LOD	30.0	0.0500	0.0150
LOQ	100	0.0500	0.0500
20 \times LOQ	200	0.500	1.00

3. Add 10 mL of 1:1 (v/v) hexane:toluene to the separatory funnel.
4. Shake for approximately 15 seconds and allow the phases to separate. Methanol may be added drop-wise to break up any emulsions.
5. Drain the lower aqueous layer into a beaker.
6. Using a glass transfer pipet, transfer the top organic layer to a 25 mL glass vial. Take care not to transfer any of the aqueous phase that may be adhering to the sides of the separatory funnel.
7. Pour the aqueous phase in the beaker back into the separatory funnel and repeat steps 3-6 above, combining the organic layer with the previous extract in the 25 mL glass vial.
8. Evaporate the extract to less than 0.5 mL under a stream of nitrogen (N-Evap). Use a water bath temperature setting of 35°C.

Silica Gel Cleanup

1. Place a small plug of glass wool in a disposable Pasteur pipette and add approximately 20 mm of deactivated silica gel. Place another small plug of glass wool on top of the deactivated silica gel.

2. Rinse the column with 2 mL of hexane. Do not collect the hexane as it comes through the silica column.
3. Place a 4 mL glass vial beneath the silica column to collect the column eluate.
4. Add the concentrated extract from step 8 above to the column and allow it to drain to the top of the bed.
5. Rinse the 25 mL vial with 1 mL toluene and add it to the column. Allow the column to drain to the top of the bed.
6. Repeat the previous step using 1 mL of toluene.
7. Repeat the previous step using 0.5 mL of toluene.
8. Once the column stops dripping, force any remaining solvent from the silica column with a pipet bulb.
9. Adjust the eluate final volume in the vial to 3 mL (as compared to a reference vial) with toluene and vortex mix. Store sample solutions refrigerated when not needed in the laboratory.
10. Transfer a portion of the sample solution to an autosampler vial and cap. Analyze by GC/EI-MS.

Instrumental Conditions – GC/EI-MS Analysis

The following GC/EI-MS parameters are used to determine the concentration of nitrapyrin in sample extract solutions. The parameters may be modified to achieve adequate chromatographic resolution and/or detector sensitivity. The actual parameters used are documented with each GC/EI-MS sample sequence.

GC System: Agilent 7890A

Injector: Agilent Model 7693 with 150 sample capacity autosampler tray

Detector: Agilent Model 5975C mass spectrometer

Column: Agilent DB-5MS (30 m x 0.25 mm ID x 1.00 μ m film)

A splitless column inlet liner and the specified column should be installed into the split/splitless injection port according to manufacturer's recommended procedures.

Oven Method:

Column	130 °C for 1.0 minute
	130 °C to 210 °C at 10 °C/minute
	210 °C to 280 °C at 15 °C/minute
	280 °C for 2.0 minute

Carrier Gas:	Helium
Constant Flow:	1 mL/min

Injection: Splitless
Injector Temperature: 180 °C
Injection Volume: 2 µL
Purge flow: 40.0 mL/min

Detection Method: Positive-ion Electron-impact ionization mode
Source Temperature: 250 °C
Quad Temperature: 180 °C

Ions Monitored:

	Ion	Dwell
Quantitation:	<i>m/z</i> 194	50 msec
Confirmation (primary):	<i>m/z</i> 196	50 msec
Confirmation (secondary):	<i>m/z</i> 198	50 msec

Determination of 6-Chloropicolinic Acid (6-CPA) in Water by Liquid Chromatography with Tandem Mass Spectrometry Detection

Note: the method described herein was modified from EPL-BAS Method No. 205G907B Draft 11/25/14 to account for differences in differences in laboratory equipment or instrumentation. Changes to both methods provided by the sponsor are described in the "Changes to Method" section.

Method Summary:

6-CPA is extracted from water samples (20 mL) by C18 solid-phase extraction (SPE). 6-CPA is eluted from the SPE cartridge with acetonitrile. The acetonitrile eluate is evaporated to dryness. The dried residues are dissolved in 25/75 (v/v) methanol/DI water and submitted for liquid chromatography-tandem mass spectrometry analysis. The limit of quantitation for 6-CPA in water is 0.05 µg/L. The limit of detection is 0.015 µg/L.

Reference Substance and Reference Substance Solutions

Common Name	6-Chloropicolinic acid (6-CPA)
CAS Number	4684-94-0
Purity (%)	>95
Storage Conditions	Ambient

Volumes given may be adjusted as necessary as long as proportionality is maintained. Additional working and calibration standard concentrations may be used as needed.

Stock Standard Solution, 1,000 µg/mL

A stock solution of 6-CPA is prepared at a target concentration of 1,000 µg/mL. Weigh 0.010 g of the neat reference substance into a 10 mL volumetric flask and dilute to volume with methanol.

Working Standard Solution, 10 µg/mL

Pipet 0.1 mL of the 1,000 µg/mL stock solution into a 10 mL volumetric flask and bring to volume with methanol. Maintain refrigerated up to 365 days when not needed in the laboratory.

Working/Fortification Standard Solution, 1 µg/mL (1,000 ng/mL)

Pipet 1 mL of the 1,000 µg/mL stock solution into a 10 mL volumetric flask and bring to volume with methanol. Maintain refrigerated up to 365 days when not needed in the laboratory.

Fortification Standard Solution, 0.1 µg/mL (100 ng/mL)

Pipet 0.1 mL of the 10 µg/mL working solution into a 10 mL volumetric flask and bring to volume with methanol. Maintain refrigerated up to 365 days when not needed in the laboratory.

Fortification Standard Solution, 0.03 µg/mL (30 ng/mL)

Pipet 0.3 mL of the 1 µg/mL working solution into a 10 mL volumetric flask and bring to volume with methanol. Maintain refrigerated up to 365 days when not needed in the laboratory.

Calibration Standard Solution, 100 ng/mL

Pipet 0.1 mL of the 10 µg/mL working solution into a 10 mL volumetric flask and bring to volume with 25/75 (v/v) methanol/DI water. Maintain refrigerated up to 365 days when not needed in the laboratory.

Calibration Standard Solution, 50 ng/mL

Pipet 0.05 mL of the 10 µg/mL working solution into a 10 mL volumetric flask and bring to volume with 25/75 (v/v) methanol/DI water. Maintain refrigerated up to 365 days when not needed in the laboratory.

Calibration Standard Solution, 10 ng/mL

Pipet 0.1 mL of the 1 µg/mL working solution into a 10 mL volumetric flask and bring to volume with 25/75 (v/v) methanol/DI water. Maintain refrigerated up to 365 days when not needed in the laboratory.

Calibration Standard Solution, 5 ng/mL

Pipet 0.05 mL of the 1 µg/mL working solution into a 10 mL volumetric flask and bring to volume with 25/75 (v/v) methanol/DI water. Maintain refrigerated up to 365 days when not needed in the laboratory.

Calibration Standard Solution, 2 ng/mL

Pipet 0.2 mL of the 0.1 µg/mL working solution into a 10 mL volumetric flask and bring to volume with 25/75 (v/v) methanol/DI water. Maintain refrigerated up to 365 days when not needed in the laboratory.

Calibration Standard Solution, 1 ng/mL

Pipet 0.1 mL of the 0.1 µg/mL working solution into a 10 mL volumetric flask and bring to volume with 25/75 (v/v) methanol/DI water. Maintain refrigerated up to 365 days when not needed in the laboratory.

Calibration Standard Solution, 0.5 ng/mL

Pipet 0.05 mL of the 0.1 µg/mL working solution into a 10 mL volumetric flask and bring to volume with 25/75 (v/v) methanol/DI water. Maintain refrigerated up to 365 days when not needed in the laboratory

Reagents and Reagent Solutions

A. Reagents

Acetic acid, glacial, reagent grade

Acetonitrile, HPLC grade

Methanol, HPLC grade

Milli-Q filtered water

Hydrochloric acid, concentrated (36.5-38.0%), reagent grade (used to prepare 1N and 0.1N HCl below)

1.0N Hydrochloric acid, prepared internally

0.1N Hydrochloric acid, prepared internally

B. Reagent Solutions

1N Hydrochloric Acid

Fill a 1000 mL volumetric flask most of the way with Milli-Q water . Add 42.5 mL of concentrated hydrochloric acid and bring to volume with water. Mix well. The solution is stored ambient for up to 365 days.

0.1N Hydrochloric Acid

Fill a 1000 mL volumetric flask partially with Milli-Q water. Add 100 mL of 1N hydrochloric acid, bring to volume with water, and mix well. The solution is stored ambient for up to 365 days.

25/75 (v/v) Methanol/DI Water

Combine 250 mL of methanol with 750 mL of DI water. The solution is stored ambient for up to 365 days.

0.01% (v/v) Acetic Acid

Pipet 0.1 mL of glacial acetic acid into 1,000 mL of DI water and mix. The solution is stored ambient for up to 180 days.

Laboratory Equipment

Balance, analytical, capable of weighing to the nearest 0.01 mg

Scintillation vials, 40 mL, glass with PTFE-lined screw caps
Culture tubes, 15 mL, glass with PTFE-lined screw caps
SPE extraction manifold with vacuum pressure gauge pump
SPE Columns, Bakerbond Octadecyl (C18), 6 mL/1 g, cat. No. 7020-07
Pipettes, air-displacement, with disposable tips
Nitrogen evaporation apparatus, N-Evap, Organomation
Vortex mixer
Volumetric flasks and pipets
Autosampler vials with caps, low volume recovery

I. Methods

Sample Extraction

1. Dispense 20 mL of water sample into a 40 mL glass vial.
2. Fortify laboratory QC samples (20 mL) with the appropriate spiking solution.
 - i. *LOD Fortification (0.015 µg/L)*
Pipet 0.01 mL (10.0 µL) of the 30 ng/mL fortification standard solution into the water sample.
 - ii. *LOQ Fortification (0.05 µg/L)*
Pipet 0.01 mL (10.0 µL) of the 100 ng/mL fortification standard solution into the water sample.
 - iii. *20x LOQ Fortification (1 µg/L)*
Pipet 0.02 mL (20.0 µL) of the 1,000 ng/mL fortification standard solution into the water sample.
3. Add 2.0 mL of 1N hydrochloric acid. Cap and mix.

B. Solid-Phase Extraction

1. Place a C18 SPE column on an SPE vacuum manifold.
2. Rinse the column with 5.0 mL of acetonitrile and discard the rinse solvent.
3. Condition the SPE column with 5.0 mL of 0.1 N HCl. Do not allow the column bed to dry. Discard the conditioning solution.
4. Transfer the entire sample from Step IV.A.3 to the SPE column and draw through at approximately 1-2 mL/minute using a vacuum if gravity flow is inadequate to achieve the approximate flow rate of 1-2 mL/minute. Discard the column effluent.
5. Add 2.0 mL of 0.1 N HCl to the SPE column and draw through at approximately 1-2 mL/minute. Discard the column effluent.
6. Use vacuum (-10 inches Hg) to dry the SPE column for 40 minutes. **Do not use a vacuum pressure less than -10 inches Hg and do not exceed the 40 minute drying time or analyte loss may occur.**

7. Elute 6-CPA into a clean 15 mL glass culture tube with 5.0 mL of acetonitrile using a flow rate of approximately 1 mL/minute.
8. Once the column stops dripping, apply vacuum to collect the remaining eluent from the SPE column.
9. Evaporate the acetonitrile eluate just to dryness under a stream of nitrogen (N Evap). A water bath temperature setting of 35°C may be used. **Do not exceed the 35°C heating block temperature and do not over dry or analyte loss will occur. Evaporating the solvent to less than 100 µL then manually taking to dryness is recommended.**
10. Redissolve the dried sample residue by pipetting 0.5 mL of 25/75 (v/v) methanol/DI water into the culture tube.
11. Vortex mix for about 10 seconds at a speed adequate to thoroughly wash the sides of the culture tube.
12. Transfer the sample solution to a low volume recovery autosampler vial and cap. Store refrigerated if LC-MS/MS analysis cannot be initiated the same day.

LC-MS/MS Analysis

The following LC-MS/MS parameters are used to determine 6-CPA concentrations. The parameters may be modified to achieve adequate chromatographic resolution or detector sensitivity and selectivity. The actual parameters used are documented with each analysis sequence.

LC System: Agilent 1260 SL HPLC system (including a vacuum solvent degasser, binary HPLC pump, and column oven) and CTC Analytics HTC-Pal Autosampler,
MS Detector: Applied Biosystems MDS Sciex API 5500 linear ion trap MS/MS system with TurboIonspray (ESI) source

HPLC Conditions:

Column	Waters XSelect® HSS T3 3.5 µm, 3.0 x 100 mm
Column Temperature	40°C
Injection Volume	10 µL for batches 14-0606, 15-0003, 14-0618, 15-0002, 14-0625 and 15-0001 5 µL for batches 15-0016, 15-0027 and 15-0028
Mobile Phase	A: 0.01% (v/v) Acetic Acid B: Methanol

HPLC Gradient Program:

Time (minutes)	%A	%B	Curve	Flow (mL/minute)
0.00	95	5	1	0.40
0.10	95	5	1	0.40

4.00	30	70	6	0.40
4.50	30	70	1	0.40
4.51	5	95	1	0.40
6.50	5	95	1	0.40
6.51	95	5	1	0.40
8.50	95	5	1	0.40

Divert flow to waste from 0 to 0.5 minutes. Direct flow to the MS detector from 0.5 to 8.5 minutes. Divert flow to waste again at end of run.

Ion Source Conditions:

Source temperature	350 °C
Gas supply (GS 1):	40
Gas supply (GS 2)	60
Curtain gas:	30
CAD gas:	Medium
ESI Polarity:	Negative
Entrance potential:	-10 V
Ionspray Voltage	-4000 V

MS/MS Conditions for 6-CPA:

Quantification MS/MS transition	155.8/111.9
Collision energy (CE):	-13 V
Cell exit potential (CXP):	-13 V
Dwell time:	100 ms
Declustering potential (DP):	-50 V
Confirmation MS/MS transition	157.9/114.0
Collision energy (CE):	-13 V
Cell exit potential (CXP):	-13 V
Dwell time:	100 ms
Declustering potential (DP):	-50 V

Note that samples analyzed before 13 January 2015, utilized slightly different HPLC conditions and samples analyzed before 08 January 2015 has slightly different MS conditions. These conditions are described below:

HPLC Conditions:

Column	Water HSS T3 1.8 µm, 2.1 x 100 mm
Column Temperature	Same as above
Injection Volume	Same as above
Mobile Phase	A: Same as above B: Same as above

HPLC Gradient Program:

Time (minutes)	%A	%B	Curve	Flow (mL/minute)
-------------------	----	----	-------	---------------------

0.00	95	5	1	0.20
0.10	95	5	1	0.20
4.00	30	70	6	0.20
4.50	30	70	1	0.20
4.51	5	95	1	0.20
6.50	5	95	1	0.20
6.51	95	5	1	0.20
8.50	95	5	1	0.20

Flow diversion same as above.

Ion Source Conditions:

Source temperature	Same as above
Gas supply (GS 1):	Same as above
Gas supply (GS 2)	Same as above
Curtain gas:	Same as above
CAD gas:	Same as above
ESI Polarity:	Same as above
Entrance potential:	Same as above
Ionspray Voltage	Same as above

MS/MS Conditions for 6-CPA:

Quantification MS/MS transition	Same as above
Collision energy (CE):	Same as above
Cell exit potential (CXP):	Same as above
Dwell time:	Same as above
Declustering potential (DP):	Same as above
Confirmation MS/MS transition	Same as above
Collision energy (CE):	-25 V
Cell exit potential (CXP):	Same as above
Dwell time:	Same as above
Declustering potential (DP):	Same as above

The instrument data system is used to construct a linear regression external standard calibration curve of the calibrant peak areas (y-axis) vs. the standard concentrations in ng/mL (x-axis). A weighting factor of $1/x$ is used. The curve origin is not included in the linear regression equation.