Test Material:	Flufenacet				
MRID:	48897601				
Title:	Bethem, R.A., Peterson, R.G., and Leimkuehler, W.M. Determination of FOE 5043 and Four Metabolites in Soil By High Performance Liquid Chromatography Electrospray Tandem Mass Spectrometry (LC-ESI/MS/MS).				
MRID:	48897602				
Title:		wski, J. Independent Laboratory Validation of MFOE2 for the Analysis of Soil for FOE 5043 and dues			
EPA PC Code:	121903				
OCSPP Guideline:	850.6100				
For CDM Smith					
Primary Reviewer: R	ichard Lester	Signature: Richard R Sester Date: 5/14/14			
Secondary Reviewer:	Dan Hunt	Signature: $\int \mathcal{D} \mathcal{D} \mathcal{D} \mathcal{D} \mathcal{D} \mathcal{D} \mathcal{D} \mathcal{D}$			
QC/QA Manager: Joa	an Gaidos	Signature: Jour St. Date: 5/14/14			

#### Analytical method for flufenacet in soil

Reports:	<ul> <li>ECM: MRID 48897601. Bethem, R.A., Peterson, R.G., and Leimkuehler,</li> <li>W.M. 1995. Determination of FOE 5043 and Four Metabolites in Soil By</li> <li>High Performance Liquid Chromatography Electrospray Tandem Mass</li> <li>Spectrometry (LC-ESI/MS/MS). Submitting Laboratory: Bayer</li> <li>Corporation, Kansas City, Missouri, Report No. 106964. Performing</li> <li>Laboratory: Alta Analytical, El Dorado Hills, California, Study No.</li> <li>AMFOE2. September 6, 1995. 20 pages.</li> <li>ILV: Bruns, G. and Hoshowski, J. 1998. Independent Laboratory</li> <li>Validation of Analytical Method AMFOE2 for the Analysis of Soil for FOE</li> <li>5043 and Four Metabolite Residues. Submitting Laboratory: Bayer</li> <li>Corporation, Kansas City, Missouri, Report No. 108145, Study No.</li> <li>F3112106. Performing Laboratory: Enviro-Test Laboratories, Edmonton,</li> <li>Alberta, Canada, Document No. 98BAY04.REP. March 5, 1998. 90 pages.</li> </ul>					
<b>Document No.:</b>	MRIDs 48897601 and 48897602					
Guideline:	850.6100					
Statements:	<ul> <li>ECM: The study was not conducted in compliance with FIFRA GLP standards (p. 3). A signed and dated Data Confidentiality statement was provided (p. 2). Quality Assurance and Certificate of Authenticity statements were not provided from the performing laboratory.</li> <li>ILV: The study was conducted in accordance with FIFRA GLP standards (p. 3). Signed and dated Data Confidentiality, GLP Compliance, and Quality Assurance statements were provided (pp. 2-4, 8-9). A Certificate of Authenticity statement was not provided from the performing laboratory.</li> </ul>					
Classification:	This analytical method is classified as <b>Supplemental</b> . The ECM did not provide performance data at $10 \times \text{LOQ}$ . The soil used in the ECM and ILV was not characterized. Representative chromatograms were not provided for reagent blanks, matrix blanks or spiked samples in the ECM.					
PC Code:	121903					
<b>Reviewer:</b>	Gabriel RothmanSignature:Environmental Scientist, USEPADate: 9/24/2015					

### All page citations refer to MRID 48897601 unless otherwise noted.

#### **Executive Summary**

This analytical method, Alta Analytical Study No. AMFOE2, is designed for the quantitative determination of FOE 5043 (flufenacet) and metabolites FOE 5043 sulfonic acid, FOE 5043 alcohol, FOE 5043 oxalate, and FOE 5043 thiadone in soil using LC-ESI/MS/MS (see Table 1). The method is quantitative for the analytes at the stated LOQ of 10  $\mu$ g/kg. Lowest toxicological levels of concern in soil were not reported. No major issues were discovered by the independent laboratory.

	MR	ID						Limit of
Analyte(s) by Pesticide	Environmental Chemistry Method	Independent Laboratory Validation	EPA Review	Matrix	Method Date	Registrant	Analysis	Quantitation (LOQ)
FOE 5043 (flufenacet), FOE 5043 sulfonic acid, FOE 5043 alcohol, FOE 5043 oxalate, FOE 5043 thiadone	48897601	48897602		Soil	09/06/95	Bayer Corporation	LC- ESI//MS	10 µg/kg

## **Table 1. Analytical Method Summary**

# I. Principle of the Method

The method was developed for the analysis of FOE 5043 and four metabolites (FOE 5043 sulfonic acid, FOE 5043 alcohol, FOE 5043 oxalate, and FOE 5043 thiadone) in soil (p. 6). Soil samples (10 g) are extracted with 20 mL of 0.1N HC1:ACN (1:1) by shaking on an orbital shaker for 1 hour (*ca.* 220 rpm; pp. 10-11). After centrifuging for 10 minutes, a 10 mL aliquot of the extract is transferred to a conical tube. After the addition of 100  $\mu$ L of an internal standard solution and *ca.* 1 mL of methanol, the solution is concentrated to *ca.* 4.9 mL using nitrogen and a water bath at 25-30°C. The concentrate is brought up to 5 mL with 0.1% formic acid and a portion of the extract is syringe filtered ( $\leq 0.45 \mu$ m).

Extracts were analyzed by LC-ESI/MS/MS (5g/5mL) using an Inertsil ODS-2 column (50 x 3 mm, 5  $\mu$ ; Tables I-II, pp. 15-16).The analytes are determined by reversed phase liquid chromatography, using 0.1% formic acid and acetonitrile with 0.1% formic acid. The standards and sample extracts are injected onto a base deactivated reversed phase column connected to a triple quadrupole mass spectrometer. An electrospray (ESI) atmospheric pressure ionization (API) inlet is used. The sulfonic acid, thiadone, and oxalate metabolites are analyzed by negative ionization MS/MS and the alcohol and parent are analyzed by positive ionization MS/MS using a second injection.

Quantitation is performed using the area response factors of the native compounds relative to their stable isotope internal standards. A calibration check standard (CCS) is analyzed at the onset and completion of every analytical sample set. The response factor from each calibration check standard is compared to the average response factor from a triplicate 4 point calibration curve.

# **II. Recovery Findings**

Mean recoveries and relative standard deviations (RSD) were within guideline requirements (mean 70-120%; RSD  $\leq$ 20%) for all analytes in soil. No validation was performed at the 10×LOQ fortification level in the ECM.

Analyte	Fortification Level (µg/kg)	Number of Tests	Recovery Range (%) <sup>1</sup>	Mean Recovery (%) <sup>1</sup>	Standard Deviation (%) <sup>2</sup>	Relative Standard Deviation (%)
FOE 5043	10 (LOQ)	7	90-109	101	5.8	5.7
FOE 5043 sulfonic acid	10 (LOQ)	7	89-110	99	8.1	8.2
FOE 5043 oxalate	10 (LOQ)	7	83-101	91	5.5	6.0
FOE 5043 thiadone	10 (LOQ)	7	98-109	104	4.1	3.9
FOE 5043 alcohol	10 (LOQ)	7	97-112	104	5.4	5.2

Data were obtained from Tables III-IV, pp. 17-18 of the ECM study report.

<sup>1</sup> Recoveries are based on internal standard quantitation as specified in the method (p. 12). The ECM separately reports recoveries based on external standard quantitation.

<sup>2</sup> Standard deviations of measured concentrations are reported in the ECM. See reviewer Excel file for standard deviations of percent recoveries.

Analyte	Fortification Level (µg/kg)	Number	•	Mean Recovery (%)	Standard Deviation (%) <sup>1</sup>	Relative Standard Deviation (%)
FOE 5043	11.5 (LOQ)	5	90-97	93	2.3	2.5
FUE 3043	115	5	83-91	87	3.9	4.4
FOE 5043 sulfonic acid	10.5 (LOQ)	5	96-105	99	3.7	3.7
	105	5	81-86	84	2.0	2.4
FOE 5043 oxalate	10.2 (LOQ)	5	78-90	85	4.7	5.6
	102	5	74-79	77	2.6	3.4
FOE 5043 thiadone	10.9 (LOQ)	5	94-103	96	3.8	4.0
	109	5	80-90	87	3.9	4.5
FOE 5043 alcohol	10.4 (LOQ)	5	82-95	90	5.0	5.6
	104	5	82-86	84	2.2	2.6

 Table 3. Independent Validation Method Recoveries for Analytes in Soil

Data were obtained from MRID 48897602, pp. 27.

<sup>1</sup> Standard deviations of measured concentrations are reported in the ILV. See reviewer Excel file for standard deviations of percent recoveries.

## **III. Method Characteristics**

The LOD was calculated in the ECM as the students' t-value  $\times$  the standard deviation, citing 40 CFR Part 136, Appendix B (p. 17). The LOQ was calculated as 10  $\times$  the standard deviation, using scientifically accepted calculation procedures. Reproducibility is not satisfactory because no LODs or LOQs were calculated in the ILV. Observed standard deviations in the ILV, however, were similar to those in the ECM.

	FOE 5043	FOE 5043 sulfonic acid	FOE 5043 oxalate	FOE 5043 thiadone	FOE 5043 alcohol
Limit of Quantitation (LOQ)	5.8 µg/kg	8.1 µg/kg	5.5 µg/kg	4.1 µg/kg	5.4 µg/kg
Limit of Detection (LOD)	1.8 µg/kg	2.5 µg/kg	1.7 µg/kg	1.3 µg/kg	1.7 µg/kg
Linearity (calibration curve $r^2$ and concentration range)	$r^2 = 1.00$ 5-200 µg/L	$r^2 = 0.998$ 5-200 µg/L	$r^2 = 1.00$ 5-200 µg/L	$r^2 = 0.999$ 5-200 µg/L	$r^2 = 1.00$ 5-200 µg/L
Repeatable	Yes	Yes	Yes	Yes	Yes
Reproducible	No	No	No	No	No
Specific	Yes	Yes	Yes	Yes	Yes

### **Table 4. Method Characteristics**

The LOQ and LOD are those reported in the ECM (Tables III and IV, pp. 17-18). Linearity data are those reported in the ILV (Trial II, pp. 19, 26).

### **IV. Method Deficiencies and Reviewer's Comments**

- 1. The ECM did not provide performance data at  $10 \times LOQ$ .
- 2. LODs and LOQs were not calculated in the ILV.
- 3. The soil used in both the ECM and the ILV was not characterized. Only the source of the soil used in the ILV was reported (Wisconsin; MRID 48897602, p. 14).
- 4. The ECM only provided representative chromatograms for standards and calibration curves were not provided (Figures 1-2, pp. 19-20). Representative chromatograms should be provided for reagent blanks, matrix blanks, standard curves, and spiked samples at all spiking levels for all analytes.
- 5. The lowest toxicological level of concern in soil was not reported. The established LOQ must be less than toxicological levels of concern.
- 6. The ILV was validated with the second trial. Low response of the thiadone internal standard from the initial QC check calibration curve used for Trial 1 resulted in low recoveries for the 100 ng/g fortification (mean of 66%, with a RSD of 29%; MRID 48897602; pp. 13, 25). Following communication with the Study Monitor, a new standard for thiadone was provided and new standard solutions were prepared for use in the second trial.
- 7. The ECM provided calculation steps (pp. 13-14), but did not provide example calculations showing how the raw data were converted to a final concentration.

- 8. A communication log between the submitting laboratory (Bayer) and the ILV (Enviro-Test Laboratories) was provided in Appendix III of MRID 48897602 (p. 61). A brief summary of the discussion of the phone conversations following the Trial 1 results was provided on page 13 of the ILV.
- 9. It was noted in the ILV that the LC/MS/MS system sensitivity was significantly greater than what appeared in the analytical method and in order to obtain similar sensitivities to those in the method, the injection volume was changed from 50 to 20 μL for the negative ion analysis (MRID 48897602, p. 23). Also a conversion in the mass spectrometry and interface parameters was required due to the use of a Turbo Ion Spray API interfaced with a PE Sciex LC/MS/MS system versus a heated capillary API interfaced with a Finnigan LC/MS/MS which was employed in the original method.

## V. References

Attachment 1: Chemical Na Code Name/ Synonym Flufenacet	ames and Structures Chemical Name	Chemical Structure
(FOE 5043)	IUPAC: 4'-Fluoro-N-isopropyl-2-[5-(trifluoromethyl)-1,3,4-thiadiazol- 2-yloxy]acetanilide	F
	CAS: N-(4-fluorophenyl)-N-(1-methylethyl)-2-[[5-(trifluoromethyl)- 1,3,4-thiadiazol-2-yl]oxy]acetamide	N S O
	CAS #: 142459-58-3	О <u>С</u> H <sub>2</sub> С N <u>С</u> С С H <sub>3</sub> С С С H <sub>3</sub>
	Formula: C <sub>14</sub> H <sub>13</sub> F <sub>4</sub> N <sub>3</sub> O <sub>2</sub> S MW: 363.3 g/mol SMILES: CC(C)N(c1ccc(cc1)F)C(=O)COc2nnc(s2)C(F)(F)F	
FOE 5043 Sulfonic Acid, Sodium Salt, Monohydrate	IUPAC: 4-Fluoro-N-methylethylanilinesulfoacetamide,sodiumsalt, monohydrate	
	Formula: C11H15FNO5SNa MW: 315.29 g/mol SMILES: [H]O([H])[Na+][O-]S(=O)(=O)CC(=O)N(c1ccc(cc1)F)C(C)C	F F
FOE 5043 Alcohol	IUPAC: N-(4-fluorophenyl)-2-hydroxy-N-(1-methylethyl)acetamide	
	Formula: C <sub>11</sub> H <sub>14</sub> FNO <sub>2</sub> MW: 211.2 g/mol SMILES: CC(C)N(c1ccc(cc1)F)C(=O)CO	Р ОН
	IUPAC: N-(4-fluorophenyl)-N-methylethyl-aminooxoacetic acid	
FOE 5043 Oxalate	Formula: C <sub>11</sub> H <sub>12</sub> FNO <sub>3</sub> MW: 225.2 g/mol SMILES: CC(C)N(c1ccc(cc1)F)C(=O)C(=O)O	Р ОН
	CAS: 5-Trifluoromethyl-1,3,4-thiadiazol-2(3H)ol	NN
FOE 5043 Thiadone	CAS #: 84352-75-0	HO
	Formula: C <sub>3</sub> HF <sub>3</sub> N <sub>2</sub> OS MW: 170.1 g/mol SMILES: c1(nnc(s1)O)C(F)(F)F	S F