




ACR #: 91.3
EFFECTIVE: April 29, 1991
SUPERCEDES: NEW

DETERMINATION OF RESIDUES OF DE-498 IN SOIL
BY
CAPILLARY GAS CHROMATOGRAPHY/MASS SPECTROMETRY

E. L. Olberding, D. R. Foster, B. J. Harnick, and J. L. Balcer
Formulation and Environmental Chemistry
DowElanco
Midland, Michigan 48641-1706

1. Scope

This method is applicable for the quantitative determination of DE-498 (N-(2,6-difluorophenyl)-5-methyl-1,2,4-triazolo-[1,5a]-pyrimidine-2-sulfonamide) in soil at a validated lower level of quantitation of 2.5 ppb.



N-(2,6-difluorophenyl)-5-methyl-1,2,4-triazolo-[1,5a]-
pyrimidine-2-sulfonamide (DE-498)

2. Principle

DE-498 residues are extracted from soil using a 90% acetone/10% 0.1 N hydrochloric acid solution. Following evaporation of the acetone, the sample is diluted with 0.005 N hydrochloric acid and purified using a Cis solid-phase extraction (S-P-E). The eluent from the S-P-E is evaporated to dryness, and the residue reconstituted with acetonitrile. The sample is then derivatized with methyl iodide to form the N-methyl derivative. The derivatized sample solution is evaporated to dryness, reconstituted with toluene containing N-d₃-methyl DE-498 as an internal standard, and analyzed by capillary gas chromatography/mass spectrometry.

3. Safety Precautions

- a. Each analyst should be acquainted with the potential hazards of the reagents, products, and solvents used in this method before commencing laboratory work. SOURCES OF INFORMATION INCLUDE: MATERIAL SAFETY DATA SHEETS, LITERATURE, AND OTHER RELATED DATA. Safety information on non-DowElanco products should be requested from the supplier. Disposal of reagents, reactants, and solvents must be in compliance with local, state, and federal laws and regulations.
- b. Acetone, acetonitrile, methanol, methyl-t-butyl ether, toluene, and triethylamine are flammable and should be used in well-ventilated areas away from ignition sources.

4. Equipment

- a. Gas chromatograph, Model 5890A, Hewlett-Packard, Avondale, PA 19311.
- b. Automatic sampler, Model 7673A, Hewlett-Packard, Avondale, PA 19311.
- c. Mass selective detector, Model 5971A, Hewlett-Packard, Palo Alto, CA 94304.
- d. Mass spectrometer data system, Model 59970, Hewlett-Packard, Palo Alto, CA 94304.
- e. Balance, analytical, Model AE200, Mettler Instrument Corp., Hightstown, NJ 08520.
- f. Balance, pan, Model BB2440, Mettler Instrument Corp.
- g. Centrifuge, with head to accommodate 10-dram vials, Model Centra-8, International Equipment Company, Needham Heights, MA 02194.
- h. Desiccator, 250 mm I.D. with Drierite adsorbent, Catalog Number 08-595E, Fisher Scientific, Pittsburg, PA 15219.
- i. Evaporator, N-Evap, Model 111, Organomation Associates, Inc., South Berlin, MA 01549.
Set at a water bath temperature of 40°C and a nitrogen flow rate of 200 mL/min.
- j. Oven, Model OV-490A-2, Blue M Electric Company, Blue Island, IL 60406.
- k. Shaker, variable-speed reciprocating with box carrier, Model 6000, Eberbach Corp., Ann Arbor, MI 48106.
- l. Ultrasonic bath, Model 1200, Branson Cleaning Equipment Company, Shelton, CT 06484.

- m. Vacuum manifold box, Model spe-21, J. T. Baker Chemical Company, Phillipsburg, NJ 08865.
- n. Vial crimper, Part Number 8710-0979, Hewlett-Packard, Avondale, PA 19311.
- o. Vortex mixer, Model G-560, Scientific Industries, Inc., Bohemia, NY 11716.

5. Glassware and Materials

- a. Column, capillary gas chromatography, Durabond-17 liquid phase, 10 m x 0.18 mm i.d., 0.3 μ m film thickness, Catalog Number 121-1713, J&W Scientific, Folsom, CA 95630.
- b. Column inlet liner, deactivated, Catalog Number 5181-3315, Hewlett-Packard, Avondale, PA 19311.
- c. Column, Cis S-P-E, Catalog Number 7020-07, J. T. Baker Chemical Company. (Note 16.b)
- d. Cylinder, graduated, 2000 mL, Catalog Number 131-9058, National Scientific Company, Lawrenceville, Georgia 30245.
- e. Dish, weighing, Catalog Number 08-732, Fisher Scientific.
- f. Gas, helium, 99.995% purity, Scott Specialty Gases, Troy, MI 48063.
- g. Gas, nitrogen, technical grade, Scott Specialty Gases.
- h. Moisture trap, Catalog Number 7971, Chrompack, Inc., Raritan, NJ 08869. (Note 16.c.)
- i. Charcoal scrubber, Catalog Number 7972, Chrompack, Inc. (Note 16.c.)
- j. Oxygen trap, Catalog Number 7970, Chrompack, Inc. (Note 16.c.)
- k. Syringes, 10, 50, and 500 μ L, Model 700 Series, Hamilton Company, Reno, NV 89520.
- l. Vials, 2 dram, with poly(tetrafluoroethylene)-lined screw caps, Catalog Number B7800-3, National Scientific Company.
- m. Vials, 10 dram, with poly(tetrafluoroethylene)-lined screw caps, Catalog Number B7800-6, National Scientific Company.
- n. Vials, autosampler, 2 mL, Catalog Number C4011-2, National Scientific Company.
- o. Vial seals, Catalog Number C4011-1A, National Scientific Company.

6. Reagents and Chemicals

- a. Acetone, acetonitrile, methanol, methyl-t-butyl ether, and toluene (Optima Grade), Fisher Scientific.
- b. Hydrochloric acid, 0.1 N, reagent grade, certified concentration, Fisher Scientific.
- c. Sodium chloride, ACS reagent grade, Fisher Scientific.

d. Hydrochloric acid, 0.005 N.

Prepare by diluting 50 mL of 0.1 N hydrochloric acid to volume in a 1000-mL volumetric flask with distilled/deionized water.

e. Sodium chloride, 5% (w/v).

Prepare by dissolving 50 grams of sodium chloride in distilled/deionized water in a 1000-mL volumetric flask. Adjust to volume with distilled/deionized water.

f. 90% acetone/10% 0.1 N hydrochloric acid solution.

Prepare by pouring 200 mL of 0.1 N hydrochloric acid into a 2000-mL graduated cylinder. Add 1500 mL of acetone, swirl the cylinder, and allow to equilibrate to room temperature. Adjust to volume with acetone.

g. Methyl iodide, minimum 99.5% purity, Catalog Number 28,956-6, Aldrich Chemical Company, Milwaukee, WI 53233.

h. Methyl iodide, stable-isotope labeled, $^{12}\text{CD}_3\text{I}$, Catalog Number 29,675-9, Aldrich Chemical Company.

i. Triethylamine, minimum 99% purity, Catalog Number 13,206-3, Aldrich Chemical Company.

j. Water, distilled/deionized, Corning MEGA-PURE Still, Model MP-12A, Corning Glass Works, Science Products Division, Corning, NY 14831.

k. Standard

N-(2,6-difluorophenyl)-5-methyl-1,2,4-triazolo-[1,5a]-pyrimidine-2-sulfonamide (DE-498), analytical standard.

Obtain from Sample Coordinator, DowElanco, P.O. Box 1706, Midland, Michigan 48641-1706.

7. Preparation of Standards**a. Preparation of Calibration Standards/Spiking Solutions**

- (1) Dissolve 0.1000 gram of DE-498 analytical standard in acetone in a 100-mL volumetric flask. Dilute to volume to obtain a 1000 $\mu\text{g/mL}$ stock solution.
- (2) Dilute 5 mL of the above 1000 $\mu\text{g/mL}$ solution to 1000 mL with acetone in a 1000-mL volumetric flask to obtain a 5.0 $\mu\text{g/mL}$ (5.0 $\text{ng}/\mu\text{L}$) initial solution.
- (3) Solutions for spiking soil samples are prepared by diluting the initial solution from Section 7.a.(2) above with acetone as follows:

<u>Allquot of Initial Soln.</u>	<u>Final Soln. Volume</u>	<u>Spiking Soln. Final Conc.</u>	<u>Equivalent Sample Conc.</u>
mL	mL	ng/mL	ppb
5.0	2000	12.5	2.5
10.0	2000	25.0	5.0
20.0	2000	50.0	10.0
50.0	2000	125.0	25.0
10.0	200	250.0	50.0
20.0	200	500.0	100.0
50.0	200	1250.0	250.0

- (4) Solutions for calibration standards are prepared by pipeting 1.0 mL of the DE-498 standards in Section 7.a.(3) above into 2-dram vials and derivatizing according to the procedure described in Section 9, steps k through w.

b. Preparation of Internal Standard Solution

- (1) Pipet 2.0 mL of the 1000 $\mu\text{g/mL}$ DE-498 stock solution from Section 7.a.(1) into a 2-dram vial.
- (2) Evaporate the solution to dryness using an N-Evap evaporator.
- (3) Add 1.0 mL of acetonitrile, cap the vial, and sonicate for 5-10 seconds.
- (4) Add 50 μL of triethylamine and 50 μL of stable-isotope labeled methyl iodide (Section 6.h), cap the vial, and sonicate for 5-10 seconds.
- (5) Allow the mixture to react with the methyl iodide for 30 minutes at room temperature.
- (6) Evaporate the solution to dryness using an N-Evap evaporator.
- (7) Add 1.0 mL of a 5% sodium chloride solution, cap the vial, and sonicate for 5-10 seconds.

- (8) Add 5.0 mL of methyl-t-butyl ether, cap the vial, and vortex the sample for 5-10 seconds.
- (9) Centrifuge the vial for 5 minutes at 2500 rpm.
- (10) Carefully transfer the methyl-t-butyl ether layer to a clean 10-dram vial.
- (11) Repeat Steps 8-9 three additional times, combining the methyl-t-butyl ether layers in the 10-dram vial.
- (12) Evaporate the solution to dryness using an N-Evap evaporator.
- (13) Add 20 mL of acetone, cap the vial, and sonicate for 5-10 seconds.
- (14) Transfer the acetone to a 200-mL volumetric flask.
- (15) Rinse the 10-dram vial again with 20 mL of acetone, and transfer the acetone to the 200-mL volumetric flask.
- (16) Dilute the solution to volume with acetone. This solution contains 10.0 µg/mL N-d3-methyl DE-498.
- (17) Dilute 10.0 mL of the above 10.0 µg/mL solution to 1000 mL with toluene in a 1000-mL volumetric flask to obtain a 0.100 µg/mL (0.100 ng/µL) solution.

8. Gas Chromatography/Mass Spectrometry

a. Column

Install the splitless liner (5.b) and the capillary column (5.a) on the split/splitless injection port of the GC/MS following the manufacturer's recommended procedure.

b. Typical operating conditions for the determination of DE-498 by capillary MS:

Instrumentation: Hewlett-Packard Model 5890A Gas Chromatograph / Model 5971A Mass Selective Detector

Column: J&W Scientific fused silica capillary Durabond-17 liquid phase
10 m x 0.18 mm i.d.
0.30 µm film thickness

Temperatures:

Column	120°C for 1.1 minutes
	120°C to 325°C at 20°C/minute
	325°C for 5.65 minutes
Injector	320°C
Interface	310°C

EFFECTIVE: April 29, 1991

ACR #: 91.3

Carrier Gas: helium
Head Pressure: 100 kPa
Linear Velocity: 25 cm/sec
Injection Mode: splitless
Purge Delay: 1.0 minutes
Splitter Flow: 50 mL/min
Septum Purge: 1.0 mL/min
Injection Volume: 2 µL
Ions Monitored: N-methyl DE-498
m/z 134 (base peak ion)
(M⁺-205; see Figure 1)
m/z 142 (M⁺-197; see Figure 1)
N-d3-methyl DE-498 (internal standard)
m/z 145 (M⁺-197; see Figure 2)

Electron Multiplier: 1800 volts

- c. A typical calibration curve is shown in Figure 3.
- d. Typical chromatograms of a standard, control sample, and a 2.5 ppb recovery sample are shown in Figures 4-6, respectively.

9. Recovery of DE-498 from Soil

- a. Weigh 5.0-gram portions of control soil into a series of 10-dram vials.
- b. For preparing fortified samples, use part of the samples as controls and fortify the remaining samples by adding 1.0-mL aliquots of the appropriate spiking solutions (Section 7.a.(3)) in acetone to obtain concentrations ranging from 2.5 to 250 ppb.
- c. Add 25 mL of a 90% acetone/10% 0.1 N hydrochloric acid extracting solution.
- d. Cap the vial and sonicate the sample for 30-45 seconds.
- e. Shake the sample for a minimum of 2 hours on a reciprocating shaker at approximately 180 excursions/minute.
- f. Centrifuge the sample container for 10 minutes at 2500 rpm.
- g. Transfer the acetone/hydrochloric acid solution to a clean 10-dram vial.
- h. Evaporate the acetone using an N-Evap evaporator.

Page 7 of 19

62A ASOILS AM 10

22 TOTAL PAGES

1. Add 15.0 mL of 0.005 N hydrochloric acid, cap the vial, and sonicate the sample for 10-15 seconds.
- J. The sample is then purified using the following S-P-E procedure (Note 16.b):
 - (1) Place a Cis S-P-E column on the vacuum manifold box.
 - (2) Rinse the S-P-E column with 5 mL of methanol.
 - (3) Condition the S-P-E column with 5 mL of 0.005 N hydrochloric acid. (Do not allow the column bed to dry.)
 - (4) Transfer the sample solution from Step 9.1 to the S-P-E column and, with the aid of vacuum, slowly pull the sample through the column. Without allowing the column bed to dry, wash the sample vial with a 10 mL aliquot of 0.005 N hydrochloric acid and transfer the wash to the S-P-E column.
 - (5) Thoroughly dry the S-P-E column by drawing air through it for approximately 45 minutes.
 - (6) Remove the S-P-E column from the vacuum manifold box and elute the DE-498 by passing 2.5 mL of methanol through the S-P-E column, collecting the eluent in a 2-dram vial. (Note 16.d)
- k. Evaporate the solution to dryness using an N-Evap evaporator.
 1. Add 500 µL of acetonitrile, cap the vial, and sonicate for 5-10 seconds.
 - m. Add 10 µL of triethylamine and 10 µL of methyl iodide, cap the vial, and sonicate for 5-10 seconds.
 - n. Allow the sample to react with the methyl iodide for 30 minutes at room temperature.
 - o. Evaporate the solution to dryness using an N-Evap evaporator.
 - p. Add 1.0 mL of a 5X sodium chloride solution, cap the vial, and sonicate for 5-10 seconds.
 - q. Add 5.0 mL of methyl-t-butyl ether, cap the vial, and vortex the sample for 5-10 seconds.
 - r. Centrifuge the vial for 5 minutes at 2500 rpm.
 - s. Carefully transfer the methyl-t-butyl ether layer to a 2-dram vial.
 - t. Evaporate the solution to dryness using an N-Evap evaporator.
 - u. Add 1.0 mL of toluene containing the N-d3-methyl DE-498 internal standard, cap the vial, and sonicate for 5-10 seconds.
 - v. Centrifuge the vial for 5 minutes at 2600 rpm.

EFFECTIVE: April 29, 1991

ACR #: 91.3

- w. Transfer the solution to a 2-ml autosampler vial. Seal the vial with a cap and crimper.
- x. Analyze the sample by capillary gas chromatography/mass spectrometry.

10. Determination of Percent Recovery of DE-498

- a. Inject the calibration standards described in Section 7.a.(4) and determine the peak areas at m/z 134 and m/z 142 for methylated DE-498 and at m/z 145 for d_3 -methylated DE-498.

For each standard calculate the DE-498 confirmation ratio. The average confirmation ratio for all of the calibration standards will be used to confirm the presence of DE-498 in the soil samples.

For example, using the data from Figure 4:

$$\text{Confirmation Ratio} = \frac{\text{peak area at } m/z \text{ 142}}{\text{peak area at } m/z \text{ 134}}$$

$$\text{Confirmation Ratio} = \frac{42425}{76598}$$

$$\text{Confirmation Ratio} = 0.55387$$

Positive confirmation of the presence of DE-498 is indicated when the confirmation ratio for the samples is in the range of $\pm 10\%$ of the average found for the standards.

- b. Prepare a standard curve by plotting the equivalent DE-498 concentration on the abscissa (x-axis) and the m/z 134/145 peak area ratio on the ordinate (y-axis) as shown in Figure 3. Using regression analysis, determine the equation for the curve with respect to the abscissa.

For example, using power regression with the data from Figure 3

$$Y = \text{constant} \cdot X^{\text{(exponent)}}$$

$$X = \left(\frac{Y}{\text{constant}} \right)^{1/\text{exponent}}$$

$$\text{DE-498 Conc. (ppb)} = \left(\frac{m/z \text{ 134/145 peak area ratio}}{\text{constant}} \right)^{1/\text{exponent}}$$

$$\text{DE-498 Conc. (ppb)} = \left(\frac{m/z \text{ 134/145 peak area ratio}}{0.089690} \right)^{1/0.963883}$$

- c. Determine the net concentration in each recovery sample by first subtracting the average DE-498 peak area ratio in the control samples from that of the recovery sample. Substitute the peak area ratio obtained into the above equation and solve for the concentration.

$$\text{DE-498 Conc. (ppb)} = \left(\frac{\text{net } m/z \text{ 134/145 peak area ratio}}{0.089690} \right) 1/0.963885$$

$$\text{DE-498 Conc. (ppb)} = \left(\frac{0.210707 - 0.000000}{0.089690} \right) 1/0.963885$$

$$\text{DE-498 Conc.} = 2.43 \text{ ppb}$$

- d. Determine the percent recovery by dividing the net concentration of each recovery sample by the theoretical concentration added.

$$\text{Recovery} = \frac{\text{Concentration Found}}{\text{Concentration Added}} \times 100\%$$

$$\text{Recovery} = \frac{2.43 \text{ ppb}}{2.50 \text{ ppb}} \times 100\%$$

$$\text{Recovery} = 97\%$$

11. Determination of DE-498 in Soil Samples

- Prepare control, recovery, and treated samples as described in Section 9.
- Prepare a standard curve and determine the DE-498 concentration in the recovery samples as described in Section 10.
- Determine the concentration in each treated sample by substituting the DE-498 m/z 134/145 peak area ratio obtained into the equation for the standard curve and solve for the concentration.

For example, using the data from Figure 6:

$$\text{DE-498 Conc. (ppb)} = \left(\frac{m/z \text{ 134/145 peak area ratio}}{\text{constant}} \right) 1/\text{exponent}$$

$$\text{DE-498 Conc. (ppb)} = \left(\frac{0.210707}{0.089690} \right) 1/0.963885$$

$$\text{DE-498 Conc.} = 2.43 \text{ ppb}$$

13