

INTRODUCTION

A method for determining Iprodione, RP-30228 and RP-32490 residues in field sediment is presented. The methodology for extraction, purification and detection of Iprodione and its degradates was developed for freshwater and seawater sediments.

PRINCIPLE AND APPLICATION

Iprodione and its degradates are extracted from field sediment by acidifying with hydrochloric acid and partitioning on a shaker table with 50% dichloromethane : 25% acetone : 25% ethyl acetate. The extract is filtered through a glass fiber membrane, dried through sodium sulfate, rotary evaporated to approximately 0.5 mL, evaporated to dryness under nitrogen and the residues dissolved in acetonitrile. The acetonitrile solution is partitioned twice with an equal volume of hexane, the acetonitrile is rotary evaporated to 0.5 mL and subsequently evaporated to dryness under nitrogen. The extracted residue is dissolved in 50% hexane : 50% dichloromethane. This solution is then submitted to florisil column chromatography and the Iprodione and RP-30228 residues are collected separately from the RP-32490 fraction. The solvent in each of the samples is evaporated to dryness, the residues are dissolved in a requisite volume of acetonitrile:water and analyzed by HPLC. Iprodione and RP-30228 are separated isocratically from coextracted artifacts using two HPLC columns in series and visualized with UV detection. The second degradate (RP-32490) is separated using a single column system. Linear regression analysis of Iprodione, RP-30228 and RP-32490 peak heights for samples and reference standards permits calculation of residue sample concentrations.

ANALYTICAL METHOD

Reagents

Acetonitrile, Burdick and Jackson, HPLC grade, UV cutoff @ 188 nm
Water, Barnstead NANOpure, HPLC grade
Sodium sulfate, Mallinckrodt, Analytical Reagent
Hydrochloric acid, Mallinckrodt, AR Select

Springborn Laboratories, Inc.

301

Report No. 91-1-3519

- Dichloromethane, Fisher Scientific, Optima
- Ethyl acetate, Burdick and Jackson, Residue grade
- Hexane, Burdick and Jackson, HPLC grade
- Toluene, Burdick and Jackson, UV cutoff @ 263 nm
- Acetone, Burdick and Jackson, Reagent grade
- Florisil, Fisher, 60-100 mesh, activated at 120 °C overnight
- Glass wool, washed with dichloromethane
- Iprodione, Lot No. EA 2002 SD8, 99.9% active ingredient, supplied by Rhone-Poulenc
- RP-30228, Lot No. EA 2025 RFI, 100% active ingredient, supplied by Rhone-Poulenc
- RP-32490, Lot No. EA 2026 RFI, 100% active ingredient, supplied by Rhone-Poulenc

Equipment

- Balance, Ohaus Galaxy 160, four-place analytical balance
- Balance, Ohaus MB 200, for moisture determination
- Nalgene bottles, Nalgene, HDPE, 500 mL
- Flasks, volumetric, assorted sizes
- Filter flasks, 1000 mL
- Buchner funnels, 11.0 cm diameter
- Glass fiber filters, Whatman, 934-AH, 11.0 cm diameter
- Separatory funnels, 250 mL
- Roundbottom flasks, assorted sizes
- Chromatographic columns, glass, 270 mm (length) x 10 mm (ID)
- Pipets, volumetric, assorted sizes
- Serum bottles, Wheaton, assorted sizes, with Teflon-lined lids and metal crimp caps
- Syringes, Hamilton, assorted sizes
- Rotary evaporator, Buchii Model R110, with vacuum pump, 40 °C water bath
- Shaker table, Eberbach, 0-500 rpm, or equivalent

Springborn Laboratories, Inc.

Springborn Laboratories, Inc.

302

Detailed Procedure

I. Preparation of Stock Solution

A. Iprodione (0.10 mg/mL)

1. Weigh 10.0 milligrams (A.I.) of Iprodione on an analytical balance.
2. Transfer the Iprodione to a 100-mL volumetric flask and dissolve to the mark with acetone.
3. Transfer the stock solution to a 100-mL amber serum vial and seal with a Teflon-lined crimp cap.
4. Store this stock solution in a refrigerator maintained at 4 °C.

B. RP-30228 and RP-32490 (0.10 mg/mL)

1. Weigh, dissolve, seal and refrigerate as previously described for Iprodione.

II. Quality Control Sample Fortification

A. Quality Control

1. Rinse all glassware with 50% dichloromethane : 25% acetone : 25% ethyl acetate prior to fortification.
2. To each 500 mL Nalgene bottle add approximately 50 grams (wet weight) of field sediment; see section IIIC.
3. To each sediment sample add 2.0 mL of 1.0 M HCl. Stir the sample in order to ensure homogeneity.
4. For a method validation or quality control sample, fortify each sample with Iprodione, RP-30228 and RP-32490 by volumetric addition of dilutions of the primary stock solutions.

NOTE: The fortification levels produced in the control field sediment samples for the method validation/recovery were 100, 50.0, and 20.0 ppb (three replicates at each level). An additional three field sediment samples were left unfortified and utilized as control samples.

Springborn Laboratories, Inc.

BEST AVAILABLE COPY

Springborn Laboratories, Inc.

30.3

III. Extraction

A. Field Sediment

1. To each Nalgene bottle, add approximately 50 grams (wet weight) of homogenized field sediment. Record the sediment weight to the nearest 0.01 grams.
NOTE: At this stage, a secondary subsample should be taken and weighed in order to determine the percent soil content; see section IIIC.
2. To each sediment sample add 2.0 mL of 1.0 M HCl. Stir the sample in order to ensure homogeneity.
3. Add 300 mL of 50% dichloromethane : 25% ethyl acetate : 25% acetone to the 50 gram sample.
4. Place the samples on a shaker table; adjust the rpm to 300 and allow the samples to shake for 50-55 minutes.
5. Filter the sample through a Whatman 934-AH 11.0 cm glass fiber filter. Add 25 mL of the extraction solvent to the Nalgene bottle, cap and shake for approximately 15 seconds; pour the rinse through the sediment (on the filter paper) and collect the eluent in the filter flask. Repeat this 25 mL rinse two additional times.
6. Dry the organic extract through sodium sulfate and collect the eluent in a 1000 mL roundbottom flask. Rinse the filter flask with 25 mL of the extraction solvent. Pour this rinse through the sodium sulfate and collect the eluent in the roundbottom flask. Repeat this rinse one additional time; combine rinses in the roundbottom flask.
7. Evaporate the sample to approximately 0.5 mL on a rotary evaporator.
8. Evaporate the remaining solvent under a gentle stream of nitrogen.
9. To the flask add 50 mL of acetonitrile and transfer the dissolved residues into a 250 mL separatory funnel. Rinse the flask with 10-15 mL of acetonitrile and transfer the solution into the separatory funnel. Repeat this rinse one additional time.

Springborn Laboratories, Inc

BEST AVAILABLE COPY

10. To the separatory funnel add 50 mL of hexane. Shake the separatory funnel for 1-2 minutes. Drain the lower acetonitrile layer into a 250 mL beaker; discard the hexane layer. Return the acetonitrile extract (with rinsing) to the separatory funnel. Repeat the 50 mL hexane back-partition; drain the acetonitrile into the original 1000 mL roundbottom flask. Discard the hexane layer.
11. Rotary evaporate the acetonitrile to approximately 0.5 mL and evaporate the remaining acetonitrile under a gentle stream of nitrogen.
12. To the flask add 10 mL of 50% dichloromethane : 50% hexane. Swirl the solution in order to dissolve the extracted residues. Proceed to Section III B, Column Chromatography.

B. Column Chromatography

1. Place (and tamp) a plug of glass wool in a glass chromatographic column (270 mm length x 10 mm ID) fitted with a teflon stopcock.
2. Add approximately 1 inch of sodium sulfate to the column.
3. To the column add enough 100% activated Florisil to form a 3 inch column (excluding the one inch of sodium sulfate). Add an additional 1 inch of sodium sulfate to the top of the Florisil column.
4. To the column add 25 mL of acetone (as a rinse); drain the column completely dry and discard the eluent.

Note: This step is necessary in order to clean the florisil of UV absorbing components. Add 25-30 mL of 50% dichloromethane : 50% hexane to the column and drain the solvent just into the upper sodium sulfate layer. Do not allow the florisil to go dry after this stage.

5. Transfer the dissolved residues from section A (step 8) to the column and open the stopcock fully. Rinse the roundbottom flask with 10 mL of 50% dichloromethane : 50% hexane and transfer the rinse onto the column. Drain the solvent into the upper sodium sulfate layer. Repeat this 10 mL

305

- rinse two additional times and drain into the sodium sulfate. Close the stopcock. Discard all rinses.
6. Add 20 mL of dichloromethane to the column and drain into the sodium sulfate. Discard the eluent.
 7. Place a 50 mL roundbottom flask under the column. Add 20 mL of 5% acetone : 95% dichloromethane to the column and drain into the sodium sulfate. Collect this fraction which contains Iprodione and RP-30228.
 8. To the column add 20 mL of 15% acetone : 85 % dichloromethane; drain into the sodium sulfate. Collect this fraction separately which contains RP-32490.
 9. Rotary evaporate the column eluents (separately) to approximately 0.5 mL. Evaporate the solution to dryness under a gentle stream of nitrogen.
 10. Dissolve the residues in a requisite volume of 50% acetonitrile:50% HPLC grade water. Proceed to section IV for HPLC analysis of Iprodione, RP-30228 and,RP-32490 residues.

C. Percent Soil Determination

1. Weigh approximately 10 grams of the sediment sample into an aluminum weigh-dish.
2. Place the sample into an Ohaus MB 200 moisture determination balance and close the cover.
3. Bake the sample at 205 °C for 10 minutes.
4. Document the percent solids read from the LED.

IV. High Pressure Liquid Chromatography

A. Method:

Fractions containing Iprodione and RP-30228 are analyzed isocratically on a two-column system. Fractions containing RP-32490 are analyzed isocratically on a one column system. Both analyses utilize UV detection at 200 nm.

Springborn Laboratories, Inc.

BEST AVAILABLE COPY

Springborn Laboratories, Inc.

306

Report No. 91-3-3699

Page 9 of 27

Report No. 91-1-3619

Instrumental conditions for the analysis of lprodione and RP-30228 are as follows:

Instrument: Waters Model 510 liquid chromatograph solvent pump equipped with Waters Intelligent Sample Processor Model 710 B, ABI Model 783 Variable Wavelength Detector and Hewlett-Packard Model 3396A integrator.

Column: Phenomenex Ultramex 5 μ m ODS, 250 mm (length) x 4.6 mm ID; 2 columns in series

Mobile Phase: 80% acetonitrile:20% HPLC grade water

Mobile Phase Flowrate: 1.5 mL/minute

Pressure: ca 2500 psi

Chartspeed: 0.3 cm/minute

Injection Volume: 125 μ L

Wavelength: 200 nm

Sensitivity: 0.010 AUFS

Rise Time: 0.1 seconds

Attenuation: 2⁶

Threshold: 6

Peak Width: 0.10 seconds

Area Reject: 500

Instrumental conditions for the analysis of RP-32490 are as follows:

Instrument: Waters Model 510 liquid chromatograph solvent pump equipped with Waters Intelligent Sample Processor Model 710 B, ABI Model 783 Variable Wavelength Detector and Hewlett-Packard Model 3396A integrator.

Column: Phenomenex Ultramex 5 μ m ODS, 250 mm (length) x 4.6 mm ID

Mobile Phase: 60% acetonitrile:40% HPLC grade water

Mobile Phase Flowrate: 1.2 mL/minute

Springborn Laboratories, Inc.

BEST AVAILABLE COPY

Report No. 91-1-3619

Pressure:	ca 1000 psi
Chartspeed:	0.3 cm/minute
Injection Volume:	100 μ L
Wavelength:	200 nm
Sensitivity:	0.010 AUFS
Rise Time:	0.1 seconds
Attenuation:	2 ⁶
Threshold:	6
Peak Width:	0.15 seconds
Area Reject:	500

B. Analysis

1. Prepare standard solutions containing both lprodione and RP-30228 and separately RP-32490. Standard solution concentrations used for the recovery study were 125, 250, 500, 750 and 1000 μ g/L.
 2. Inject 125 μ L (or 100 μ L for RP-32490) of the 125 μ g/L standard solution. Adjust the attenuation so that the peak signal results in at least a fifteen percent deflection from the baseline.
 3. Inject 125 μ L of each of the mixed standards, document the peak heights, and determine the correlation coefficient of the line. Proceed to step 4 if the correlation coefficient is greater than or equal to 0.985.
 4. Inject 125 μ L of several samples.
 5. Identify each analyte by its retention time and document the respective peak heights.
 6. After each set of samples, reinject 125 μ L of each of the mixed standards and document the peak heights.
 7. Construct a standard curve for each analyte (using all standard results) by plotting peak height observed versus the concentration (μ g/L) of the standard injected.
- B. The standard linear regression analysis for lprodione, RP-30228 and RP-32490 is used to determine the concentration in each sample.

Springborn Laboratories, Inc.

Springborn Laboratories, Inc.

308

9. In order to determine the analytical result for each sample, the following equation is used:

$$\text{Analytical Result (ppb)} = A \times \text{D.F.}$$

where:

Analytical Result = concentration of Iprodione, RP-30228 or RP-32490

A = concentration ($\mu\text{g/L}$) of sample from the regression analysis

D.F. = dilution factor, ratio of the final volume (mL) of the extracted sample to the initial mass (g) of sample