


~~This method cancels and replaces PPRAM 2A dated November 1960.~~

1. SCOPE

The analytical procedures described are suitable for the determination of residues of paraquat (I, 1,1'-dimethyl-4,4'-bipyridinium ion) in soil.



The limit of determination of the method is 0.05 $\mu\text{g kg}^{-1}$.

2. METHOD SUMMARY

The sample is boiled in 6N sulphuric acid solution. The filtered digest is percolated through a column of cation-exchange resin which retains the paraquat and some of the natural soil constituents. The column is washed with dilute hydrochloric acid, 2.5% ammonium chloride solution and water; the paraquat is eluted with saturated ammonium chloride solution. A portion of the column eluent is treated with sodium dithionite in alkali. This reduces paraquat to a free radical the light absorption of which is measured with a spectrophotometer.


3. REAGENTS

- Sodium chloride solution - saturated
- Ammonium chloride solution - 2.5% (w/v) and saturated
- Sulphuric acid - concentrated (13M)
- Hydrochloric acid solution (2N)
- Cation-exchange resin: Duolite C225 (SRC 14) chromatographic grade resin, 52-100 mesh, 0.68-0.85 water regain, sodium form.
- Standard paraquat solutions:

- Stock solution (250 $\mu\text{g/ml}$ of paraquat)

Dissolve 0.0864 g pure paraquat dichloride, $\text{C}_{14}\text{H}_{14}\text{N}_2\text{Cl}_2$ (mol. wt., 257.2; 72.46 cation), in saturated ammonium chloride solution and make up to 250 ml with saturated ammonium chloride solution. Paraquat salts are hygroscopic; they should be dried at 100°C for 5 hours and cooled in a desiccator before use.

BEST AVAILABLE COPY



(ii) Working solutions

Make dilutions of the stock solution to give a range of working solutions. Concentrations of these working solutions will vary depending on the level of residues in the samples analysed.

These solutions are stable under normal laboratory conditions provided that they are not exposed to sunlight for long periods.

- (g) Sodium dithionite solution, 0.2% (w/v) in 0.1M sodium hydroxide. This solution should, if possible be used immediately, and must not be used more than half an hour after preparation. When preparing this solution do not shake vigorously.

Solid sodium dithionite is unstable in the presence of moisture, and should therefore be stored in a tightly sealed plastic container.

- (h) Octan-2-ol

4. SAFETY CONSIDERS

The following information is included as an indication to the analyst of the nature and hazards of the reagents used in this procedure. If in doubt, consult the appropriate safety manual (eg ICI Laboratory Safety Manual) containing recommendations and procedures for handling chemicals, and a monograph such as 'Hazards in the Chemical Laboratory', edited by G D Muir, The Chemical Society, London.


(a) HYDROCHLORIC ACID

Irritant vapour
Corrosive - causes burns
Avoid breathing vapour.
Prevent contact with eyes and skin
(TLV 7 mg m⁻³ as HCl).

(b) SULPHURIC ACID - concentrated

Corrosive - causes burns
Prevent contact with skin and eyes
Do not put water into container
(TLV 1 mg m⁻³)

BEST AVAILABLE COPY



(c) **PARAQUAT**

Toxic by ingestion

Harmful dust.

Avoid contact with eyes, skin and mouth. Avoid breathing dust. Wash hands and exposed skin before meals and after work.

Exposure to paraquat should be regarded as a dire emergency and action taken immediately. Details of remedial action/antidotes should be available in the laboratory.

(d) **OCEAN-2-OL**

Harmful vapour

Harmful if taken internally

Highly flammable

Avoid breathing vapour or contact with skin and eyes
(TLV 2.60 mg m⁻³)

5. **APPARATUS**

- (a) Equipment which can be used for the initial preparation of samples is, metal trays for drying soils.
- (b) Heating mantles for the extraction procedure (e.g. Electrothermal Heating Unit available from Electrothermal Engineering Ltd, London E7 UK).
- (c) Boiling flasks - 500 ml round bottom flasks with 234 necks fitted with water cooled reflux condensers.
- (d) Glass columns for chromatography of 1.0 cm i.d. and 50 cm long (25 ml burettes are suitable).
- (e) Scanning Spectrophotometer eg, Perkin Elmer Lambda 5 UV/VIS Spectrophotometer.
- (f) Sieve, 2-3 mm mesh size.


6. **PROCEDURE**

6.1 **Extraction and Chromatographic Separation**

- (a) Thoroughly mix the soil sample and weigh about 100 g into a shallow metal tray or glass petri dish. Dry the soil sample at room temperature for several days.
- (b) Grind the dry soil sample to pass through a 2-3 mm sieve. Weigh an aliquot (25 g) of the ground soil into a 500 ml round bottomed flask and add water (65 ml) followed by 10M sulphuric acid (35 ml).

Note: **CAUTION!** The sulphuric acid must be added to the aqueous solution with care.

BEST AVAILABLE COPY



- (c) After adding anti-bumping granules and octan-2-ol (1 ml) which acts as an anti-foaming agent, place the 500 ml flask containing the sample on a heating mantle. Attach a water-cooled reflux condenser and heat to boiling. Swirl the flask contents occasionally to minimize local overheating and charring until the solution is boiling steadily.

It is important to avoid local overheating and charring and this can be done by ensuring that only the bottom quarter of the flask makes contact with the heating mantle (which should be oversized).

- (d) Boil under reflux for 5 hours and allow to cool. The solution can be left overnight at this stage.
- (e) While the samples are being refluxed the ion-exchange columns are prepared as follows: Wash 5.0 g of resin with water into a burette (25 ml) containing a glass wool plug placed near the stopcock. Pass successively through the column at the rate of 5 ml/min. saturated sodium chloride solution (20 ml) and water (50 ml). Prepare a separate column for each sample.
- (f) Wash the reflux condenser attached to the boiling flask with water (50 ml) into the cooled contents of the 500 ml round-bottomed flask from 5.1 (d) above. Dilute the flask contents to near 500 ml by the addition of water and filter the solution by suction through 2 Whatman No.5 filter papers. Soak the filter pad dry and wash the filter twice with water (100 ml), allowing the first portion to be sucked through before adding the second.
- (g) Transfer the filtrate to a 1 litre separating funnel. Allow the solution to percolate through a prepared resin column from 5.1 (e) above at a flow rate of 5-10 ml/min.
- (h) Remove the funnel and wash the column at a flow rate of 3-4 ml/min successively with water (25 ml), 2N-hydrochloric acid (100 ml), water (25 ml), 2.5% (w/v) ammonium chloride solution (50 ml) and water (25 ml). (The process can be left overnight provided the resin column has been covered with water).
- (i) Elute the paraquat from the column with saturated ammonium chloride solution at a flow rate of about 1 ml/min. Collect the first 50 ml of the eluent in a 50 ml volumetric flask and mix.

5.3 The recovery of the paraquat from the resin column will be adversely affected if the flow rate of the eluent exceeds 1.0 ml/min.

BEST AVAILABLE COPY

6.2

Determination

- (a) Pipette an aliquot (10.0 ml) of the eluent into a 15 ml glass-stoppered test-tube. Add by pipette 0.2% (w/v) sodium dithionite solution (2.0 ml) and mix by gently inverting the tube once.
- (b) Within 5 minutes of adding the sodium dithionite, use a recording spectrophotometer to record the spectrum of the solution in a 4 cm path length cell over the range 360-430 nm. against a reference solution prepared from saturated ammonium chloride (10.0 ml) and sodium dithionite (2.0 ml).
- (c) Draw a baseline as a tangent to the curve from the valley in the region of 390 nm. Measure the height of the peak above the baseline at 396 nm.
- (d) Draw a calibration curve relating the peak height at 396 nm (nm above the baseline) to the concentration of paraquat in $\mu\text{g/ml}$.

N.B When using a spectrophotometer with a derivative function, operation in the 2nd derivative mode will give an enhanced response to paraquat. See Appendix 1 for comparison of spectra and methods of measuring peak heights.

7.

CALCULATION

Read off from the prepared calibration curve, using the peak height measured at 396 nm, the concentration ($\mu\text{g/ml}$) of paraquat in the eluent.

N.B It has been found useful to use the linear regression function on an electronic calculator to produce the calibration curve.

Then paraquat concentration in sample ($\mu\text{g/kg}$)

$$= \frac{\text{Volume of eluent, ml from column} \times \text{Concentration in eluent } (\mu\text{g/ml})}{\text{Weight of sample, (g)}}$$

To correct the paraquat concentration for the experimental recovery multiply by 100 and divide by the % of recovery figure.

BEST AVAILABLE COPY