

## Introduction

### Objective of Study

The purpose of this project is to analyze water samples collected from Arkansas and California for orthosulfamuron and its major metabolites during a method validation using a liquid chromatography-mass spectrometry/mass spectrometry (LC-MS/MS) system. The target limit of quantitation (LOQ) is 0.17 ng/mL.

## Section 1 Analytical Summary

### Section 1.1 Method Summary

The residue method was developed and validated (Appendix 1) for orthosulfamuron and its major metabolites in rice paddy source water obtained from the orthosulfamuron aquatic field dissipation test sites located in Arkansas and California (Waterborne Study Number WEI 642.03). The room temperature water sample was buffered with 10 mM disodium hydrogen phosphate if the pH range was not between 8 and 9. After mixing the water sample, a 10 mL aliquot was transferred into a plastic centrifuge tube and fortified if necessary. The sample was shaken at ~ 180 excursions per minute for 5 minutes. A ~ 1.5 mL aliquot was filtered through a syringe filter into an autosampler vial, capped and analyzed by LC/MS/MS.

### Section 1.2 Method Validation and Limit of Quantitation

The method was successfully validated by analysis of blank untreated control (UTC) water samples fortified with 0.170 and 1.70 ng/mL of each analyte. These data support a method limit of quantitation (LOQ) of 0.170 ng/mL for each analyte.

### Section 1.3 Determination of Limit of Detection


The limit of detection (LOD) (Appendix 1) is the lowest analyte concentration for a measurement statistically different from that of a blank UTC sample. The LOD for orthosulfamuron and its major metabolites in water were calculated by multiplying **SD**, the standard deviation of the analyte recovery measurements at the target LOQ, by  $t_{0.99}$ , the appropriate one-tailed Student's *t* statistic. See Section 4 for the LOD calculations for orthosulfamuron and its major metabolites. The LOD was verified by analysis of UTC water samples fortified at the target LOQ level of 0.170 ng/mL each of orthosulfamuron and its major metabolites. The equation for calculating the LOD is given below.

$$\text{LOD} = (t_{0.99} \times \text{SD})$$


### Section 1.4 Analytical Standards Used in this Study

Solutions of the analytical reference standards were prepared in acetonitrile, or in a mixture of acetonitrile and water with 10 mM disodium hydrogen phosphate, and were corrected for purity during initial standard preparation. Standard solutions and linearity curve solutions were stored in a laboratory freezer (BLDG2WF1) at an average temperature of  $\leq -18^{\circ}\text{C}$ . The standards for orthosulfamuron and its metabolites were obtained from ABC Laboratories, Inc.


Common Name: DB Amine  
 Molecular Weight: 164.20 g/mol  
 Appearance: off-white solid  
 Lot Number: D16M  
 Purity: 99.7%  
 Expiration Date: 2 years from date of C of A  
 Storage Conditions: Room Temperature  
 Chemical Name (IUPAC): 2-Amino-N,N-dimethylbenzamide  
 CAS#: 6526-66-5  
 Structure:




Common Name: O-Desmethyl-Orthosulfamuron  
 Molecular Weight: 410.41 g/mol  
 Appearance: white solid  
 Lot Number: 81481-65  
 Purity: 97.6  
 Expiration Date: 2 years from date of C of A  
 Storage Conditions: ≤-20°C  
 Chemical Name (IUPAC): 1-(4-hydroxy-6-methoxypyrimidin-2-yl)-3-[2-(dimethoxycarbonyl)phenylsulfamoyl]-urea  
 Structure:




Common Name: DBS Acid Ammonium Salt  
 Molecular Weight: 261.30 g/mol  
 Appearance: white solid  
 Lot Number: 81481-1-39-2  
 Purity: 97.8  
 Expiration Date: 2 years from date of C of A  
 Storage Conditions: ≤-20°C  
 Chemical Name (IUPAC): (2-dimethylcarbamoylphenyl)-sulfamic acid ammonium salt  
 Structure:




Common Name: DOP Amine  
 Molecular Weight: 155.15 g/mol  
 Appearance: white solid  
 Lot Number: S61826V  
 Purity: 99.8  
 Expiration Date: 2 years from date of C of A  
 Storage Conditions: Room Temperature  
 Chemical Name (IUPAC): (2-amino-4,6-dimethoxy)-pyrimidine  
 CAS#: 36315-01-2  
 Structure:




Common Name: DBS Amide  
 Molecular Weight: 243.28 g/mol  
 Appearance: pale-orange solid  
 Lot Number: 81481-1-16  
 Purity: 91.9  
 Expiration Date: 2 years from date of C of A  
 Storage Conditions: ≤-20°C  
 Chemical Name (IUPAC): (2-dimethylcarbamoylphenyl)-sulfanamide  
 Structure:




Common Name: DOP Urea  
 Molecular Weight: 198.18 g/mol  
 Appearance: off-white Solid  
 Lot Number: 81481-1-33-4  
 Purity: 99.8  
 Expiration Date: 2 years from date of C of A  
 Storage Conditions: ≤-20°C  
 Chemical Name (IUPAC): (4,6-dimethoxy-2-pyrimidinyl)-urea  
 CAS#: 151331-81-6  
 Structure:



Common Name: O-Desmethyl DOP urea  
Molecular Weight: 184.15 g/mol  
Appearance: white solid  
Lot Number: 81481-1-72-1  
Purity: 97.7  
Expiration Date: 2 years from date of C of A  
Storage Conditions: ≤-20°C  
Chemical Name (IUPAC): N-(4-hydroxy-6-methoxy-2-pyrimidin-2-yl)-urea  
CAS#: 888225-63-6  
Structure:



Common Name: Orthosulfamuron  
Molecular Weight: 424.43 g/mol  
Appearance: white solid  
Lot Number: 81481-1-41-1  
Purity: 99.1  
Expiration Date: 2 years from date of C of A  
Storage Conditions: ≤-20°C  
Chemical Name (IUPAC): 1-(4,6-Dimethoxy-2-pyrimidinyl)-3-[2-(dimethylcarbamoyl)phenylsulfamoyl]-urea  
CAS#: 213464-77-8  
Structure:



### Section 1.5 Instrument Detector Response – Range and Linearity

Calibration standards and samples were analyzed using LC/MS/MS. Calibration curves and recovery values were calculated using Analyst 1.5.1 or higher data handling software using linear regression with 1/x weighting.

The standards were fit to the linear equation:  $Y = MX + B$  with 1/x weighting

where: X is the concentration of the reference standard  
M is the calibration line slope  
B is the calibration line intercept  
Y is the native peak area

The relative response of the LC/MS/MS to orthosulfamuron and its metabolites was linear over a range of 0.10 ng/mL to 2.5 ng/mL using matrix-matched standards. The correlation coefficients (r value in Analyst file) of the linearity curves were all ≥0.99 (see worksheets in Section 5). Matrix-matched standards were necessary due to the difference in detector response for a solvent standard compared to a matrix sample. An enhancement effect was observed when analyzing matrix samples and that necessitated the use of matrix-matched standards to have comparable responses to accurately calculate the residues. See Appendix 2 for an example of the instrument conditions that were used.

**Appendix 1****ANALYTICAL METHOD VALIDATION FOR THE DETERMINATION OF  
ORTHOSULFAMURON AND ITS MAJOR METABOLITES IN WATER USING  
LC-MS/MS**

Prepared by: Walter R. Vandaveer, Research Scientist  
SynTech Research Laboratory Services

Date: 04/13/2015

**1.0 OBJECTIVE**

The purpose of this project is to analyze water samples collected from Arkansas and California for orthosulfamuron and its major metabolites during a method validation using a liquid chromatography-mass spectrometry/mass spectroscopy (LC-MS/MS) system. The target limit of quantitation (LOQ) is 0.17 ng/mL.

**2.0 EXPERIMENTAL DESIGN**

- 2.1 Test solutions from the study will be analyzed.
- 2.2 The analytical method will be developed and validated before the sample analysis.
- 2.3 O-Desmethyl DOP urea, Lot Number 81481-1-72-1, 97.7% purity, expiration date 11/26/2016; DBS acid ammonium salt, Lot Number 81481-1-39-2, 97.8% purity, expiration date 11/26/2016; DBS amide, Lot Number 81481-1-16, 91.9% purity, expiration date 11/26/2016; DOP amine, Lot Number S61826V, 99.8% purity, expiration date 11/26/2016; DBS amine, Lot Number D16M, 99.7% purity, expiration date 11/26/2016; O-Desmethyl Orthosulfamuron, Lot Number 81481-65, 97.6% purity, expiration date 11/26/2016; DOP urea, Lot Number 81481-1-33-4, 99.8% purity, expiration date 11/26/2016; Orthosulfamuron, Lot Number 81481-1-41-1, 99.1% purity, expiration date 11/26/2016 will be supplied by ABC Laboratories, Inc.
- 2.4 Method validation will be completed prior to the analysis of study samples. This will consist of a solvent blank, two controls, and 5 spiked controls at the target LOQ (0.17 ng/mL) and 5 spiked controls at 10 x the target LOQ (1.7 ng/mL) with recovery determined using a calibration curve. All procedures developed during the method validation phase of the project will be employed during the study sample analysis. A recovery spike for orthosulfamuron and its major metabolites will be required to monitor methodology and individual performance. The acceptable recoveries are 70 - 120%. Results will be considered acceptable if the average recovery value (of the two or more concurrent recovery samples) falls within this range. If the lab spike falls outside of this range, the results may be rejected.
- 2.5 All analytical results will be reported to the study director.

**3.0 METHODOLOGY****3.1 Preparation of Calibration Standards**

3.1.1 Weigh ~ 1 mg of each standard into a 50 mL volumetric flask, dilute to volume with the appropriate solvent, and mix well by inversion (see the following table) to create a 20 µg/mL standard stock solution.

Compound Name	Lot No.	Amount (mg)	Purity (%)	Solvent	Volume (mL)	Target Concentration (µg/mL)
O-Desmethyl DOP urea	81481-1-72-1	1.03	97.7	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	50	20
DBS acid ammonium salt	81481-1-39-2	1.03	97.8	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	50	20
DBS amide	81481-1-16	1.09	91.9	ACN	50	20
DOP amine	S61826V	1.0	99.8	ACN	50	20
DBS amine	D16M	1.0	99.7	ACN	50	20
O-Desmethyl Orthosulfamuron	81481-65	1.03	97.6	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	50	20
DOP urea	81481-1-33-4	1.0	99.8	ACN	50	20
Orthosulfamuron	81481-1-41-1	1.01	99.1	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	50	20

3.1.2 Pipette 1.0 mL of each ~ 20 µg/mL standard stock solution into a 100 mL volumetric flask, dilute to volume with the appropriate solvent, and mix well by inversion (see the following table) to create a 0.2 µg/mL standard stock solution.

Compound Name	Lot No.	Amount (mL)	Solvent	Volume (mL)	Target Concentration (µg/mL)
O-Desmethyl DOP urea	81481-1-72-1	1.0	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	100	0.20
DBS acid ammonium salt	81481-1-39-2	1.0	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	100	0.20
DBS amide	81481-1-16	1.0	ACN	100	0.20
DOP amine	S61826V	1.0	ACN	100	0.20
DBS amine	D16M	1.0	ACN	100	0.20
O-Desmethyl Orthosulfamuron	81481-65	1.0	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	100	0.20
DOP urea	81481-1-33-4	1.0	ACN	100	0.20
Orthosulfamuron	81481-1-41-1	1.0	ACN/H <sub>2</sub> O with 10mM Na <sub>2</sub> HPO <sub>4</sub> (50/50, %/%)	100	0.20

- 3.1.3 Pipette ~ 3 mL of each ~ 0.2 µg/mL standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 25 ng/mL mixed standard stock solution. The exact volume of standard stock solution will need to be calculated based on the actual standard stock concentration for each analyte.
- 3.1.4 Pipette 5.0 mL of the 25 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 5 ng/mL mixed standard stock solution.
- 3.1.5 Pipette 0.5 mL of the 5 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 0.1 ng/mL calibration standard solution.
- 3.1.6 Pipette 0.75 mL of the 5 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 0.15 ng/mL calibration standard solution.
- 3.1.7 Pipette 1.0 mL of the 5 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 0.20 ng/mL calibration standard solution.

- 3.1.8 Pipette 2.5 mL of the 5 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 0.50 ng/mL calibration standard solution.
- 3.1.9 Pipette 1.0 mL of the 25 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 1.0 ng/mL calibration standard solution.
- 3.1.10 Pipette 1.5 mL of the 25 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 1.5 ng/mL calibration standard solution.
- 3.1.11 Pipette 2.0 mL of the 25 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 2.0 ng/mL calibration standard solution.
- 3.1.12 Pipette 2.5 mL of the 25 ng/mL mixed standard stock solution into a 25 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 2.5 ng/mL calibration standard solution.

### 3.2 Preparation of Test Solutions

- 3.2.1 If necessary, thaw the water sample(s) thoroughly. The day before the analysis, remove the sample from the freezer and place it in a refrigerator overnight. The day of analysis, remove the sample from the refrigerator and place it on a benchtop and allow to warm to room temperature.
- 3.2.2 Determine the pH of the water sample using a calibrated pH probe.
- 3.2.3 Desired pH range is 8 to 9 to assist with stabilizing the analytes in the water sample.
- 3.2.4 If necessary, buffer the water sample to above pH 8 using 10 mM disodium hydrogen phosphate: for example add approximately 1.42 g into 1 L of water and shake vigorously, this will yield a buffered water sample with 10 mM disodium hydrogen phosphate.  
Observation: water sample goes from clear to a bit cloudy after the addition of disodium hydrogen phosphate
- 3.2.5 Verify the pH using a calibrated pH probe.
- 3.2.6 Label 12 x 50 mL plastic centrifuge tubes: 2 controls, 5 for LOQ, and 5 for 10 x LOQ
- 3.2.7 Shake the water sample and then add 10 mL of buffered water sample into each of the 12 labeled centrifuge tubes.
- 3.2.8 The 2 tubes labeled as control are left unfortified.

- 3.2.9 To prepare 7 fortified recovery QC samples at the target LOQ (0.17 ng/mL), pipette 352  $\mu$ L of the 5 ng/mL mixed standard stock solution into each labeled centrifuge tube containing 10 mL of buffered water sample, for a total volume of 10.352 mL, and mix well by inversion.
- 3.2.10 To prepare 5 fortified recovery QC samples at 10 x the target LOQ (1.7 ng/mL), pipette 730  $\mu$ L of the 25 ng/mL mixed standard stock solution into the labeled centrifuge tube containing 10 mL of buffered water sample, for a total volume of 10.730 mL, and mix well by inversion.
- 3.2.11 Shake all samples (12 tubes) at ~ 180 excursions per minute for 5 minutes.
- 3.2.12 Approximately 1.5 mL of each sample will be aliquoted and filtered through a 0.45  $\mu$ m syringe filter into an autosampler vial, cap vial and proceed to LC-MS/MS analysis.
- 3.2.13 Analyze all samples using the parameters shown in Section 3.4.

### 3.3 Calculation of Orthosulfamuron and its Major Metabolites Concentrations and Recoveries

Generate a calibration curve with each set of samples analyzed by plotting the average response of the peak areas versus the concentration of each standard solution.

Calibration curve equation:  $y = mx + b$ , where  $b$  = intercept,  $m$  = slope of the curve,  $x$  = sample concentration,  $y$  = response of each analytes peak area.

$$\text{Calculated Conc. (ng/mL)} = \frac{[(\text{Peak Area Response}) - (\text{intercept})]}{\text{Slope of Calibration Curve}} \times \text{DF}$$

Where: DF = Dilution Factor

$$\text{Percent Recovery} = \frac{\text{Analyte Conc., ng/mL}}{\text{Nominal Conc., ng/mL}} \times 100$$

### 3.4 Instrument Conditions

Note: The parameters provided below are suggested parameters for the method. These parameters may be adjusted to compensate for column-to-column and system-to-system variability.



<b>LC-MS/MS System Conditions for Analysis of Orthosulfamuron and Metabolites</b>				
Liquid Chromatograph (LC): Shimadzu LC-20AD HPLC Pumps and CBM 20A Controller				
Autosampler: Shimadzu SIL-20A/HT				
Valco 2-position switching valve (optional)				
Shimadzu CTO-20A Column Oven				
Column: Phenomenex Luna C18(2), 150 mm x 2.0 mm ID, 5 µm particle size				
Guard Column: C18, 4 mm x 2.0 mm ID				
Mass Spectrometer (MS): Sciex Triple Quad LC/MS/MS System				
Software: Analyst Software (Version 1.6.1)				
<b>LC Conditions</b>				
Oven Temperature:	35°C			
Initial Flow:	0.250 mL/min			
Injection Volume:	25 µL			
Mobile Phase:	Organic = 10 mM Ammonium Formate in water/ACN with 0.1% formic acid (10:90, v/v)			
	Aqueous = 10 mM Ammonium Formate in water with 0.1% formic acid			
	Rinse = ACN/Water (1:1, v/v)			
Gradient Method:				
	<u>Time (min)</u>	<u>%Aqueous</u>	<u>%Organic</u>	<u>Flow (mL/min)</u>
	0.00	100	0	0.250
	2.00	100	0	0.250
	7.50	10	90	0.250
	10.00	10	90	0.250
	10.10	100	0	0.250
	15.00	100	0	0.250
Approx. Retention Times (min):				
O-Desmethyl DOP urea	-	6.6		
DBS acid ammonium salt	-	6.8		
DBS amide	-	7.3		
DOP amine	-	7.7		
DBS amine	-	7.5		
O-Desmethyl Orthosulfamuron	-	7.5		
DOP urea	-	7.7		
Orthosulfamuron	-	9.0		
<b>Autosampler Conditions</b>				
Rinsing Volume:	1000 µL			
Needle Stroke:	52 mm			
Rinsing Speed:	35 µL/sec			
Sampling Speed:	2.0 µL/sec			
Purge Time:	25.0 min			
Rinse Dip Time:	0 sec			
Rinse Mode:	Before and after aspiration			
Control Vial Needle Stroke:	52 mm			
Rinse Method:	Rinse Port Only			

Rinse Time:	1 sec
<b>Mass Spectrometer Conditions for LC-MS/MS Analysis</b>	
<b>Period 1:</b>	
Analyte:	O-Desmethyl DOP urea
Scan Type:	MRM
Ion Source:	Turbospray
Polarity:	Positive
Q1 Mass (amu):	185.000
Q3 Mass (amu):	142.100
Dwell (msec):	20.00
Period Time (min):	0.0 to 6.60
Declustering Potential (DP):	94.0 V
Entrance Potential (EP):	7.0 V
Collision Energy (CE):	23.0 V
Collision Cell Exit Potential (CXP):	6.0 V
Curtain Gas (CUR):	35.0
Temperature (TEM):	550°C
Ion Source Gas 1 (GS1):	35.0
Ion Source Gas 2 (GS2):	55.0
Collision Gas (CAD):	7.0
Ionspray Voltage (IS):	3000 V
<b>Period 2:</b>	
Analyte:	DBS acid ammonium salt
Scan Type:	MRM
Ion Source:	Turbospray
Polarity:	Negative
Q1 Mass (amu):	243.165
Q3 Mass (amu):	80.104
Dwell (msec):	20.00
Period Time (min):	6.60 to 6.85
Declustering Potential (DP):	-40.0 V
Entrance Potential (EP):	-5.0 V
Collision Energy (CE):	-21.0 V
Collision Cell Exit Potential (CXP):	-12.0 V
Curtain Gas (CUR):	35.0
Temperature (TEM):	550°C
Ion Source Gas 1 (GS1):	35.0
Ion Source Gas 2 (GS2):	55.0
Collision Gas (CAD):	7.0
Ionspray Voltage (IS):	-4400 V

Period 3:						
Analyte	DBS amide	DOP amine	DBS amine	O-Desmethyl Orthosulfamuron	DOP urea	Orthosulfamuron
Scan Type	MRM	MRM	MRM	MRM	MRM	MRM
Ion Source	Turbospray	Turbospray	Turbospray	Turbospray	Turbospray	Turbospray
Polarity	Positive	Positive	Positive	Positive	Positive	Positive
Q1 Mass (amu)	244.100	156.017	164.740	411.074	199.200	425.089
Q3 Mass (amu)	165.000	99.700	119.800	227.000	156.000	199.000
Dwell (msec)	20.00	20.00	20.00	20.00	20.00	20.00
Period Time (min)	6.85 to 15.0	6.85 to 15.0	6.85 to 15.0	6.85 to 15.0	6.85 to 15.0	6.85 to 15.0
Declustering Potential (DP)	39.0 V	50.0 V	55.0 V	51.0 V	60.0 V	56.0 V
Entrance Potential (EP)	6.0 V	6.0 V	7.0 V	3.0 V	5.0 V	4.0 V
Collision Energy (CE)	13.0 V	21.0 V	13.0 V	17.0 V	19.0 V	15.0 V
Collision Cell Exit Potential (CXP)	6.0 V	50.0 V	4.0 V	28.0 V	10.0 V	22.0 V
Curtain Gas (CUR)	35.0	35.0	35.0	35.0	35.0	35.0
Temperature (TEM)	550.0°C	550.0°C	550.0°C	550.0°C	550.0°C	550.0°C
Ion Source Gas 1 (GS1)	35.0	35.0	35.0	35.0	35.0	35.0
Ion Source Gas 2 (GS2)	55.0	55.0	55.0	55.0	55.0	55.0
Collision Gas (CAD)	7.0	7.0	7.0	7.0	7.0	7.0
Ionspray Voltage (IS)	3000 V	3000 V	3000 V	3000 V	3000 V	3000 V

Divert Valve Program: 0.0 to 5.5 min – To Waste  
5.5 to 10.5 min – To Instrument  
10.5 to 15.0 min – To Waste

Note: The indicated LC-MS/MS parameters are guidelines and should be optimized for the instrument and column actually used. Instrument parameters and mobile phase compositions may be adjusted to improve separation from interfering peaks.

### 3.5 Revisions

None

Method Amendment

---

**Modification Number: 1**  
**Study Number: 292SRLS14R01**

---

**Method Title:** Analytical Method Validation for the Determination of Orthosulfamuron and its Major Metabolites in Water using LC-MS/MS

**Reason for the Change/addition and impact on the Method:**

The method modifications address typographical errors, using matrix matched standards, and different spiking concentrations for recovery QC samples. No adverse impact on the method.

**Location and Description of Changes:**

**2.0 EXPERIMENTAL DESIGN**

2.4 Clarify that there shall be at least 7 spiked controls at the target LOQ and at least 5 spiked controls at 10 x the target LOQ

**3.0 METHODOLOGY**

**3.1 Preparation of Calibration Standards**

3.1.2 Pipette 1.0 mL of each ~ 20 µg/mL standard stock solution into one 10 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 2.0 µg/mL (2000 ng/mL) mixed standard stock solution.

3.1.3 Pipette 0.1 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 4.0 ng/mL spiking matrix matched calibration standard solution.

3.1.4 Pipette 0.15 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 6.0 ng/mL spiking matrix matched calibration standard solution.

3.1.5 Pipette 0.2 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 8.0 ng/mL spiking matrix matched calibration standard solution.

3.1.6 Pipette 0.5 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 20 ng/mL spiking matrix matched calibration standard solution.

3.1.7 Pipette 1.0 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and

### Method Amendment

mix well by inversion to create a 40 ng/mL spiking matrix matched calibration standard solution.

3.1.8 Pipette 1.5 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 60 ng/mL spiking matrix matched calibration standard solution.

3.1.9 Pipette 2.0 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 80 ng/mL spiking matrix matched calibration standard solution.

3.1.10 Pipette 2.5 mL of the 2000 ng/mL mixed standard stock solution into a 50 mL volumetric flask, dilute to volume with ACN/H<sub>2</sub>O with 10 mM Na<sub>2</sub>HPO<sub>4</sub> (50/50, %/%), and mix well by inversion to create a 100 ng/mL spiking matrix matched calibration standard solution.

3.1.11 Pipette 25 µL of the 4.0 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 0.1 ng/mL matrix matched calibration standard solution.

3.1.12 Pipette 25 µL of the 6.0 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 0.15 ng/mL matrix matched calibration standard solution.

3.1.13 Pipette 25 µL of the 8.0 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 0.20 ng/mL matrix matched calibration standard solution.

3.1.14 Pipette 25 µL of the 20 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 0.50 ng/mL matrix matched calibration standard solution.

3.1.15 Pipette 25 µL of the 40 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 1.0 ng/mL matrix matched calibration standard solution.

3.1.16 Pipette 25 µL of the 60 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 1.5 ng/mL matrix matched calibration standard solution.

3.1.17 Pipette 25 µL of the 80 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 2.0 ng/mL matrix matched calibration standard solution.

3.1.18 Pipette 25 µL of the 100 ng/mL spiking matrix matched calibration standard solution into an HPLC vial that contains 1 mL of buffered and filtered matrix. Mix solution well by inversion to create a 2.5 ng/mL matrix matched calibration standard solution.

## Method Amendment

### 3.2 Preparation of Test Solutions

3.2.6 Label the appropriate number of 15 mL plastic centrifuge tubes: 2 controls, at least 7 for target LOQ, and at least 5 for 10 x target LOQ

3.2.7 Shake the water sample and then add 10 mL of buffered water sample into each of the labeled centrifuge tubes.

3.2.9 To prepare at least 7 fortified recovery QC samples at the target LOQ (0.17 ng/mL), pipette 170  $\mu$ L of the 10 ng/mL mixed standard solution into each labeled centrifuge tube containing 10 mL of buffered water sample and mix well by inversion.

3.2.10 To prepare at least 5 fortified recovery QC samples at 10 x the target LOQ (1.7 ng/mL), pipette 170  $\mu$ L of the 100 ng/mL mixed standard solution into each labeled centrifuge tube containing 10 mL of buffered water sample and mix well by inversion.

3.2.11 Shake all samples at ~ 180 excursions per minute for 5 minutes.

### 3.5 Revisions

See Method Amendment, Modification Number 1

---