


## 2 Study Objective

The study objective was to perform an independent laboratory validation (ILV) on an analytical method for the determination of residues of Thiophanate-methyl and its metabolite Carbendazim in water. The primary method validation has been performed at PTRL Europe, Ulm, Germany and described in study no. P 2681 G according to guideline SANCO/825/00 rev. 8.1.

## 3 Material and Methods


### 3.1 Test and Reference Item(s)

Common Name:	<b>Thiophanate-methyl</b>
CAS-Registry-No.:	23564-05-8
Structural Formula:	
Molecular Formula:	$C_{12}H_{14}N_4O_4S_2$
Molecular Mass:	342.4 g/mol
Supplier:	Labor Dr. Ehrenstorfer, Germany
Batch No.:	30529
CIP (test facility) code:	07016
Purity:	99.0 %
Date of Certification:	25 October 2013
Date / Quantity of Receipt:	19 February 2014 / 250 mg
Appearance:	colourless solid
Recommended Storage Conditions:	Room temperature
Storage at CIP:	Refrigerated (nominally 4 °C)
Date of Expiry:	25 October 2017

Common Name: **Carbendazim**

CAS-Registry-No.: 10605-21-7

Structural Formula:



Molecular Formula:

$C_9H_9N_3O_2$

Molecular Mass:

191.2 g/mol

Supplier:

Labor Dr. Ehrenstorfer, Germany

Batch No.:

21129

CIP (test facility) code:

07024

Purity:

99.0 %

Date of Certification:

30 January 2013

Date / Quantity of Receipt:

05 March 2014 / 250 mg

Appearance:

colourless solid

Recommended

Storage Conditions:

Room temperature

Storage at CIP:

Refrigerated (nominally 4 °C)

Date of Expiry:

30 November 2018

All specifications of purity and composition of the reference item(s) were provided by the supplier. The certificate of analysis is depicted in Appendix 5.

## 3.2 Analytical Procedure

### 3.2.1 Specimen Origin

Drinking water sample was taken from the public water supply of the village of Schömburg, Kreis Calw, Germany.

The parameters of the used water sample are given in the following table:

Table 1: Typical parameters of the used water sample

	<b>Drinking water</b>
Origin	Schömburg, Kreis Calw, Germany
pH	8.03
conductivity ( $\mu\text{S}/\text{cm}$ )	170
TOC (mg/L)	0.56
SAK 254 nm (1/m)	1.32
Total hardness (mmol/L)*	0.87 (4.9 dGH)

\*One millimole of calcium (either  $\text{Ca}^{2+}$  or  $\text{CaCO}_3$ ) per litre of water corresponds to a hardness of 100.09 ppm or 5.608 dGH, since the molar mass of calcium carbonate is 100.09 g/mol.


All parameters were analyzed according to DIN EN ISO/IEC 17025:2005 (non-GLP).

### 3.2.2 Specimen Storage and Preparation

The water specimens were stored refrigerated ( $< 8\text{ }^\circ\text{C}$ ) in the dark. Before use, the sample was homogenized by shaking.

### 3.3 Outline of the Method

Water specimens are analysed for residues of Thiophanate-methyl and Carbendazim using the analytical method reported in study P 2681 G (JOOß, S., PTRL Europe, Ulm, Germany, 2012):



Representative chromatograms are presented in Appendix 4.

### 3.4 Equipment and Apparatus





- Adjustable pipettes
- Hamilton syringes
- HPLC autosampler vials, 1.8 mL
- HPLC with MS/MS detection (Dionex Ultimate 3000 with AB Sciex API 5500 QTRAP)

### 3.5 Reagents and Materials

- |                            |           |                  |
|----------------------------|-----------|------------------|
| • Acetonitrile, HPLC grade | Promochem | No. SO-9128-B025 |
| • Methanol, HPLC grade     | Promochem | No. SO-9260-B025 |
| • Formic acid, p.a.        | Merck     | No. 1.00264.1000 |
| • Acetic acid, p.a.        | Merck     | No. 1.00063.1000 |
| • Acetonitrile, LCMS grade | Promochem | No. SO-9340-B010 |
| • Water, LCMS grade        | Promochem | No. SO-9368-B010 |

### 3.5.1 Hazards or Precautions

For conduct of this method the German guidelines for laboratories „Working Safely in Laboratories, Basic Principles and Guidelines“ or comparable guidelines in other countries are to be observed. The following chemicals are used, which are classified by the hazardous material regulations according to Regulation (EC) No 1272/2008 [EU-GHS/CLP].

Acetonitrile:	
Methanol:	
Acetic acid	
Formic acid	

The pertinent safety instructions must be observed when working with all compound mentioned in this method (e.g. hazard (H) and precautionary (P) statements).

### 3.6 Preparation of Water Samples

Aliquots of the water samples were transferred to HPLC vials by means of adjustable pipettes. In case of recovery samples, 10  $\mu$ L of the respective standard solution, prepared in acetonitrile and 10  $\mu$ L 10% formic acid<sub>(aq)</sub> (final concentration 0.1 % formic acid) were added. The vials were capped tightly, shaken well and used for analysis by HPLC-MS/MS.

Time required for one sample set (12 samples including preparation of the analytical sequence and calibration standards) is about 2 hours for one person.

### 3.7 HPLC-MS/MS Conditions

**HPLC-MS/MS:** Dionex Ultimate 3000 with AB Sciex API 5500 QTRAP  
**Column:** Agilent Zorbax Eclipse XDB-C8 150 mm × 4.6 mm, 5.0 μm  
 (Part No. 993967-906)  
**Mobile phase:** A: Water+0.1% Formic acid  
 B: Acetonitrile +0.1% Formic acid

Time [min]	A [%]	B [%]	Flow [μL/min]	Gradient
0.00	80	20	600	linear
3.00	0	100	600	
6.00	0	100	600	
6.10	80	20	600	
10.00	80	20	600	

Column temp.:	35 °C
Interface:	ESI
Source polarity:	Positive
Ion Source:	Turbo Spray
Curtain gas	15 units
Collision gas:	Medium
Ion spray voltage:	5.5 kV
Interface temperature:	600 °C
GS1:	60 units
GS2:	0 units
Injection volume:	100 μL
Retention time:	Carbendazim ~ 3.6 min Thiophanate-methyl ~ 5.1 min
Split:	Approx. 1:5
Valve:	0-1 min to waste; 1-6 min to MS; 6-10 min to waste
Quantification:	Peak areas of the fragment ions, external standards



## Mass spectrometer parameters:

Analyte Monitored	Ions Monitored MRM	DP [V]	EP [V]	CE [V]	CXP [V]	Dwell Time [ms]
Carbendazim	192 → 160	148	10	15	16	200
	192 → 132	148	10	43	26	200
Thiophanate-methyl	343 → 151	142	10	26	48	300
	343 → 192	39	10	22	32	500
Ion Mode:		ESI (positive)				

Representative mass spectra are presented in Appendix 1.

The mass spectrometer was operated in MS/MS-Single Reaction Monitoring (MRM) positive ion mode. For Thiophanate-methyl, MRM 343 → 151 was suitable for quantification, MRM 343 → 192 was suitable as confirmation method.

For Carbendazim, MRM 192 → 160 was suitable for quantification, MRM 192 → 132 was suitable as confirmation method.

At the beginning of the sample sequence, a calibration curve was constructed by injecting standard solutions in the range from 0.01 µg/L to 1 µg/L (five points, single injection each). Injections of samples were interspersed with injections of standards to provide a continuous check on the instrument calibration.



### 3.8 Calculation of the Residues

External standard solutions comparable to the concentration expected in specimens were injected after each two or three samples.

The residues (R) in  $\mu\text{g/L}$  are calculated using external standards with one-point calibration according to the following equation:

$$R = \frac{A_x \cdot C_{St, \text{nominal}}}{A_{St}}$$

where

$A_x$ : Peak area of the analyte

$C_{St, \text{nominal}}$ : Nominal concentration of the analyte in external bracketing standard solutions, in  $\mu\text{g/L}$

$A_{St}$ : Peak area of the analyte in external standard solution (mean from bracketing standards)

The response factor (see Appendix 3) is calculated as follows:

$$\text{Response Factor} = \text{Peak Area [counts]} / \text{Concentration } [\mu\text{g/L}]$$

Recoveries were calculated by the following equation:

$$\text{Rec} = \frac{R_{\text{found}}}{R_{\text{fortified}}} \cdot 100 \%$$

Rec Recovery [%]

$R_{\text{found}}$  Analyte determined [ $\mu\text{g/L}$ ]

$R_{\text{fortified}}$  Fortification level [ $\mu\text{g/L}$ ]

Example calculation:

Recovery 0.5  $\mu\text{g/L}$  in drinking water (Thiophanate-methyl; quantifier):

Sample ID: RDW 0.5-1

$A_x$ : 60627

$C_{St, \text{nominal}}$ : 0.5  $\mu\text{g/L}$

$A_{St}$ : 63408 counts

$$R = \frac{60627 \text{ counts} \cdot 0.5 \mu\text{g/L}}{63408 \text{ counts}} = 0.478 \mu\text{g/L} \quad \text{Rec} = \frac{0.478 \mu\text{g/L}}{0.5 \mu\text{g/L}} \cdot 100 \% = 96\%$$



### 3.9 Standard Solutions

#### 3.9.1 Stock Solution

Stock solutions containing 1000 mg/L of Thiophanate-methyl (S1000) and 200 mg/L of Carbendazim (S200) were prepared in acetonitrile:methanol 7/3 (v/v) + 0.1% acetic acid as follows:

S1000: 20.2 mg Thiophanate-methyl (purity 99.0 %) in 20 mL

S200: 10.1 mg Carbendazim (purity 99.0 %) in 50 mL

Standards mixtures of Thiophanate-methyl and Carbendazim were prepared in acetonitrile and water containing 0.1% formic acid to achieve solutions in the concentration range needed for this study. Preparation of these working solutions is shown in the following table:

Table 2: Preparation of working solutions

Fortification Solution	Solution used for preparation	Volume used* [mL]	Final volume [mL]	Solvent	Concentration obtained [mg/L]
SM 50	S1000 Thiophanate-methyl S200 Carbendazim	1 5	20	Acetonitrile + 0.1% Formic acid	50 of each
SM 5	SM 50	1	10	Acetonitrile + 0.1% Formic acid	5 of each
SM 0.5	SM 5	1	10	Acetonitrile + 0.1% Formic acid	0.5 of each
SM 0.05	SM 0.5	1	10	Acetonitrile + 0.1% Formic acid	0.05 of each
SM 0.005	SM 0.05	1	10	Acetonitrile + 0.1% Formic acid	0.005 of each
SM 0.01	SM 0.5	0.2	10	Water + 0.1% Formic acid	0.01 of each

\* Volume used = Volume withdrawn from the solution indicated in the 2<sup>nd</sup> column.

All solutions in acetonitrile were stored deep frozen ( $\leq -18^{\circ}\text{C}$ ) in the dark. The solution in water was stored refrigerated ( $< 8^{\circ}\text{C}$ ) in the dark.

### 3.9.2 HPLC Standard Solutions

Chromatographic external calibration solutions were prepared by diluting working solutions in LCMS water using adjustable pipettes and Hamilton syringes. Preparation of calibration standards is described in Table 3.

A series of matrix-matched calibration standards in drinking water was prepared to determine the matrix effects of the water sample under investigation. Preparation of matrix-matched standards is shown in Table 4.

Table 3: Preparation of calibration standards in LCMS water

Standard Solution	Standard Solution used for preparation	Volume Used* [μL]	Solvent added [μL]	Solvent	10% formic acid <sub>(aq)</sub> added [μL]	Concentration obtained [μg/L]
Std 1 μg/L	SM 0.01	100	900	Water (LCMS grade)	10	1
Std 0.5 μg/L	SM 0.01	50	950		10	0.5
Std 0.1 μg/L	Std 1 μg/L	100	900		10	0.1
Std 0.05 μg/L	Std 0.5 μg/L	100	900		10	0.05
Std 0.01 μg/L	Std 0.1 μg/L	100	900		10	0.01

\* Volume used = Volume withdrawn from the solution indicated in the 2<sup>nd</sup> column.

Table 4: Preparation of matrix-matched standards in drinking water

Matrix-matched Standard Solution	Standard Solution used for preparation	Volume Used* [μL]	Solvent added [μL]	Solvent	10% formic acid <sub>(aq)</sub> added [μL]	Concentration obtained [μg/L]
MStd 1 μg/L DW	SM 0.01	100	900	Drinking Water	10	1
MStd 0.5 μg/L DW	SM 0.01	50	950		10	0.5
MStd 0.2 μg/L DW	MStd 1 μg/L DW	100	900		10	0.1
MStd 0.1 μg/L DW	MStd 0.5 μg/L DW	100	900		10	0.05
MStd 0.01 μg/L DW	MStd 0.1 μg/L DW	100	900		10	0.01

\* Volume used = Volume withdrawn from the solution indicated in the 2<sup>nd</sup> column.

## 4 Fortifications

Control (untreated) specimens of drinking water were fortified with the previously prepared working solutions as follows:

Table 5: Fortification of water samples

Matrix	Fortification level [µg/L]	Specimen volume [µL]	Fortification solution	Added [µL]	10% Formic acid <sub>(aq)</sub> added [µL]
Drinking water	0.05	1000	SM 0.005	10	10
	0.5	1000	SM 0.05	10	10

## 5 Deviations from the Study Plan

The study was performed according to the study plan dated 07 March 2014.  
 This report reflects the conduct of this study.