

4. DESCRIPTION OF THE ANALYTICAL METHOD

4.1. Reference Substances, Standard Solutions, and Fortification Solutions

The reference substances (analytical reference standards) used in this validation were IR9792/F9990, 3-hydroxy-IR9792/F9990, 1-carboxy-IR9792/F9990 (mixture of diastereomers), and pyrazole carboxamide. The details are described below:

| | |
|---------------------------|---------------------|
| Compound Name: | IR9792/F9990 |
| Reference Number (Lot #): | PL13-0195 |
| Percent purity: | 98.6% |
| CAS Number: | Not registered |
| FMC ID: | 510142 |
| Storage Conditions: | Ambient |
| Expiration date: | July 01, 2018 |
| Received date: | August 20, 2014 |

| | |
|---------------------------|-------------------------------|
| Compound Name: | 3-hydroxy-IR9792/F9990 |
| Reference Number (Lot #): | 30440/43 |
| Percent purity: | 97.2% |
| CAS Number: | Not registered |
| FMC ID: | 510152 |
| Storage Conditions: | Freezer |
| Expiration date: | April 30, 2017 |
| Received date: | April 27, 2015 |

| | |
|---------------------------|--|
| Compound Name: | 1-carboxy-IR9792/F9990 (mixture of diastereomers) |
| Reference Number (Lot #): | 55261-4-8 |
| Percent purity: | 97.8% |
| Diastereomeric ratio | 1.72:1 |
| CAS Number: | Not registered |
| FMC ID: | 510169 |
| Storage Conditions: | Freezer |
| Expiration date: | March 31, 2017 |
| Received date: | April 27, 2015 |

| | |
|---------------------------|-----------------------------|
| Compound Name: | Pyrazole carboxamide |
| Reference Number (Lot #): | R14151 |
| Percent purity: | 99.5% |
| CAS Number: | 925689-10-7 |
| FMC ID: | 510151 |
| Storage Conditions: | Freezer |
| Expiration date: | June 30, 2017 |
| Received date: | June 24, 2015 |

4.1.1. IR9792/F9990 Stock Standard Solution

A stock solution containing IR9792/F9990 was prepared by weighing approximately 0.010 g of IR9792/F9990 reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of IR9792/F9990.

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

4.1.1. 3-hydroxy-IR9792/F9990 Stock Standard Solution

A stock solution containing 3-hydroxy-IR9792/F9990 was prepared by weighing approximately 0.010 g of 3-hydroxy-IR9792/F9990 reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of 3-hydroxy-IR9792/F9990.

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

4.1.2. 1-carboxy-IR9792/F9990 Stock Standard Solution

A stock solution containing the 1-carboxy-IR9792/F9990 diastereomers was prepared by weighing approximately 0.010 g of 1-carboxy-IR9792/F9990 reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of total 1-carboxy-IR9792/F9990. Based on the diastereomeric ratio, the stock solution contained approximately 632 µg/mL and 368 µg/mL for diastereomers 1 and 2, respectively.

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

4.1.3. Pyrazole carboxamide Stock Standard Solution

A stock solution containing pyrazole carboxamide was prepared by weighing approximately 0.010 g of pyrazole carboxamide reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of pyrazole carboxamide.

For each compound, the following intermediate standards and fortification solutions were made individually. This was done so all four compounds could be used to spike the soil samples individually for a storage stability study, once the method was validated.

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

4.1.4. 10 µg/mL Intermediate Standard

An intermediate standard containing 10 µg/mL of compound was prepared by pipetting approximately 100 µL of the compound stock standard solution into a 10 mL class A volumetric flask and bringing it to volume with methanol.

The intermediate standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

4.1.5. 1 µg/mL Intermediate Standard / 10X LOQ Fortification Solution

An intermediate standard containing 1 µg/mL of compound was prepared by pipetting approximately 5 mL of the compound 10 µg/mL intermediate standard solution into a 50 mL class A volumetric flask and bringing it to volume with methanol.

The intermediate standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

4.1.6. 100 ng/mL Intermediate Standard / LOQ Fortification Solution

An intermediate standard containing 100 ng/mL of compound was prepared by pipetting approximately 1 mL of the 1 µg/mL intermediate standard solution into a 10 mL class A volumetric flask and bringing it to volume with methanol.

The intermediate standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

4.1.7. Working Standards (Calibration Standards)

Sets of working standards were prepared in volumetric flasks, and brought up to volume using methanol:water (8:2) (v:v).

The working standards were stored chilled at approximately 1-6 °C when not being used in the laboratory.

Example Preparation:

| Target Concentration (ng/mL) | Volume of Intermediate Standard (mL) | Concentration of Intermediate Standard | Final Volume (mL) |
|------------------------------|--------------------------------------|--|-------------------|
| 0.1 | 1.0 | 1 ng/mL Working Std | 10.0 |
| 1.0 | 1.0 | 10 ng/mL Working Std | 10.0 |
| 2.5 | 0.25 | 100 ng/mL Int Std | 10.0 |
| 5.0 | 0.5 | 100 ng/mL Int Std | 10.0 |
| 10.0 | 1.0 | 100 ng/mL Int Std | 10.0 |

4.2. Determination of IR9792/F9990, 3-hydroxy-IR9792/F9990, 1-carboxy-IR9792/F9990 (diastereomer 1), 1-carboxy-IR9792/F9990 (diastereomer 2), and pyrazole carboxamide in Soil by Ultra Performance Liquid Chromatography-Tandem Mass Spectrometry

4.2.1. Reagents and Solutions

Formic acid
 Ammonium Acetate
 Methanol, Pesticide Grade
 Methanol, HPLC Grade
 Acetone, pesticide grade

Hydrochloric acid
MilliQ Deionized (DI) Water, 18.2 Ω

Extraction solution 1: Acetone:water (9:1) (v:v). Using a graduated cylinder, add 900 mL of pesticide grade acetone to a 1L solution bottle. Next, using a graduated cylinder, add 100 mL of deionized water to the same 1L solution bottle. Shake well to mix. The solution was stored at room temperature for up to 30 days after preparation.

Extraction solution 2: Acetone:water (1:1) (v:v). Using a graduated cylinder, add 500 mL of pesticide grade acetone to a 1L solution bottle. Next, using a graduated cylinder, add 500 mL of deionized water to the same 1L solution bottle. Shake well to mix. The solution was stored at room temperature for up to 30 days after preparation.

Extraction solution 3: Acetone:0.5N HCl (1:1) (v:v). Using a graduated cylinder, add 500 mL of pesticide grade acetone to a 1L solution bottle. Next, using a graduated cylinder, add 500 mL of 0.5 N HCl to the same 1L solution bottle. Shake well to mix. The solution was stored at room temperature for up to 30 days after preparation. Note: 0.5N HCl was prepared by adding approximately 20.83 mL of concentrated HCl to 500 mL of deionized water.

Mobile Phase A: 10mM ammonium acetate and 0.2% formic acid in water. Using a graduated cylinder, 1000 mL of DI water was measured and transferred to a mobile phase bottle. 2.0 mL of formic acid and 0.7708g of ammonium acetate was added to the bottle, and mixed well. The solution was stored at room temperature for up to 30 days after preparation.

Mobile phase B: 0.2% formic acid in methanol. Using a graduated cylinder, 1000 mL of methanol was measured and transferred to a mobile phase bottle. 2.0 mL formic acid was added to the bottle and mixed well. The solution was stored at room temperature for up to 30 days after preparation.

4.2.2. Equipment and Instruments

Freezer(s), capable of maintaining $-20\text{ }^{\circ}\text{C} \pm 10\text{ }^{\circ}\text{C}$

Refrigerator(s), capable of maintaining $1\text{-}6\text{ }^{\circ}\text{C}$

Balance, Analytical, capable of weighing to the nearest 0.01 mg

Balance, Analytical, capable of weighing to the nearest 0.1 mg

Balance, Top Loading, capable of weighing to the nearest 0.1g

Sonicator

Vortex

Adjustable pipettes, 200 μL , 1000 μL , 5000 μL

Graduated cylinders

Wrist action shaker

Centrifuge capable of 4000 rpm

Autosampler vials and caps, 2-mL, screw cap

UHPLC column,

UHPLC system, Dionex Ultimate 3000

Mass Spectrometer, ABSciex Qtrap 4000

UHPLC system, Shimadzu Nexera XR

Mass Spectrometer, Sciex 6500+ QQQ

UPLC software, Dionex Chromeleon, bundled with Analyst version 1.6.2

MS/MS software, ABSciex Analyst, version 1.6.2
Data analysis software, ABSciex MultiQuant, version 2.0.2
Microsoft Office Excel 2007

4.2.3. Sample Extraction

- Weigh 5.0 ± 0.1 g of each sample into a 50 mL centrifuge tube.
- Using untreated sample, make fortifications at the LOQ and 10X LOQ.
- Add 250 μ L of the LOQ fortification and 250 μ L of the 10 X LOQ fortification solutions to their respective samples.
- Add 25 mL of acetone:water (9:1) to each centrifuge tube and place on a wrist action shaker for approximately 60 minutes.
- Centrifuge samples at 4000 rpm for approximately 10 minutes.
- Decant supernatant into a 100 mL volumetric flask. (A funnel with filter paper may be used).
- Repeat extraction by adding 25 mL of acetone:water (1:1) to the sample.
- Sonicate the sample at approximately 40 °C for 60 minutes.
- Shake on wrist-action shaker for approximately 60 minutes.
- Centrifuge samples at 4000 rpm for approximately 10 minutes.
- Decant supernatant and combine with first supernatant in the 100 mL volumetric flask.
- Repeat extraction by adding 25 mL of acetone:0.5N HCl (1:1) to the sample.
- Sonicate the sample at approximately 40 °C for 60 minutes.
- Shake on wrist-action shaker for approximately 60 minutes.
- Centrifuge samples at 4000 rpm for approximately 10 minutes.
- Decant supernatant and combine with the other two supernatants in the 100 mL volumetric flask.
- Bring up to the 100 mL volume mark with acetone. Shake well to mix.
- Pour the mixed extract into a 50 mL centrifuge tube and fill to the top. The rest of the sample can be discarded.
- Centrifuge the sample at 4000 rpm for approximately 10 minutes.
- Using a transfer pipette, transfer a portion of the supernatant to a 15 mL centrifuge tube, up to the 5 mL graduation mark.
- Place the sample under nitrogen evaporation in order to evaporate off the acetone. Once this is done, approximately 1 to 2 mL of water should be left.
- Remove from nitrogen evaporation, and bring back up to the 5 mL graduation mark with methanol.
- Shake well to mix.
- Transfer a portion to an HPLC vial for LC-MS/MS analysis.

4.2.4. Instrument Analysis

UPLC-MS/MS was used to quantify residues of IR9792/F9990 its metabolites, 3-hydroxy-IR9792/F9990, 1-carboxy-IR9792/F9990 diastereomers, and pyrazole carboxamide in soil. Two different sets of instruments were utilized over the course of the study. At first, a Dionex HPLC with a Sciex QTrap 4000 mass spectrometer, and later a Shimadzu HPLC with a Sciex 6500+ mass spectrometer. Both sets of parameters are below:

| | |
|---------------|----------------------------------|
| UHPLC system: | Dionex Ultimate 3000 |
| MS Detector: | Sciex Triple Quadrapole API 4000 |

UHPLC column: Phenomenex Kinetex 2.6 μ C18 100A, 50 x 4.6 mm
Guard column: SecurityGuard Ultra Cartridge UHPLC C18 for 4.6 mm
Column temp: 20°C
Injection volume: 10 μ L
Run Time: 24.0 minutes
Mobile phase A: 10mM ammonium acetate and 0.2% formic acid in water
Mobile phase B: 0.2% formic acid in methanol
Flow rate: 0.600 mL/minute
Approx. Retention Times: pyrazole carboxamide: 2.7 min (positive polarity)
1-carboxy-IR9792/F9990 Diastereomer 1: 11.80 min (positive polarity)
1-carboxy-IR9792/F9990 Diastereomer 2: 12.63 min (positive polarity)
3-hydroxy-IR9792/F9990: 17.63 min (negative polarity)
IR9792/F9990: 18.78 min (positive polarity)

Gradient program:

| Time (min) | %A | % B |
|------------|----|-----|
| 0.0 | 90 | 10 |
| 3.5 | 90 | 10 |
| 6.0 | 60 | 40 |
| 22.0 | 20 | 80 |
| 23.0 | 20 | 80 |
| 23.1 | 90 | 10 |
| 24.0 | 90 | 10 |

Curtain Gas: 10.00
Collision Gas: High
IonSpray Voltage: 5500, -4500
Temperature: 450 °C
Ion Source Gas 1: 60
Ion Source Gas 2: 60

| Compound | Parent Ion (da) | Daughter Ion (da) | Dwell time (msec) | Declustering Potential (V) | Entrance Potential (V) | Collision Energy (V) | Collision Cell Exit Potential (V) | Polarity |
|--|-----------------|-------------------|-------------------|----------------------------|------------------------|----------------------|-----------------------------------|----------|
| Pyrazole carboxamide | 176.059 | 136 | 500 | 41 | 10 | 23 | 8 | Positive |
| Pyrazole carboxamide* | 176.059 | 156 | 500 | 41 | 10 | 13 | 10 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 1) | 382.1 | 336.200 | 500 | 76 | 10 | 21 | 10 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 1)* | 382.1 | 296.100 | 500 | 76 | 10 | 35 | 20 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 2) | 382.1 | 336.201** | 500 | 76 | 10 | 21 | 10 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 2)* | 382.1 | 296.101** | 500 | 76 | 10 | 35 | 20 | Positive |
| 3-hydroxy-IR9792/F9990 | 366.235 | 174.8 | 1000 | -110 | -10 | -30 | -11 | Negative |
| 3-hydroxy-IR9792/F9990* | 366.235 | 130.7 | 1000 | -110 | -10 | -40 | -7 | Negative |
| IR9792/F9990 | 352.185 | 332.1 | 1000 | 106 | 10 | 23 | 8 | Positive |
| IR9792/F9990* | 352.185 | 256.1 | 1000 | 106 | 10 | 41 | 16 | Positive |

*Confirmatory ion transition.

**For the diastereomers, AB Sciex Analyst program required that the daughter ions not have the exact same molecular weight.

| | |
|-------------------|--|
| UHPLC system: | Shimadzu Nexera XR HPLC |
| MS Detector: | Sciex Triple Quadrupole 6500+ |
| UHPLC column: | Phenomenex Kinetex 2.6 μ C18 100A, 50 x 4.6 mm |
| Guard column: | SecurityGuard Ultra Cartridge UHPLC C18 for 4.6 mm |
| Column temp: | 20°C |
| Injection volume: | 10 μ L |
| Run Time: | 24.0 minutes |

Mobile phase A: 10mM ammonium acetate and 0.2% formic acid in water
Mobile phase B: 0.2% formic acid in methanol
Flow rate: 0.600 mL/minute
Approx. Retention Times: pyrazole carboxamide: 2.8 min (positive polarity)
1-carboxy-IR9792/F9990 Diastereomer 1: 11.8 min (positive polarity)
1-carboxy-IR9792/F9990 Diastereomer 2: 12.6 min (positive polarity)
3-hydroxy-IR9792/F9990: 17.8 min (negative polarity)
IR9792/F9990: 18.9 min (positive polarity)

Gradient program:

| Time (min) | %A | % B |
|------------|----|-----|
| 0.0 | 90 | 10 |
| 3.5 | 90 | 10 |
| 6.0 | 60 | 40 |
| 22.0 | 20 | 80 |
| 23.0 | 20 | 80 |
| 23.1 | 90 | 10 |
| 24.0 | 90 | 10 |

Curtain Gas: 20.00
Collision Gas: 7.0
IonSpray Voltage: 5500, -4500
Temperature: 450 °C
Ion Source Gas 1: 75.0
Ion Source Gas 2: 75.0

| Compound | Parent Ion (da) | Daughter Ion (da) | Dwell time (msec) | Declustering Potential (V) | Entrance Potential (V) | Collision Energy (V) | Collision Cell Exit Potential (V) | Polarity |
|--|-----------------|-------------------|-------------------|----------------------------|------------------------|----------------------|-----------------------------------|----------|
| Pyrazole carboxamide | 176.0 | 136.0 | 350 | 26 | 10 | 21 | 16 | Positive |
| Pyrazole carboxamide* | 176.0 | 156.0 | 350 | 26 | 10 | 13 | 18 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 1) | 382.1 | 336.1 | 350 | 101 | 10 | 19 | 20 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 1)* | 382.1 | 296.1 | 350 | 101 | 10 | 33 | 20 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 2) | 382.1 | 336.1 | 350 | 101 | 10 | 19 | 20 | Positive |
| 1-carboxy-IR9792/F9990 (diastereomer 2)* | 382.1 | 296.1 | 350 | 101 | 10 | 33 | 20 | Positive |
| 3-hydroxy-IR9792/F9990 | 366.1 | 175.0 | 1000 | -130 | -10 | -28 | -11 | Negative |
| 3-hydroxy-IR9792/F9990 * | 366.1 | 131.0 | 1000 | -130 | -10 | -36 | -13 | Negative |
| IR9792/F9990 | 352.1 | 256.1 | 1000 | 86 | 10 | 39 | 16 | Positive |
| IR9792/F9990 * | 352.1 | 312.1 | 1000 | 86 | 10 | 25 | 22 | Positive |

*Confirmatory ion transition.

4.5. Calculations

4.5.1. Linear Regression (IR9792/F9990, 3-hydroxy-IR9792/F9990, 1-carboxy-IR9792/F9990 diastereomers, and pyrazole carboxamide).

The following formulas were used to calculate IR9792/F9990, 3-hydroxy-IR9792/F9990, 1-carboxy-IR9792/F9990 diastereomers, and pyrazole carboxamide fortification recoveries. The amount found (ppb) of IR9792/F9990, 3-hydroxy-IR9792/F9990, 1-carboxy-IR9792/F9990 diastereomers, and pyrazole carboxamide residue was calculated from the linear regression equation with 1/x weighting.

The linear calibration equation is:

$$y = mx + b$$

Where:

y: area count from LC-MS/MS chromatogram
m: slope of the linear regression line
b: y-intercept of the linear regression line
x: sample concentration in ng/mL

4.5.2. Amount Found (ppb)

Below is an example of how the amount of residue found is determined for IR9792/F9990. This same equation is used for all compounds in this validation. The amount found (ppb) IR9792/F9990 was calculated by solving the linear equation for x, concentration, using the detector response (area) substituted in for y, and multiplying the result by the dilution factor.

$$x = \frac{y-b}{m} * \text{dilution factor}$$

$$\text{Regression line: } y = 2553330 (x) + 18546.41437$$

| | |
|-------------------|------------------|
| Example: | IR9792/F9990 |
| Lab Sample ID: | T1-SC-4-SPA-6-R1 |
| Sample Peak Area: | 27940977 |
| Dilution factor: | 20 |

$$x = \frac{27940977 - (-18546.41437)}{2553330} * 20 = 219.004 \text{ ppb}$$

Analyst MultiQuant software (ABSciex, Version 2.0.2) was used to calculate linear, 1/x weighted regression calibration curves and amount found (ppb).

4.5.3. Dilution Factor

The dilution factor for IR9792/F9990, 3-hydroxy-IR9792/F9990, 1-carboxy-IR9792/F9990 (diastereomer 1), 1-carboxy-IR9792/F9990 (diastereomer 2), and pyrazole carboxamide samples are calculated with the equation below:

$$\text{dilution factor} = \frac{\text{total dilution volume (mL)}}{\text{sample weight (g)}}$$

Example: IR9792/F9990
 Lab Sample ID: LOQ
 Sample Weight: 5.0 g
 Total Dilution Vol: 100 mL

$$\text{dilution factor} = \frac{100 \text{ mL}}{5.0 \text{ g}} = 20$$

4.5.4. Fortification Level IR9792/F9990 (ppb added)

Below are examples of the fortifications used for the validation.

$$\text{Fortification level (ppm)} = \frac{\text{Vol. Spiking Soln. (mL)} \times \text{Spiking Soln. Conc. } \left(\frac{\mu\text{g}}{\text{mL}}\right)}{\text{Sample weight (g)}}$$

Example: IR9792/F9990
 Lab Sample ID: Interval 4 LOQ 1
 Volume Spiking Soln: 0.250 mL
 Identity of Spiking Soln: LOQ Fortification Solution
 Spiking Soln. Conc.: 0.107474 $\mu\text{g/mL}$
 Sample Weight: 5.0 g

$$\text{Fortification level} = 0.250 \text{ mL} * 0.107474 \mu\text{g/mL} / 5.0 \text{ g} = 0.0053737 \text{ ppm (5.3737 ppb)}$$

Example: IR9792/F9990
 Lab Sample ID: Interval 4 10X LOQ 1
 Volume Spiking Soln: 0.250 mL
 Identity of Spiking Soln: 10X LOQ Fortification Solution
 Spiking Soln. Conc.: 1.07474 $\mu\text{g/mL}$
 Sample Weight: 5.0 g

$$\text{Fortification level} = 0.250 \text{ mL} * 1.07474 \mu\text{g/mL} / 5.0 \text{ g} = 0.053737 \text{ ppm (53.737ppb)}$$

Summary Title

Method Validation of F9990, 3-hydroxy-F9990, F9990-1-carboxylate and
3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide in soil

Author

Christopher Skaggs

Sponsor

FMC Corporation
Agricultural Solutions
701 Princeton South Corporate Center
Ewing, NJ 08628
United States of America

Performing Laboratory

SGS North America, Inc.
GLP Laboratory
241 34th Avenue
Brookings, SD 57006
United States of America

Total Number of Pages

50

INTRODUCTION

The following is a summary of the steps taken to validate an analytical method for F9990 and its soil degradates, 3-hydroxy-F9990, F9990-1-carboxylate diastereomers, and 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide in soil. Two types of soil were used, a sandy loam from Nebraska, USA, and a silt loam from New York, USA. This validation is to be used as part of the larger ongoing studies being performed by SGS for FMC. A description of the analytical methods, along with instrumental parameters and raw data from the datasets are found below.

DESCRIPTION OF THE ANALYTICAL METHOD

Reference Substances, Standard Solutions, and Fortification Solutions

The reference substances (analytical reference standards) used in this validation were F9990, 3-hydroxy-F9990, F9990-1-carboxylate, and 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide. The details are described below:

| | |
|---------------------------|-----------------|
| Compound Name: | F9990 |
| Reference Number (Lot #): | PL13-0195 |
| Percent purity: | 98.6% |
| CAS Number: | Not registered |
| Storage Conditions: | Ambient |
| Expiration date: | July 01, 2018 |
| Received date: | August 20, 2014 |

| | |
|---------------------------|------------------------|
| Compound Name: | 3-hydroxy-F9990 |
| Reference Number (Lot #): | 30440/43 |
| Percent purity: | 97.2% |
| CAS Number: | Not registered |
| Storage Conditions: | Freezer |
| Expiration date: | April 30, 2017 |
| Received date: | April 27, 2015 |

| | |
|---------------------------|---|
| Compound Name: | F9990-1-carboxylate (mixture of diastereomers) |
| Reference Number (Lot #): | 55261-4-8 |
| Percent purity: | 97.8% |
| Diastereomeric ratio | 1.72:1 |
| CAS Number: | Not registered |
| Storage Conditions: | Freezer |
| Expiration date: | March 31, 2017 |
| Received date: | April 27, 2015 |

| | |
|---------------------------|--|
| Compound Name: | 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide |
| Reference Number (Lot #): | R14151 |
| Percent purity: | 99.5% |
| CAS Number: | 925689-10-7 |
| Storage Conditions: | Freezer |
| Expiration date: | June 30, 2017 |
| Received date: | June 24, 2015 |

F9990 Stock Standard Solution

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

A stock solution containing F9990 was prepared by weighing approximately 0.010 g of F9990 reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of F9990.

3-hydroxy-F9990 Stock Standard Solution

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

A stock solution containing 3-hydroxy-F9990 was prepared by weighing approximately 0.010 g of 3-hydroxy-F9990 reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of 3-hydroxy-F9990.

F9990-1-carboxylate Stock Standard Solution

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

A stock solution containing the F9990-1-carboxylate diastereomers was prepared by weighing approximately 0.010 g of F9990-1-carboxylate reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of F9990-1-carboxylate.

3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide Stock Standard Solution

The stock standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

A stock solution containing 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide was prepared by weighing approximately 0.010 g of 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide reference substance into a 10 mL volumetric flask. The flask was diluted to volume with methanol, capped, and mixed thoroughly. The stock solution contained approximately 1000 µg/mL of 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide.

For each compound, the following intermediate standards and fortification solutions were made individually. This was done so all four compounds could be used to spike the soil samples individually for a storage stability study, once the method was validated.

10 µg/mL Intermediate Standard

The intermediate standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

An intermediate standard containing 10 µg/mL of compound was prepared by pipetting approximately 100 µL of the compound stock standard solution into a 10 mL class A volumetric flask and bringing it to volume with methanol.

1 µg/mL Intermediate Standard / 10X LOQ Fortification Solution

The intermediate standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

An intermediate standard containing 1 µg/mL of compound was prepared by pipetting approximately 5 mL of the compound 10 µg/mL intermediate standard solution into a 50 mL class A volumetric flask and bringing it to volume with methanol.

100 ng/mL Intermediate Standard / LOQ Fortification Solution

The intermediate standard was stored chilled at approximately 1-6 °C when not being used in the laboratory.

An intermediate standard containing 100 ng/mL of compound was prepared by pipetting approximately 1 mL of the 1 µg/mL intermediate standard solution into a 10 mL class A volumetric flask and bringing it to volume with methanol.

Working Standards (Calibration Standards)

The working standards were stored chilled at approximately 1-6 °C when not being used in the laboratory.

Sets of working standards were prepared in volumetric flasks, and brought up to volume using methanol:water (8:2) (v:v).

Example Preparation:

| Target Concentration (ng/mL) | Volume of Intermediate Standard (mL) | Concentration of Intermediate Standard | Final Volume (mL) |
|------------------------------|--------------------------------------|--|-------------------|
| 0.1 | 1.0 | 1 ng/mL Working Std | 10.0 |
| 1.0 | 1.0 | 10 ng/mL Working Std | 10.0 |
| 2.5 | 0.25 | 100 ng/mL Int Std | 10.0 |
| 5.0 | 0.5 | 100 ng/mL Int Std | 10.0 |
| 10.0 | 1.0 | 100 ng/mL Int Std | 10.0 |

Determination of F9990, 3-hydroxy-F9990, F9990-1-carboxylate diastereomers and 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide in Soil by Ultra Performance Liquid Chromatography-Tandem Mass Spectrometry

Reagents and Solutions

Formic acid
Ammonium Acetate
Methanol, Pesticide Grade
Methanol, HPLC Grade
Acetone, pesticide grade
Hydrochloric acid
MilliQ Deionized (DI) Water, 18.2 Ω

Extraction solution 1: Acetone:water (9:1) (v:v). Using a graduated cylinder, add 900 mL of pesticide grade acetone to a 1L solution bottle. Next, using a graduated cylinder, add 100 mL of deionized water to the same 1L solution bottle. Shake well to mix. The solution was stored at room temperature for up to 30 days after preparation.

Extraction solution 2: Acetone:water (5:5) (v:v). Using a graduated cylinder, add 500 mL of pesticide grade acetone to a 1L solution bottle. Next, using a graduated cylinder, add 500 mL of deionized water to the same 1L solution bottle. Shake well to mix. The solution was stored at room temperature for up to 30 days after preparation.

Extraction solution 3: Acetone:0.5N HCl (5:5) (v:v). Using a graduated cylinder, add 500 mL of pesticide grade acetone to a 1L solution bottle. Next, using a graduated cylinder, add 500 mL of 0.5 N HCl to the same 1L solution bottle. Shake well to mix. The solution was stored at room temperature for up to 30 days after preparation. Note: 0.5N HCl was prepared by adding approximately 20.83 mL of concentrated HCl to 500 mL of deionized water.

Mobile Phase A: 10mM ammonium acetate and 0.2% formic acid in water. Using a graduated cylinder, 1000 mL of DI water was measured and transferred to a mobile phase bottle. 2.0 mL of formic acid and 0.7708g of ammonium acetate was added to the bottle, and mixed well. The solution was stored at room temperature for up to 30 days after preparation.

Mobile phase B: 0.2% formic acid in methanol. Using a graduated cylinder, 1000 mL of methanol was measured and transferred to a mobile phase bottle. 2.0 mL formic acid was added to the bottle and mixed well. The solution was stored at room temperature for up to 30 days after preparation.

Equipment and Instruments

Freezer(s), capable of maintaining $-20\text{ }^{\circ}\text{C} \pm 10\text{ }^{\circ}\text{C}$
Refrigerator(s), capable of maintaining $1\text{-}6\text{ }^{\circ}\text{C}$
Balance, Analytical, capable of weighing to the nearest 0.01 mg
Balance, Analytical, capable of weighing to the nearest 0.1 mg
Balance, Top Loading, capable of weighing to the nearest 0.1g
Sonicator
Vortex
Adjustable pipettes, 200 μL , 1000 μL , 5000 μL

Graduated cylinders
Wrist action shaker
Centrifuge capable of 4000 rpm
Autosampler vials and caps, 2-mL, screw cap
UHPLC column,
UHPLC system, Dionex Ultimate 3000
Mass Spectrometer, ABSciex Qtrap 4000
UPLC software, Dionex Chromeleon, bundled with Analyst version 1.5.2
MS/MS software, ABSciex Analyst, version 1.5.2
Data analysis software, ABSciex MultiQuant, version 2.0.2
Microsoft Office Excel 2007

Sample Extraction of F9990, 3-hydroxy-F9990, F9990-1-carboxylate, and 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide

- Weigh 5.0 ± 0.1 g of each sample into a 50 mL centrifuge tube.
- Using untreated samples, make fortifications at the LOQ and 10X LOQ.
- Add 250 μ L of the LOQ fortification and 250 μ L of the 10 X LOQ fortification solutions to their respective samples.
- Add 25 mL of acetone:water (9:1) to each centrifuge tube and place on a wrist action shaker for approximately 60 minutes.
- Centrifuge samples at 4000 rpm for approximately 10 minutes.
- Decant supernatant into a 100 mL volumetric flask. (A funnel with filter paper may be used).
- Repeat extraction by adding 25 mL of acetone:water (5:5) to the sample.
- Sonicate the sample at approximately 40 °C for 60 minutes.
- Shake on wrist-action shaker for approximately 60 minutes.
- Centrifuge samples at 4000 rpm for approximately 10 minutes.
- Decant supernatant and combine with first supernatant in the 100 mL volumetric flask.
- Repeat extraction by adding 25 mL of acetone:0.5N HCl (5:5) to the sample.
- Sonicate the sample at approximately 40 °C for 60 minutes.
- Shake on wrist-action shaker for approximately 60 minutes.
- Centrifuge samples at 4000 rpm for approximately 10 minutes.
- Decant supernatant and combine with the other two supernatants in the 100 mL volumetric flask.
- Bring up to the 100 mL volume mark with acetone. Shake well to mix.
- Pour the mixed extract into a 50 mL centrifuge tube and fill to the top. The rest of the sample can be discarded.
- Centrifuge the sample at 4000 rpm for approximately 10 minutes.
- Using a transfer pipette, transfer a portion of the supernatant to a 15 mL centrifuge tube, up to the 5 mL graduation mark.
- Place the sample under nitrogen evaporation in order to evaporate off the acetone. Once this is done, approximately 1.5 mL of water should be left.
- Remove from nitrogen evaporation, and bring back up to the 5 mL graduation mark with methanol.
- Shake well to mix
- Transfer a portion to an HPLC vial for LC-MS/MS analysis.

Instrument Analysis

UPLC-MS/MS was used to quantify residues of F9990 its metabolites, 3-hydroxy-F9990, F9990-1-carboxylate diastereomers, and 3-(difluoromethyl)-1H-pyrazole-4-carboxamide in soil.

UHPLC system: Dionex Ultimate 3000
 MS Detector: Triple Quadrapole API 4000
 UHPLC column: Phenomenex Kinetex 2.6 μ C18 100A 5 μ m, 50 x 4.6 mm
 Guard column: SecurityGuard Ultra Cartridge UHPLC C18 for 4.6 mm
 Column temp: 20°C
 Injection volume: 10 μ L
 Run Time: 24.0 minutes
 Mobile phase A: 10mM ammonium acetate and 0.2% formic acid in water
 Mobile phase B: 0.2% formic acid in methanol
 Flow rate: 0.600 mL/minute
 Approx. Retention Times: 3-(difluoromethyl): 2.7 min (positive polarity)
 F9990-1-carboxylate Diastereomer 1: 11.80 min (positive polarity)
 F9990-1-carboxylate Diastereomer 2: 12.63 min (positive polarity)
 3-hydroxy-F9990: 17.63 min (negative polarity)
 F9990: 18.78 min (positive polarity)

Gradient program:

| Time (min) | %A | % B |
|------------|----|-----|
| 0.0 | 90 | 10 |
| 3.5 | 90 | 10 |
| 6.0 | 60 | 40 |
| 22.0 | 20 | 80 |
| 23.0 | 20 | 80 |
| 23.1 | 90 | 10 |
| 24.0 | 90 | 10 |

Curtain Gas: 10.00
 Collision Gas: High
 IonSpray Voltage: 5500, -4500
 Temperature: 450 °C
 Ion Source Gas 1: 60
 Ion Source Gas 2: 60

| Compound | Parent Ion (da) | Daughter Ion (da) | Dwell time (msec) | Declustering Potential (V) | Entrance Potential (V) | Collision Energy (V) | Collision Cell Exit Potential (V) | Polarity |
|--|-----------------|-------------------|-------------------|----------------------------|------------------------|----------------------|-----------------------------------|----------|
| 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide | 176.059 | 136 | 500 | 41 | 10 | 23 | 8 | Positive |
| 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide* | 176.059 | 156 | 500 | 41 | 10 | 13 | 10 | Positive |
| F9990-1-carboxylate (diastereomer 1) | 382.1 | 336.200 | 500 | 76 | 10 | 21 | 10 | Positive |
| F9990-1-carboxylate (diastereomer 1)* | 382.1 | 296.100 | 500 | 76 | 10 | 35 | 20 | Positive |
| F9990-1-carboxylate (diastereomer 2) | 382.1 | 336.201** | 500 | 76 | 10 | 21 | 10 | Positive |
| F9990-1-carboxylate (diastereomer 2)* | 382.1 | 296.101** | 500 | 76 | 10 | 35 | 20 | Positive |
| 3-hydroxy-F9990 | 366.235 | 174.8 | 1000 | -110 | -10 | -30 | -11 | Negative |
| 3-hydroxy-F9990* | 366.235 | 130.7 | 1000 | -110 | -10 | -40 | -7 | Negative |
| F9990 | 352.185 | 332.1 | 1000 | 106 | 10 | 23 | 8 | Positive |
| F9990* | 352.185 | 256.1 | 1000 | 106 | 10 | 41 | 16 | Positive |

*Confirmatory ion transition.

**For the diastereomers, AB Sciex Analyst program required that the daughter ions not have the exact same molecular weight.

Calculations

Linear Regression (F9990, 3-hydroxy-F9990, F9990-1-carboxylate diastereomers, and 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide)

The following formulas were used to calculate F9990, 3-hydroxy-F9990, F9990-1-carboxylate diastereomers, and 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide fortification recoveries.

The amount found (ppb) of F9990, 3-hydroxy-F9990, F9990-1-carboxylate diastereomers, and 3-(difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide residue was calculated from the linear regression equation with 1/x weighting.

The linear calibration equation is:

$$y = mx + b$$

The slope (m) and intercept (b) are calculated as:

$$m = (\sum w \sum wxy - \sum wx \sum wy) / D_x$$

$$b = (\sum wx^2 \sum wy - \sum wx \sum wxy) / D_x$$

Where:

$$D_x = \sum w \sum wx^2 - (\sum wx)^2$$

$$D_y = \sum w \sum wy^2 - (\sum wy)^2$$

Amount Found (ppb)

Below is an example of how the amount of residue found is determined for F9990. This same equation is used for all compounds in this validation. The amount found (ppb) F9990 was calculated by solving the linear equation for x, concentration, using the detector response (area) substituted in for y, and multiplying the result by the dilution factor.

$$x = \frac{y-b}{m} * \text{dilution factor}$$

| | |
|-------------------|----------------|
| Example: | F9990 |
| Lab Sample ID: | Nebraska LOQ 1 |
| Sample Peak Area: | 8240.95 |
| Dilution factor: | 20 |

$$x = \frac{8240.95 - (1018.48685)}{32234.90} * 20 = 4.4811 \text{ ppb}$$

Analyst MultiQuant software (ABSciex, Version 2.0.2) was used to calculate linear, 1/x weighted regression calibration curves and amount found (ppb).

Dilution Factor

The dilution factor for F9990, 3-hydroxy-F9990, F9990-1-carboxylate diastereomers, and 3-difluoromethyl)-1-methyl-1H-pyrazole-4-carboxamide samples are calculated with the equation below:

$$\text{dilution factor} = \frac{\text{total dilution volume (mL)}}{\text{sample weight (g)}}$$

Example: F9990
 Lab Sample ID: Nebraska LOQ
 Sample Weight: 5.0 g
 Total Dilution Vol: 100 mL

$$\text{dilution factor} = \frac{100 \text{ mL}}{5.0 \text{ g}} = 20$$

Fortification Level F9990 (ppb added)

Below are examples of the fortifications used for the validation.

$$\text{Fortification level (ppm)} = \frac{\text{Vol. Spiking Soln. (mL)} \times \text{Spiking Soln. Conc.} \left(\frac{\mu\text{g}}{\text{mL}}\right)}{\text{Sample weight (g)}}$$

Example: F9990
 Lab Sample ID: Nebraska LOQ 1
 Volume Spiking Soln: 0.250 mL
 Identity of Spiking Soln: LOQ Fortification Solution
 Spiking Soln. Conc.: 0.105009 $\mu\text{g/mL}$
 Sample Weight: 5.0 g

$$\text{Fortification level} = 0.250 \text{ mL} * 0.105009 \mu\text{g/mL} / 5.0 \text{ g} = 0.00525045 \text{ ppm (5.25045 ppb)}$$

Example: F9990
 Lab Sample ID: Nebraska 10X LOQ 1
 Volume Spiking Soln: 0.250 mL
 Identity of Spiking Soln: 10X LOQ Fortification Solution
 Spiking Soln. Conc.: 1.05009 $\mu\text{g/mL}$
 Sample Weight: 5.0 g

$$\text{Fortification level} = 0.250 \text{ mL} * 1.05009 \mu\text{g/mL} / 5.0 \text{ g} = 0.0525045 \text{ ppm (52.5045 ppb)}$$

