

## 1.0 INTRODUCTION

The purpose of this study was to validate an analytical method used to determine the content of dicloran in aqueous samples by liquid chromatography with tandem mass spectrometry detection (LC-MS/MS). The method was validated (9 May to 12 July 2019) to quantify the concentrations of dicloran present in recovery samples prepared in ground water and surface water. The analytical method was validated with regards to accuracy, precision, specificity, linearity, limit of quantitation (LOQ), limit of detection (LOD), method detection limit (MDL), and confirmation of analyte identification.

The method was validated in ground water and surface water by fortification with dicloran at concentrations of 0.100 (LOQ) and 1.00 (10X LOQ)  $\mu\text{g/L}$ . The samples were extracted twice with dichloromethane following acidification to pH 2 with phosphoric acid. The recovery samples were subsequently concentrated, then reconstituted and diluted with 20/80 acetonitrile/purified reagent water (v/v). The High-level (10X LOQ) recovery samples were further diluted into the calibration range with 20/80 acetonitrile/purified reagent water (v/v) and/or matrix blank final extract (see [Section 2.11](#)). All samples were analyzed using liquid chromatography with tandem mass spectrometry detection (LC-MS/MS).

The study was initiated on 27 April 2019, the day the Study Director signed the protocol, and was completed on the day the Study Director signed the final report. The experimental portion of the validation was conducted from 9 May to 12 July 2019 at Smithers Viscient, located in Wareham, Massachusetts. All original raw data, the protocol, and the final report produced during this study are stored in Smithers Viscient's archives at the above location.

## 2.0 MATERIALS AND METHODS

### 2.1 Protocol

Procedures used in this study followed those described in the Smithers Viscient protocol entitled “Environmental Chemistry Method: Validation of the Analytical Method for the Determination of Dicloran in Ground Water and Surface Water by LC-MS/MS” ([Appendix 1](#)). The study was conducted under Good Laboratory Practice (GLP) regulations and principles as described in 40 CFR 160 ([U.S. EPA, 1989](#)) and the OECD principles on GLP ([OECD, 1998](#)), and followed the SANCO/3029/99 rev 4 guidance document ([EC, 2000](#)) and OCSPP 850.6100 guideline ([U.S. EPA, 2012](#)).

### 2.2 Test and Reference Substances

#### 2.2.1 Test Substance

The test substance, dicloran technical, was received on 2 April 2019 from EPL Archives, Inc., Sterling, Virginia. The following information was provided:

Name:	Dicloran technical
Synonym(s):	BOTRAN technical; 2,6-dichloro-4-nitroaniline
Lot No.:	20130605
CAS No.:	99-30-9
Purity:	98.9% (Certificate of Analysis, <a href="#">Appendix 2</a> )
Recertification Date:	19 April 2021

Upon receipt at Smithers Viscient, dicloran technical (SMV No. 9932) was stored at room temperature in a dark, ventilated cabinet in the original container. Concentrations were adjusted for the purity of dicloran technical. This sample of dicloran technical was used to prepare recovery samples during testing.

### 2.2.2 Reference Substance

The reference substance, dichloran PESTANAL, was received on 11 April 2019 from EPL Archives, Inc., Sterling, Virginia. The following information was provided:

Name:	Dichloran PESTANAL
Synonym:	Dicloran
Batch No.:	SZBF103XV
CAS No.:	99-30-9
Purity:	99.6% (Certificate of Analysis, <a href="#">Appendix 2</a> )
Expiry Date:	13 April 2020

Upon receipt at Smithers Viscient, dichloran PESTANAL (SMV No. 9945) was stored at room temperature in a dark, ventilated cabinet in the original container. Concentrations were adjusted for the purity of dichloran PESTANAL. This sample of dichloran PESTANAL was used to prepare calibration standards during testing.

Determination of stability and characterization, verification of dicloran technical and dichloran PESTANAL identities, maintenance of records on dicloran technical and dichloran PESTANAL, and archival of a sample of dicloran technical and dichloran PESTANAL are the responsibility of the Study Sponsor.

### 2.3 Reagents

1. Acetonitrile: EMD, reagent grade
2. Dichloromethane: EMD, reagent grade
3. Phosphoric acid: EMD, reagent grade
4. Acetone: EMD, reagent grade
5. Methanol: EMD, reagent grade
6. Formic Acid: Honeywell, reagent grade
7. 0.1% formic acid in water: Fisher, reagent grade
8. 0.1% formic acid in acetonitrile: Surface Water Validation  
Honeywell, reagent grade; B & J, reagent grade  
  
Ground Water Validation  
Fisher, reagent grade

9. Purified reagent water: Prepared from a Millipore MilliQ Direct 8 water purification system (meets ASTM Type II requirements)

## 2.4 Instrumentation and Laboratory Equipment

1. Instrument: Ground Water Validation  
MDS Sciex API 5000 mass spectrometer equipped with an ESI Turbo V ion source  
Shimadzu SIL-20AHT autoinjector  
Shimadzu DGU-20A5R vacuum degassers  
Shimadzu LC-20ADXR solvent delivery pumps  
Shimadzu CTO-20AC column compartment  
Shimadzu CBM-20A communications bus  
Analyst 1.6.3 software for data acquisition
- Surface Water Validation  
MDS Sciex API 5000 mass spectrometer equipped with an ESI Turbo V ion source  
Shimadzu SIL-20AHT autoinjector  
Shimadzu DGU-20A5R vacuum degassers  
Shimadzu LC-20AD solvent delivery pumps  
Shimadzu CTO-20AC column compartment  
Shimadzu CBM-20A communications bus  
Analyst 1.6.3 software for data acquisition
2. Balance: Mettler Toledo XSE205DU
3. Centrifuge: Thermo Scientific Sorvall Legend XFR
4. Laboratory equipment: Positive displacement pipets, graduated cylinders, volumetric flasks, disposable glass pipets, stir bars, stir plate, vortex mixer, separatory funnels, round bottom flasks, pH paper, centrifuge tubes, rotary evaporator, amber bottles, clear vials with snap caps, amber vials with crimp caps, autosampler vials, and amber glass bottles with Teflon-lined caps

Other equipment or instrumentation may be used in future testing but may require optimization to achieve the desired separation and sensitivity.

## 2.5 Test Matrixes

The matrixes used during this method validation were ground water and surface water.

### **Ground water information:**

Ground water consists of unadulterated water from a 100-meter bedrock well prepared by filtering to remove any potential organic contaminants.

### **Surface water information:**

The surface water used for this method validation analysis was collected from the Taunton River (SMV Lot No. 05Feb19Wat-A) in Bridgewater, Massachusetts. The water was collected from an area of the river with approximately 30 to 60 cm of overlying water. Prior to use, the surface water was characterized by Smithers and was determined to have a pH of 6.57 and a dissolved oxygen content of 10.22 mg/L. All documentation relating to the preparation, storage, and handling is maintained by Smithers.

## 2.6 Preparation of Liquid Reagent Solutions

The volumes listed in this section were those used during the validation. For future testing, the actual volumes used may be scaled up or down as necessary.

A 20/80 acetonitrile/purified reagent water (v/v) liquid reagent solution was typically prepared by adding 200 mL of acetonitrile to 800 mL of purified reagent water. The solution was mixed well using a stir bar and stir plate for 5 minutes.

A 0.1% formic acid in purified reagent water mobile phase solution was typically prepared by adding 1.00 mL of formic acid to 1000 mL of purified reagent water. The solution was mixed well before use.

A 30/30/40 acetonitrile/methanol/purified reagent water (v/v/v) autosampler needle wash solution was typically prepared by combining 1500 mL of acetonitrile, 1500 mL of methanol, and 2000 mL of purified reagent water. The solution was mixed well before use.

## 2.7 Preparation of Stock Solutions

The volumes and masses listed in this section were those used during each separate validation. For future testing, the actual volumes and masses used may be scaled up or down as necessary.

Primary stock solutions were typically prepared as described in the table below:

Primary Stock ID	Amount Weighed (g), Net Weight	Amount Weighed (g), as Active Ingredient	Stock Solvent	Final Volume (mL)	Primary Stock Concentration (mg/L)	Primary Stock Use
<b>Test Substance</b>						
9932-1H	0.0508	0.0502	Acetonitrile	50.0	1000	Secondary stock solution
9932-1Z	0.5067	0.5011	Acetone	50.0	10,000	Secondary stock solution
9932-1AG	0.5063	0.5007	Acetone	50.0	10,000	Secondary stock solution
<b>Reference Substance</b>						
9945-1A	0.0504	0.0502	Acetonitrile	50.0	1000	Secondary stock solution
9945-1B	0.0502	0.0500	Acetonitrile	50.0	1000	Secondary stock solution
9945-1C	0.0502	0.0500	Acetonitrile	50.0	1000	Secondary stock solution
9945-1D	0.0504	0.0502	Acetonitrile	50.0	1000	Secondary stock solution

Secondary stock solutions were typically prepared as described in the table below:

Fortifying Stock ID	Fortifying Stock Concentration (mg/L)	Volume of Fortification (mL)	Final Volume (mL)	Stock Solvent	Stock ID	Stock Concentration (mg/L)	Stock Use
<b>Test Substance</b>							
9932-1H	1000	0.500	50.0	Acetonitrile	9932-1H-1	10.0	Sub-stock solutions
9932-1Z	10,000	0.0500	50.0	Acetonitrile	9932-1Z-1	10.0	Sub-stock solutions
9932-1AG	10,000	0.0500	50.0	Acetonitrile	9945-1AG-1	10.0	Sub-stock solutions
<b>Reference Substance</b>							
9945-1A	1000	0.500	50.0	Acetonitrile	9945-1A-1	10.0	Sub-stock solutions
9945-1B	1000	0.500	50.0	Acetonitrile	9945-1B-1	10.0	Sub-stock solutions
9945-1C	1000	0.500	50.0	Acetonitrile	9945-1C-1	10.0	Sub-stock solutions
9945-1D	1000	0.500	50.0	Acetonitrile	9945-1D-1	10.0	Sub-stock solutions

Sub-stock solutions were typically prepared as described in the table below:

Fortifying Stock ID	Fortifying Stock Concentration (mg/L)	Volume of Fortification (mL)	Final Volume (mL)	Stock Solvent	Stock ID	Stock Concentration (mg/L)	Stock Use
<b>Test Substance</b>							
9932-1Z-1	10.0	0.100	10.0	Acetonitrile	Tech Stk 2	0.100	LOQ recovery samples (ground water)
9932-1AG-1	10.0	1.00	10.0	Acetonitrile	Tech Stk 1	1.00	10X LOQ recovery samples (ground water)
9932-1H-1	10.0	1.00	10.0	Acetonitrile	Tech Stk 1	1.00	10X LOQ recovery samples (surface water)
		0.100	10.0	Acetonitrile	Tech Stk 2	0.100	LOQ recovery samples (surface water)
<b>Reference Substance</b>							
9945-1C-1	10.0	1.00	10.0	Acetonitrile	Ana Stk 1	1.00	Solvent-based calibration standards and matrix effects investigation samples (ground water)
9945-1D-1	10.0	1.00	10.0	Acetonitrile	Ana Stk 1	1.00	Solvent-based calibration standards for 10X samples (ground water)
9945-1A-1	10.0	1.00	10.0	Acetonitrile	Ana Stk 1	1.00	Solvent-based calibration standards and matrix effects investigation samples (surface water)
9945-1B-1	10.0	1.00	10.0	Acetonitrile	Ana Stk 1	1.00	Sub-stock solution (surface water)
Ana Stk 1	1.00	1.00	10.0	Acetonitrile	Ana Stk 2	0.100	Matrix-matched calibration standards (surface water)

All primary and secondary stock solutions were stored refrigerated (2 to 8 °C) in amber glass bottles fitted with Teflon-lined caps. Sub-stock solutions were prepared fresh on the day of use for daily use.

## 2.8 Preparation of Calibration Standards

Solvent-based calibration standards used in the quantitation of ground water and surface water samples were prepared in 20/80 acetonitrile/purified reagent water (v/v) by dosing with the 1.00 mg/L sub-stock solution to yield concentrations of 0.500, 1.00, 2.00, 3.00, 4.00, and 5.00 µg/L. Solvent-based calibration standards were used to quantify both fortified recovery samples and matrix effects samples for ground water, and matrix effects sample only for surface water.

Matrix-matched calibration standards used in the quantitation of surface water samples were prepared in matrix blank (see [Section 2.11](#)) by dosing with the 0.100 mg/L sub-stock solution to yield concentrations of 0.500, 1.00, 2.00, 3.00, 4.00, and 5.00 µg/L. Matrix-matched calibration standards were used to quantify fortified recovery samples for surface water.

## 2.9 Matrix Effect Investigation

The effects of matrix enhancement or suppression were evaluated through the assessment of matrix-matched and solvent-based calibration standards in the following manner. Calibration standards used to assess possible matrix effects were prepared in triplicate. One set was prepared in matrix blank final extract (see [Section 2.11](#)) and a second set was prepared in 20/80 acetonitrile/purified reagent water (v/v) by fortifying with the 1.00 mg/L sub-stock to yield a concentration of 1.00 µg/L. The preparation procedure for each separate matrix is outlined in the tables below.

### Ground water validation

Sample ID	Sample Type	Stock Concentration (mg/L)	Fortification Volume (mL)	Final Volume (mL)	Fortified Concentration (µg/L)
MM-Std A, B, & C	Matrix-matched calibration standard	1.00	0.0200	20.0 <sup>a</sup>	1.00
Std A, B, & C	Solvent-based calibration standard	1.00	0.0200	20.0 <sup>b</sup>	1.00

<sup>a</sup> Diluted with matrix blank final extract 12791.6319-M3

<sup>b</sup> Diluted with 20/80 acetonitrile/purified reagent water (v/v)

### Surface water validation

Sample ID	Sample Type	Stock Concentration (mg/L)	Fortification Volume (mL)	Final Volume (mL)	Fortified Concentration (µg/L)
MM-Std D, E, & F	Matrix-matched calibration standard	1.00	0.0200	20.0 <sup>a</sup>	1.00
Std D, E, & F	Solvent-based calibration standard	1.00	0.0200	20.0 <sup>b</sup>	1.00

<sup>a</sup> Diluted with matrix blank final extract 12791.6319-M2

<sup>b</sup> Diluted with 20/80 acetonitrile/purified reagent water (v/v)



## 2.10 Sample Fortification and Preparation

The recovery samples were prepared in two different matrixes (ground water and surface water) by fortification with stock solutions of dicloran at concentrations of 0.100 (LOQ) and 1.00 (10X LOQ) µg/L. Recovery samples for both matrixes were prepared separately (“de novo”) at these concentrations in separatory funnels. Five replicates were produced for each concentration level. Two samples of each matrix were left unfortified to serve as controls and were processed in the same fashion as the LOQ concentration recovery samples. In addition, one matrix blank (for the assessment of matrix effects) and one reagent blank were prepared for each sample set and processed in the same manner as the control samples. The preparation procedure for each separate matrix is outlined in the tables below.

### Ground water recovery samples

Sample ID 12791-6319-	Sample Type	Stock Concentration (mg/L)	Fortification Volume (mL)	Final Volume (mL)	Fortified Concentration (µg/L)
M3	Matrix Blank	NA <sup>a</sup>	NA	1000 <sup>b</sup>	0.00
27	Reagent Blank	NA	NA	100	0.00
28 & 29	Control	NA	NA	100	0.00
30, 31, 32, 33, & 34	LOQ	0.100	0.100	100	0.100
40, 41, 42, 42, & 44	10X LOQ	1.00	0.100	100	1.00

<sup>a</sup> NA = Not Applicable

<sup>b</sup> Volume increased for use in matrix effects assessment.

### Surface water recovery samples

Sample ID 12791-6319-	Sample Type	Stock Concentration (mg/L)	Fortification Volume (mL)	Final Volume (mL)	Fortified Concentration (µg/L)
M2	Matrix Blank	NA <sup>a</sup>	NA	1000 <sup>b</sup>	0.00
14	Reagent Blank	NA	NA	100	0.00
15 & 16	Control	NA	NA	100	0.00
17, 18, 19, 20, & 21	LOQ	0.100	0.100	100	0.100
22, 23, 24, 25, & 26	10X LOQ	1.00	0.100	100	1.00

<sup>a</sup> NA = Not Applicable

<sup>b</sup> Volume increased for use in matrix effects assessment.

## 2.11 Extraction of Samples

All samples were taken to pH 2 (measured with pH paper) with phosphoric acid prior to extraction per the table below. Samples were extracted twice with 100 mL of the extraction

solvent, dichloromethane. After the final extraction, the post-extraction water was drained from the separatory funnels and they were rinsed well with an additional 50 mL of dichloromethane. The rinsate was added to the sample extracts in round bottom flasks. The extracts were taken to a low volume (approximately 2.00 mL) by rotary evaporation using minimal heating (<35 °C). To remove any residual water from the round bottom flasks, a 100-mL aliquot of acetone was added to each round bottom flask and the extracts were again taken to low volume (approximately 5.00 mL) by rotary evaporation using minimal heating (<35 °C). The concentrated extracts were transferred to glass centrifuge tubes (or larger glass vials depending on the final volume) and taken to incipient dryness (approximately 100 µL) under a gentle stream of nitrogen at room temperature. This transfer was performed first with dichloromethane followed by acetone to better rinse the flasks. To minimize the potential for losses of the test substance, careful consideration was taken to ensure that the samples did not go dry or approach dryness at any point in the sample concentration process. Following concentration, an aliquot of 20/80 acetonitrile/purified reagent water (v/v) was added to each sample, which was vortexed for 30 seconds and sonicated for 5 minutes to mix well and aid in dissolution of the reconstituted extract. For the ground water 10X LOQ samples only, the recovery samples were further diluted with 20/80 acetonitrile/purified reagent water. For the surface water 10X LOQ samples only, the recovery samples were further diluted using matrix blank final extract. The extraction and dilution procedures for each separate matrix is outlined in the tables below.

### Ground water recovery samples

Sample ID 12791-6319-	Sample Type	Fortified Concentration (µg/L)	Sample Volume (mL)	Volume of Phosphoric Acid Added (mL)	Reconstituted Volume <sup>a</sup> (mL)	Sample Volume (mL)	Final Volume <sup>a</sup> (mL)	Dilution Factor
M3	Matrix Blank	0.00	1000	0.0200	100 <sup>b</sup>	NA <sup>c</sup>	NA	0.100
27	Reagent Blank	0.00	100	0.0200	10.0	NA	NA	0.100
28 & 29	Control	0.00	100	0.0200	10.0	NA	NA	0.100
30, 31, 32, 33, & 34	LOQ	0.100	100	0.0200	10.0	NA	NA	0.100
40, 41, 42, 43, & 44	10X LOQ	1.00	100	0.0200	10.0	3.00	10.0	0.333

<sup>a</sup> Diluted with 20/80 acetonitrile/purified reagent water (v/v)

<sup>b</sup> Volume increased for use in matrix effects assessment.

<sup>c</sup> NA = Not Applicable

**Surface water recovery samples**

Sample ID 12791-6319-	Sample Type	Fortified Concentration (µg/L)	Sample Volume (mL)	Volume of Phosphoric Acid Added (mL)	Reconstituted Volume <sup>a</sup> (mL)	Sample Volume (mL)	Final Volume <sup>b</sup> (mL)	Dilution Factor
M2	Matrix Blank	0.00	1000	0.0600	100 <sup>c</sup>	NA <sup>d</sup>	NA	0.100
14	Reagent Blank	0.00	100	0.0600	10.0	NA	NA	0.100
15 & 16	Control	0.00	100	0.0600	10.0	NA	NA	0.100
17, 18, 19, 20, & 21	LOQ	0.100	100	0.0600	10.0	NA	NA	0.100
22, 23, 24, 25, & 26	10X LOQ	1.00	100	0.0600	10.0	0.300	1.00	0.333

<sup>a</sup> Diluted with 20/80 acetonitrile/purified reagent water (v/v)

<sup>b</sup> Diluted with matrix blank final extract (12791.6319-M2). For the matrix blank sample, due to the increased sample volume, an emulsion formed during extraction. This was eliminated by centrifuging aliquots of the solvent layer at 3000 rpm for 10 minutes during the first extraction.

<sup>c</sup> Volume increased for use in matrix effects assessment.

<sup>d</sup> NA = Not Applicable

**2.12 Analysis****2.12.1 Instrumental Conditions**

The LC-MS/MS analysis was conducted utilizing the following instrumental conditions:

**LC parameters:**

Column: Phenomenex Kinetex, 2.6 µm phenyl-hexyl, 3 × 50 mm

Mobile Phase A: 0.1% formic acid in reagent grade water

Mobile Phase B: 0.1% formic acid in acetonitrile

Gradient:

Time (min.)	Flow rate (mL/min.)	Solvent A (%)	Solvent B (%)
0.01	0.500	80	20
0.50	0.500	80	20
3.00	0.500	0.0	100
4.00	0.500	0.0	100
4.10	0.500	80	20
5.50 <sup>a</sup>	0.500	80	20

Run Time: 5.5<sup>a</sup> minutes

(<sup>a</sup> Please note run time shown here is for surface water; the run time for the ground water analysis was 5.1 minutes. Both are demonstrated to be suitable for analysis of dicloran in water)

Autosampler Wash Solvent: 30/30/40 acetonitrile/methanol/reagent grade water (v/v/v)

Column Temperature: 40 °C

Sample Temperature: 10 °C

Injection Volume: 100 µL

Retention Times: approximately 3.1 minutes (ground water)  
 approximately 3.0 minutes (surface water)

**MS parameters:**

Instrument: MDS Sciex API 5000 mass spectrometer  
 Ionization Mode: Positive (+) ESI  
 Ion Spray Voltage: 5500 V  
 Scan Type: MRM  
 Dwell Time: 200 milliseconds  
 Source Temperature: 550 °C  
 Curtain Gas: 20.0  
 Ion Source – Gas 1 / Gas 2: 60.0 / 70.0  
 Collision Gas: 12.0  
 Entrance Potential: 10.0  
 Declustering Potential: 40.0  
 Resolution Q1/Q3: Unit/Unit

	Primary Transition	Confirmatory Transition
Q1/Q3 Masses (amu):	207.1/190.0	207.1/160.0
Collision Energy:	22.0	35.0
Collision Cell Exit Potential:	19.0	42.0

Other instrumentation may be used but may require optimization to achieve the desired separation and sensitivity. It is important to note that the parameters above have been established for this particular instrumentation and may not be applicable for other similar equipment that may be used.

**2.12.2 Preparation of Calibration Standard Curve**

Two sets of calibration standards were analyzed with each sample set. Calibration standards were interspersed among analysis of the recovery samples, every two to six injections. Injection of recovery samples and calibration standards onto the chromatographic system was performed by programmed automated injection.

### **2.13 Evaluation of Accuracy, Precision, Specificity, and Linearity**

The accuracy was reported in terms of percent recovery of the fortified recovery samples. Recoveries of 70.0 to 110% (for the individual mean concentrations) are acceptable. The precision was reported in terms of the relative standard deviation (RSD) for the recovery samples and retention times. RSD values less than or equal to 20% were considered acceptable for the recovery samples and RSD values less than or equal to 2% were considered acceptable for the retention times. Specificity of the method was determined by examination of the control samples for peaks at the same retention times as dicloran which might interfere with the quantitation of the analytes. Linearity of the method was determined by the coefficient of determination ( $r^2$ ), y-intercept, and slope of the regression line.

### **2.14 Limit of Quantitation (LOQ)**

The method was validated at the LOQ. This was defined as the lowest fortification level. Blank values (reagent blanks and untreated control samples) did not exceed 30% of the LOQ.

### **2.15 Limit of Detection (LOD) and Method Detection Limit (MDL)**

The LOD was calculated using three times the signal-to-noise value of the control samples. Representative calculations for the LOD can be found in [Section 3.0](#).

The MDL was defined as the lowest concentration in test samples which can be detected based on the concentration of the low calibration standard and the dilution factor of the control solutions. Representative calculations for the MDL can be found in [Section 3.0](#).

## **3.0 CALCULATIONS**

A calibration curve was constructed by plotting the analyte concentration ( $\mu\text{g/L}$ ) of the calibration standards against the peak area of the analyte in the calibration standards. The equation of the line (equation 1) was algebraically manipulated to give equation 2. The

concentration of test substance in each recovery sample was calculated using the slope and intercept from the linear regression analysis, the detector response, and the dilution factor of the recovery sample. Equations 2 and 3 were then used to calculate measured concentrations and analytical results.

$$(1) \quad y = mx + b$$

$$(2) \quad DC(x) = \frac{(y - b)}{m}$$

$$(3) \quad A = DC \times DF$$

where:

x	=	analyte concentration
y	=	detector response (peak area) from the chromatogram
b	=	y-intercept from the regression analysis
m	=	slope from the regression analysis
DC (x)	=	detected concentration ( $\mu\text{g/L}$ ) in the sample
DF	=	dilution factor (final volume of the sample divided by the original sample volume)
A	=	analytical result ( $\mu\text{g/L}$ ), concentration in the original sample

The LOD was calculated using the following equation:

$$(4) \quad LOD = ((3 \times (N_{ctl}))/Res_{PLS}) \times Con_{CLS} \times DF_{CNTL}$$

where:

$N_{ctl}$	=	mean noise in height of the control samples (or blanks)
$Res_{PLS}$	=	mean response in height of the two low calibration standards
$Con_{CLS}$	=	concentration of the low calibration standard
$DF_{CNTL}$	=	dilution factor of the control samples (smallest dilution factor used, i.e., 0.100)
LOD	=	limit of detection for the analysis

The MDL is defined as the lowest concentration that can be detected by this method in test solution samples. The MDL is calculated (equation 5) based on the concentration of the low calibration standard and the dilution factor of the control samples.

$$(5) \quad \text{MDL} = \text{MDL}_{\text{LCAL}} \times \text{DF}_{\text{CNTL}}$$

where:

- $\text{MDL}_{\text{LCAL}}$  = lowest concentration calibration standard (0.500  $\mu\text{g/L}$ )  
 $\text{DF}_{\text{CNTL}}$  = dilution factor of the control samples (smallest dilution factor used, i.e., 0.100)  
 $\text{MDL}$  = method detection limit reported for the analysis  
(0.500  $\mu\text{g/L} \times 0.100 = 0.0500 \mu\text{g/L}$ )

12791.6319

**Environmental Chemistry Method: Validation of the Analytical Method for the Determination of Dicloran in Ground Water and Surface Water by LC-MS/MS****1.0 INTRODUCTION**

The purpose of this study is to validate an analytical method used to determine the content of Dicloran in two aqueous matrices (ground water and surface water) by LC-MS/MS. The analytical method will be validated with regards to accuracy, precision, specificity, linearity, limit of quantitation (LOQ), limit of detection (LOD), method detection limit (MDL), and confirmation of identification.

**2.0 JUSTIFICATION OF THE TEST SYSTEM**

This study is conducted to support the registration of the test substance.

The method validations described in this protocol are designed to conform to EPA guideline OCSPP 850.6100: Environmental Chemistry Methods and Associated Independent Laboratory Validation and SANCO/3029/99 rev.4: Guidance for generating and reporting methods of analysis in support of pre-registration data. The study will be conducted under Good Laboratory Practices (GLP) regulations and principles as described in 40CFR160 and as accepted by OECD principles of GLP (OECD, 1998).

**3.0 TEST SUBSTANCE****3.1 Test Substances**

Upon arrival at Smithers Viscient, the test substances (and the reference substances) will be received by the Test Material Center. Records will be maintained in accordance with GLP requirements, and a Chain-of-Custody established. The condition of the external packaging of the test and reference substances will be recorded and any damage noted. The packaging will be removed, the primary storage container inspected for leakage or damage, and the condition recorded. Any damage will be reported to the Sponsor and/or manufacturer.

Each test and reference substance will be given a unique sample ID number and stored under the conditions specified by the Sponsor or manufacturer. The following information should be provided by the Study Sponsor, if applicable: test substance lot or batch number, test substance purity, water solubility (pH and temperature of solubility determination), vapor pressure, storage stability, methods of analysis of the test substance in water, MSDS, and safe handling procedures, and a verified expiration or reanalysis date.

**3.2 Test Matrices****3.2.1 Ground Water**

The ground water used in the study will be filtered well water collected on site at Smithers Viscient, Wareham, MA. This will be prepared by filtering to remove any potential organic contaminants. All documentation relating to the preparation, storage and handling will be maintained by Smithers Viscient.



12791.6319

### **3.2.2 Surface Water**

The surface water used for the method validation will be collected from river water in Massachusetts. All documentation relating to the collection, preparation, storage and handling will be maintained by Smithers Viscient.

### **3.3 Reagents**

Highly pure reagents will be used throughout the study. The actual reagent grade will be depending on the manufacturer's designation. Generally these reagents will have grades, such as high purity solvent, ACS grade, or Select. The reagents used are recorded along with test chemical information at the time of preparation.

## **4.0 VALIDATION DESIGN**

The test design will consist of two aqueous matrices (ground water and surface water) fortified with each test substance at two concentrations with five replications at the target LOQ and five replicates at 10x LOQ level for each matrix. The procedural blank will be reagent blank without matrix. The control matrix for the validation will be untreated matrix representing ground water or surface water. The validation study levels (approximate concentrations) for each test substance are:

1. Procedural blank-reagent blank	0.0 µg/L
2. Matrix blank-control matrix	0.0 µg/L
3. Control matrix fortified at LOQ	0.10 µg/L
4. Control matrix fortified at 10 x LOQ	1.0 µg/L

### **4.1 Accuracy and Precision**

The accuracy of the analytical method will be determined by applying the method to five samples at the LOQ and five samples at 10X LOQ for each test substance. Accuracy will be reported as the mean recovery at each fortification level. Mean recoveries in the range 70 – 110% of nominal concentrations of the target analyte in the fortified samples will be considered acceptable.

The precision will be calculated for the fortified samples in terms of the relative standard deviation (RSD or coefficient of variation (CV)) calculated for the retention time, peak area based quantitation (i.e., µg/L), and the observed recovery values at each fortification level (n = 5 per level). The retention time should have a RSD of less than or equal to 2%. The RSD of the peak area based quantitation (i.e., µg/L) should be less than or equal to 20% per level. The RSD of the recovery values should be less than or equal to 20% per level as well.

### **4.2 Specificity**

The specificity of the method will be determined by applying the method to the appropriate number of reagent blank (n=1) and control matrix samples (n=2). Chromatograms will be obtained for the control samples and examined for peaks that might interfere with the

12791.6319

quantitation of the analyte(s) peak of interest. Peaks attributable to the test substance(s) should be sufficiently resolved from any peaks found in the samples of control matrix to enable quantification. Blank values (including procedural blanks and untreated samples) should not exceed 30% of the LOQ. If this is exceeded, detailed justification is required. Unequivocal identification of the target analyte will be achieved by LC-MS/MS primary and confirmatory analysis.

#### **4.3 Regression Analysis**

Quantitative analysis will be achieved with the aid of a calibration curve. The calibration curve will be constructed using a minimum of five analytical standards and will extend over a range appropriate to the lowest and highest nominal concentrations of the target analyte in relevant analytical solutions  $\pm$  at least 20%.

The calibration data will be subjected to regression analysis; a plot of analyte concentration versus detector response will be included in the report along with the correlation coefficient ( $r$ ) and the equation describing the curve. The linearity of the detector response will be assessed according to the strength of the correlation coefficient: this should be  $\geq 0.995$  (or coefficient of determination,  $r^2 \geq 0.990$ ). If non-linear calibration is used an explanation will be provided.

#### **4.4 Limits of Quantitation (LOQ)**

The method will be validated at the limit of quantitation (LOQ). This will be defined as the lowest fortification level. Blank values (reagent blanks and untreated control samples) should not exceed 30% of the LOQ. If this is exceeded, it will be discussed with the Sponsor and detailed justification provided.

#### **4.5 Limits of Detection (LOD)**

The Limits of Detection (LOD) will be calculated using three times the signal-to-noise value of the control samples. The method detection limit (MDL) will be set at the lowest concentration that can be detected in sample test solutions. The value is calculated based on the concentration of the low calibration standard and the dilution factor of the control samples.

#### **4.6 Matrix Effects Determination**

Determination of LC-MS/MS matrix effects will be evaluated through the assessment of solvent-based and matrix-matched calibration standards for both primary and confirmatory transitions. Matrix effects should be evaluated at the LOQ level for each test substance. Only if experiments clearly demonstrate that matrix effects are not significant (i.e. <20%), calibration with standards in solvent may be used.

#### **4.7 Confirmatory Analyses**

Unequivocal identification of the target analytes will be achieved by LC/MS-MS using a primary quantitation ion and secondary quantitation/confirmatory ion. All of the required elements need to be met for this confirmatory method with full method validation results generated for both transitions. For triple-quad MS methods, the confirmation method would be where a confirmatory (secondary) product ion will be used for quantification. The confirmatory ion analysis will also adhere to the aforementioned method specifications (Sections 4.1 - 4.6 above)

12791.6319

#### **5.0 PROCEDURE FOR THE IDENTIFICATION OF THE TEST SYSTEM**

The test system will be defined as the fortified recovery samples. The fortified recovery samples will be labeled as defined in Section 4.0 and each sample replicate will be assigned a unique identifier. Processing of fortified recovery samples will be performed at a lab station labeled with the study number.

#### **6.0 CONTROL OF BIAS**

Bias will be effectively controlled through techniques such as, but not limited to, preparation of replicate samples and replicate analysis.

#### **7.0 RECORDS TO BE MAINTAINED**

Records to be maintained will include, but will not be limited to, correspondence and other documents relating to the interpretation and evaluation of data as well as all raw data and documentation generated as a result of the study.

#### **8.0 REPORTING**

The raw data generated at Smithers Viscient will be peer-reviewed and the final report will be reviewed by the Study Director. All values will be reported to various levels of significance depending on the accuracy of the measuring devices employed during any one process. The Quality Assurance Unit will inspect the final report to confirm that the methods, procedures, and observations are accurately and completely described, that the reported results accurately and completely reflect the raw data generated at Smithers Viscient and to confirm adherence with the study protocol. A single copy of the draft report will be submitted to the Sponsor for review. The report will be finalized according to standard operating procedures. The final report will meet the formatting requirements of EPA's PR Notice 2011-3. All reports will include, but will not be limited to reporting requirements presented in Ecological Effects Test Guidelines OCSPP 850.6100 (U.S. EPA, 2012) and SANCO/3029/99 rev.4 along with the following information:

- The report and project numbers from Smithers Viscient and Sponsor Study number (if any).
- Laboratory and site, dates of testing and personnel involved in the study, e.g. Program Coordinator (if applicable), Study Director and Principal Investigator.
- Identification of the test substance including chemical name, additional designations (e.g., trade name), chemical designation (CAS number), empirical formula, molecular structure, manufacturer, lot or batch number, degree of purity of test substance (percent test chemical) (Sponsor supplied, if available).
- A full description of the experimental design and procedures followed and a description of the test equipment used.
- The determined accuracy, precision, specificity, linearity, LOQ, LOD and MDL, and confirmation of identification.

12791.6319

- The mathematical equations and statistical methods used in generating and analyzing the data as well as calculations using these equations. Tabular and graphical representations (if appropriate) of the data.
- Description of any problems experienced and how they were resolved.
- Good Laboratory Practice (GLP) Compliance Statement signed by the Study Director.
- Date(s) of Quality Assurance reviews, and dates reported to the Study Director and management, signed by the Quality Assurance Unit.
- Location of raw data and report.
- A copy of the study protocol and study amendments, if any.

#### 9.0 PROTOCOL AMENDMENTS

All amendments to the approved protocol must be documented in writing and signed by both the Study Director and the Sponsor's contact or representative. Protocol amendments and deviations must include the reasons for the change and the predicted impact of the change on the results of the study, if any.

#### 10.0 GOOD LABORATORY PRACTICES

All test procedures, documentation, records and reports will comply with the U.S. Environmental Protection Agency's Good Laboratory Practices as set forth under the Federal Insecticide, Fungicide and Rodenticide Act (40 CFR, Part 160) and as compatible with OECD Principles of Good Laboratory Practice ([OECD, 1998](#)).

#### 11.0 REFERENCES

- European Commission, 2000. Residues: Guidance for the generating and reporting methods of analysis in support of pre-registration data requirements for Annex II (part V, section 4) and Annex III (part A, section 5) of Directive 91/414. SANCO/3029/99 rev.4.
- OECD, 1998. OECD Series on Principles of Good Laboratory Practice and Compliance Monitoring. Number 1. OECD Principles on Good Laboratory Practice (as revised in 1997). Environment Directorate Chemicals Group and Management Committee. ENV/MC/CHEM(98)17. OECD Paris. France. 41 pp.
- U.S. EPA, 1989. Federal Insecticide, Fungicide and Rodenticide Act (FIFRA); Good Laboratory Practice Standards; Final Rule (40 CFR, Part 160); FR: 8/17/89; pp. 34052. U.S. Environmental Protection Agency, Washington, D.C.
- U.S. EPA, 2011. Pesticide Registration (PR) Notice 2011-3 Standard Format for Data Submitted Under the Federal Insecticide, Fungicide, and Rodenticide Act (FIFRA) and Certain Provisions of the Federal Food, Drug, and Cosmetic Act (FFDCA). US Environmental Protection Agency Office of Pesticide Programs. November 30, 2011.

12791.6319

U.S. EPA, 2012. Office of Chemical Safety and Pollution Prevention. Ecological Effects Guideline, OCSP 850.6100. Environmental Chemistry Methods and Associated Independent Laboratory Validation. EPA 712-C-001. January 2012. U.S. Environmental Protection Agency, Washington, D.C.