



**US Environmental Protection Agency
Office of Pesticide Programs**

EPA MLB SOP MB-37-00:

**Neutralization Confirmation for Evaluating the Efficacy of Liquid Antimicrobials against
Candida auris using the OECD Quantitative Method on Hard, Non-Porous Surfaces**

Date: 03/21/17

*Please see EPA's companion interim guidance document for *Candida auris* under the Guidance tab.*

1 **Neutralization Confirmation for Evaluating the Efficacy of Liquid**
2 **Antimicrobials against *Candida auris* using the OECD Quantitative Method**
3 **on Hard, Non-Porous Surfaces**
4

5 **I. Overview**

6 A. This document describes a quantitative procedure for verifying the effectiveness
7 of the neutralizer (refer to reference A). Verify neutralization using the highest
8 concentration of test substance if there are multiple concentrations being
9 evaluated.

10 B. In brief, the test substance is first mixed with a candidate neutralizer. A diluted
11 suspension of the test organism is then added to the reaction mixture; if desired,
12 additional evaluations may be conducted using the test organism as dried
13 inoculum on a carrier. The neutralization process is deemed acceptable if the
14 criteria outlined in section II of this document are met.

15 **II. Measure of effectiveness**

16 A. For the assay to be considered valid, ensure that:

17 1. The recovered number of Colony Forming Units (CFU) in the **Titer**
18 **Control** (see section IV.D.3) using *Final Test Suspension B* yields 20-200
19 CFU per vessel.

20 B. For determining and verifying the effectiveness of the neutralizer, ensure that:

21 1. The recovered number of CFU in the **Neutralizer Toxicity Control** (see
22 section IV.D.2) is at least 50% of the **Titer Control** (see section IV.D.3).
23 A count lower than 50% indicates that the neutralizer is harmful to the
24 test organism. Note: counts higher than the **Titer Control** (e.g., 120% of
25 the **Titer Control**) are also deemed valid.

26 2. The recovered number of CFU in the **Neutralizer Effectiveness**
27 treatment (see section IV.D.1) is at least 50% of the **Titer Control**; this
28 verifies effective neutralization. Note: counts higher than the **Titer**
29 **Control** (e.g., 120% of the **Titer Control**) are also deemed valid.

30 C. For both the suspension-based assay and the carrier-based assay, meet the criteria
31 in sections II.A and II. B. If the criteria are not met, identify and verify another
32 neutralizer or mixture of neutralizers.

33 **III. Special Apparatus and Materials**

34 A. Refer to section III of EPA MLB SOP MB-35 (OECD Quantitative Method for
35 Evaluating the Efficacy of Liquid Antimicrobials against *Candida auris* on Hard,
36 Non-Porous Surfaces).

37

38 **IV. Procedure and Analysis**

39 A. Preparation and sterilization of carriers (if using the carrier-based approach)

- 40 1. Refer to section IV.A of EPA MLB SOP MB-35 (OECD Quantitative
41 Method for Evaluating the Efficacy of Liquid Antimicrobials against
42 *Candida auris* on Hard, Non-Porous Surfaces).

43 B. Preparation of test organisms

- 44 1. Refer to section IV.B of EPA MLB SOP MB-35 (OECD Quantitative
45 Method for Evaluating the Efficacy of Liquid Antimicrobials against
46 *Candida auris* on Hard, Non-Porous Surfaces).

- 47 2. Prepare *Test Suspension A (without soil load)*. Serially dilute the test
48 microbial suspension with PBS (e.g., serially dilute cultures through 10^{-4}
49 or 10^{-5}). Select appropriate dilutions of *Test Suspension A* so that after
50 the addition of the soil load, *Final Test Suspension B* will achieve an
51 average challenge of 20-200 CFU per 10 μ L. Use *Test Suspension A*
52 within 4 hours of preparation.

- 53 i. If performing the assay with dried carriers, prior testing may be
54 required to account for differences in the loss of viability of the
55 different test organisms upon drying.

- 56 3. Prepare *Final Test Suspension B (with soil load)*. Prepare the soil load:
57 vortex each component and combine 25 μ L bovine serum albumin
58 (BSA), 35 μ L yeast extract, and 100 μ L of mucin; mix well. Combine
59 340 μ L of *Test Suspension A* and 160 μ L of the soil load (SL). The test
60 microbial suspension with soil load should provide an average challenge
61 of 20-200 CFU/tube.

- 62 i. If performing the assay with dried carriers, ensure an average
63 challenge of 20-200 CFU/carrier after drying.

- 64 ii. Note: Two separate serial dilutions of *Test Suspension A* may be
65 used to prepare two different concentrations of *Final Test*
66 *Suspension B* to ensure at least one dilution with a challenge of
67 20-200 CFU. If performing the assay with dried carriers, the use
68 of two separate dilutions results in a total of 20 carriers to be
69 processed; however, the dilutions may be evaluated separately.

- 70 4. If desired, use optical density (OD @ 650 nm) to create a calibration
71 curve to estimate the number of viable organisms in *Test Suspension A*.

72 C. Carrier inoculation (for carrier-based assay)

- 73 1. Inoculate at least 13 carriers with 10 μ L of *Final Test Suspension B* (per
74 concentration of *Final Test Suspension B*) using a positive displacement
75 pipette with 10 μ L tips. Refer to the section IV.D of EPA MLB SOP

- 76 MB-35 (OECD Quantitative Method for Evaluating the Efficacy of
77 Liquid Antimicrobials against *Candida auris* on Hard, Non-Porous
78 Surfaces) for carrier drying instructions.
- 79 2. After drying, evaluate the inoculated carriers per section IV.E of this
80 document.
- 81 D. Suspension-based assay
- 82 1. ***Treatment 1: Neutralizer Effectiveness.*** Add 50 μ L of the test substance
83 to each of three vessels. At timed intervals, add 10 mL neutralizer to each
84 vessel and briefly swirl. After 10 s, gently add 10 μ L *Final Test*
85 *Suspension B* using a micropipette to each vessel and briefly vortex.
86 Proceed with section IV.F.
- 87 2. ***Treatment 2: Neutralizer Toxicity Control.*** Add 10 mL neutralizer to each
88 of three reaction vessels. At timed intervals, add 10 μ L of *Final Test*
89 *Suspension B* using a micropipette gently to each vessel and briefly vortex.
90 Proceed with section IV.F.
- 91 3. ***Treatment 3: Titer Control.*** Add 10 mL PBS to each of three reaction
92 vessels. At timed intervals, add 10 μ L of *Final Test Suspension B* with a
93 micropipette gently to each vessel and briefly vortex. Proceed with section
94 IV.F.
- 95 E. Alternative inoculated carrier-based assay
- 96 1. ***Treatment 1: Neutralizer Effectiveness.*** Add 50 μ L of the test substance
97 to each of three reaction vessels. At timed intervals, add 10 mL
98 neutralizer to each vial and briefly swirl. After 10 s, gently add one dried
99 carrier inoculated with *Final Test Suspension B* to each vessel and vortex
100 for 30 ± 2 s. Proceed with section IV.F.
- 101 2. ***Treatment 2: Neutralizer Toxicity Control.*** Add 10 mL neutralizer to each
102 of three reaction vessels. At timed intervals, add one dried carrier
103 inoculated with *Final Test Suspension B* gently to each vessel and vortex
104 for 30 ± 2 s. Proceed with section IV.F.
- 105 3. ***Treatment 3: Titer Control.*** Add 10 mL PBS to each of three reaction
106 vessels. At timed intervals, add one dried carrier inoculated with *Final*
107 *Test Suspension B* gently to each vessel and vortex for 30 ± 2 s. Proceed
108 with section IV.F.
- 109 F. Processing and recovery
- 110 1. Hold the mixtures from IV.D and IV.E for 10 ± 1 min at room temperature
111 (20 - 25°C). Steps (e.g., addition of organism and neutralizer) should be
112 conducted at timed intervals (e.g., 30 s intervals for suspension-based

- 113 assay, 1 min intervals for dried carrier-based assay) to ensure consistent
114 time of contact.
- 115 2. At the conclusion of the holding period, briefly vortex each vessel (for the
116 suspension-based assay) and pass each mixture through a separate, pre-
117 wetted 0.45 μm polyethersulfone (PES) membrane filter.
- 118 i. If performing the assay with dried carriers, vortex each vessel for
119 30 ± 2 s at the conclusion of the holding period. Use a magnet to
120 prevent carriers from falling onto the filter membrane.
- 121 3. Wash each vessel with approximately 20 mL PBS and briefly vortex; filter
122 the wash through the same filter membrane. Finish the filtering process by
123 rinsing the inside of the funnel unit with about 40 mL of PBS and filter the
124 rinsing liquid through the same filter membrane.
- 125 i. Note: Initiate filtration as soon as possible (e.g., within 30 min).
126 Two analysts are recommended to perform vortexing and filtration
127 steps to reduce holding time after vortexing.
- 128 4. Remove the membrane aseptically with sterile forceps and place it
129 carefully over the surface of the recovery medium (SDEA). Avoid trapping
130 air bubbles between the filter and the agar surface. Incubate the plates for
131 48-72 h at $30\pm 1^\circ\text{C}$.
- 132 5. Examine the plates beginning at 48 ± 4 h after incubation and conduct final
133 count at 72 ± 4 h. Calculate the average CFU for each set of test conditions.

134 V. Data Analysis and Calculations

- 135 A. For the suspension-based assay (see section IV.D), compare the average CFU of
136 the **Titer Control** with the average CFU of the **Neutralizer Toxicity Control**
137 and **Neutralizer Effectiveness** treatment. Compare results from the dried
138 carrier-based assay (see section IV.E) in the same manner.

139 VI. Attachments


- 140 A. OECD Neutralization Assay Flow Chart

141 VII. References


- 142 A. OECD Guidance Document: Quantitative Method for Evaluating Bactericidal
143 Activity of Microbicides Used on Hard, Non-Porous Surfaces (January 29, 2013).
- 144 B. EPA MLB SOP MB-35: OECD Quantitative Method for Evaluating the Efficacy
145 of Liquid Antimicrobials against *Candida auris* on Hard, Non-Porous Surfaces

OECD Neutralization Assay Flow Chart


Treatment 1



Treatment 2



Treatment 3



Alternatively, perform the assay using dried-carriers in place of the liquid suspension.