

**ESS Method 310.1:  
Ortho-Phosphorus, Dissolved  
Automated, Ascorbic Acid**

**Environmental Sciences Section  
Inorganic Chemistry Unit  
Wisconsin State Lab of Hygiene  
465 Henry Mall  
Madison, WI 53706**

**Revised October 1992**



# **ESS Method 310.1: Ortho-Phosphorus, Dissolved Automated, Ascorbic Acid**

## **1.0 Scope and Application**

- 1.1 This method may be used to determine concentrations of orthophosphate in most waters and wastewater in the range from 0.002-0.200 mg P/L. The concentration range may be extended to 0.2-2.00 mg P/L by utilizing a dilution loop.
- 1.2 Approximately 30 samples per hour can be analyzed.

## **2.0 Summary of Method**

Ammonium molybdate and antimony potassium tartrate react in an acid medium with dilute solutions of orthophosphate-phosphorus to form an antimony-phospho-molybdate complex. This complex is reduced to an intensely blue-colored complex by ascorbic acid. The color is proportional to the phosphorus concentration.

## **3.0 Sample Handling and Preservation**

Samples must be filtered through a 0.45 µm filter, cooled to 4 °C and analyzed as soon as possible.

## **4.0 Interferences**

- 4.1 Barium, lead, and silver interfere by forming a precipitate.
- 4.2 The interference from silica, which forms a pale-blue complex is small and can usually be considered negligible.
- 4.3 Arsenate is determined similarly to phosphorus and should be considered when present in concentrations higher than phosphorus.

## **5.0 Apparatus**

Technicon AutoAnalyzer II system consisting of:

- 5.1 Sampler IV with a 30/h (2:1) Cam
- 5.2 Analytical manifold (orthophosphate in seawater) with internal heating bath at 37.5°C and dilution loop
- 5.3 Proportioning pump III
- 5.4 Colorimeter equipped with 50 mm flowcells and 880 nm interference filters

5.5 Printer/Plotter

## 6.0 Reagents

- 6.1 Stock Solution A; Sulfuric acid solution, 4.9 N: Add 136 mL concentrated  $\text{H}_2\text{SO}_4$  to 800 mL Milli-Q water. Cool and dilute to 1 L with Milli-Q water.
- 6.2 Stock Solution B; Ammonium molybdate solution: Dissolve 40 g of  $((\text{NH}_4)_6 \text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O})$  in 900 mL Milli-Q water and dilute to 1 L. Store at  $4^\circ\text{C}$ .
- 6.3 Stock Solution C; Ascorbic acid: Dissolve 9 g of ascorbic acid ( $\text{C}_6\text{H}_8\text{O}_6$ ) in 400 mL Milli-Q water and dilute to 500 mL. Store at  $4^\circ\text{C}$ . Keep well stoppered. Prepare fresh monthly or as needed.
- 6.4 Stock solution D; Antimony potassium tartrate: Dissolve 3.0 g of  $(\text{K}(\text{SbO})\text{C}_4\text{H}_4\text{O}_6 \cdot \frac{1}{2}\text{H}_2\text{O})$  in 800 mL Milli-Q water and dilute to 1 L. Store at  $4^\circ\text{C}$ .
- 6.5 Combined color reagent: Combine the following solutions in order, mixing after each addition: (Prepare fresh daily)

Stock A, 6.1 (4.9 N $\text{H}_2\text{SO}_4$ )	50 mL
Stock B, 6.2 (Ammonium molybdate solution)	15 mL
Stock C, 6.3 (Ascorbic acid solution)	30 mL
Stock D, 6.4 (Antimony-tartrate solution)	5 mL

- 6.6 Water diluent solution: Add 4.0 g sodium lauryl sulfate and 5 g NaCl per L of Milli-Q water.
- 6.7 Stock phosphorus standard: Dissolve 0.4394 g of Potassium phosphate monobasic ( $\text{KH}_2\text{PO}_4$ ) (dried at  $105^\circ\text{C}$  for one hour) in 900 mL Milli-Q water. Add 2 mL of concentrated  $\text{H}_2\text{SO}_4$  and dilute to 1 L. 1.0 mL = 0.100 mg P (100 mg P/L).
- 6.8 Standard phosphorus solution 1: Dilute 100.0 mL of stock solution (6.7) to 500 mL with Milli-Q water. 1.0 mL = 0.020 mg P (20 mg P/L).
- 6.9 Standard phosphorus solution 2: Dilute 10.0 mL of stock solution (6.7) to 1 L. 1.0 mL = 0.001 mg P (1.0 mg P/L).
- 6.10 Working standard solutions:
- 6.10.1 Low Range (0.002-0.200 mg P/L): Prepare the following standards by diluting suitable volumes of standard solution 2 (6.9) to appropriate volumes with Milli-Q water:

<u>mg P/L</u>	<u>mL of standard solution 2</u>
0.005	1.0/200 mL
0.050	5.0/100 mL
0.100	50/500 mL
0.150	15/100 mL
0.200	40/200 mL



6.10.2 High Range (0.02-2.00 mg P/L): Prepare the following standards by diluting suitable volumes of standard solution 1 (6.8) to 200.0 mL with Milli-Q water:

<u>mg P/L</u>	<u>mL of standard solution 1/200.0 mL</u>
0.50	5.0
1.00	10.0
1.50	15.0
2.00	20.0

## **7.0 Procedure**

- 7.1 Set up the manifold as shown in Figure 1. For the high concentration range, use the dilution manifold (Figure 1.).
- 7.2 Allow the colorimeter, and printer to warm up for 30 minutes. Obtain a stable baseline with all reagents, feeding Milli-Q water through the sample line.
- 7.3 Load the autosampler according to the CFDA Tray Protocol.
- 7.4 Analyze according to procedures in the LIMS-CFDA Methods Manual and General AutoAnalyzer Procedures.

## **8.0 Calculations**

The phosphorus concentration is obtained directly from the LIMS plotter.

## **9.0 Precision and Accuracy**

Precision and accuracy data are available in the Inorganic Chemistry Unit Quality Assurance Manual.

## **10.0 References**

- 10.1 Methods for Chemical Analysis of Water and Wastes, U.S. Environmental Protection Agency, EPA 600/4-79-020, p 365.1, (1979).
- 10.2 Methods for Determination of Inorganic Substances in Water and Fluvial Sediments, U.S. Geological Survey Techniques of WaterResources Inv., Book #5, Ch.A1, p 514, (1985).
- 10.3 Ortho Phosphate in Water and Seawater, Industrial Method No. 155-71W, Technicon Instruments Corporation, Tarrytown, NY (1973).

Figure 1. Manifold Set Up

