

# **Standard Operating Procedure for the Determination of Dissolved O<sub>2</sub> Correction Factors**

**LG303**

---

**Revision 01, April 2000**



## TABLE OF CONTENTS

<u>Section Number</u>	<u>Subject</u>	<u>Page</u>
<b>1.0.....</b>	<b>DEFINITIONS .....</b>	<b>1</b>
<b>2.0.....</b>	<b>DETERMINE AIR TEMPERATURE CORRECTION FACTOR .....</b>	<b>1</b>
<b>3.0.....</b>	<b>CORRECTED PRESSURE .....</b>	<b>1</b>
<b>4.0.....</b>	<b>THEORETICAL O<sub>2</sub> SATURATION.....</b>	<b>1</b>
<b>5.0.....</b>	<b>ACTUAL DISSOLVED O<sub>2</sub> CONCENTRATION .....</b>	<b>2</b>
<b>6.0.....</b>	<b>STATISTICS.....</b>	<b>2</b>



## Standard Operating Procedure for the Determination of Dissolved O<sub>2</sub> Correction Factors

### 1.0 DEFINITIONS

H <sub>2</sub> O Temperature (°C)	=	$T_w$
Air Temperature (°F)	=	$T_a$
Atmospheric Pressure (mbar)	=	$P_o$
Corrected Pressure (mbar)	=	$P_{c(mbar)}$
Corrected Pressure (mm Hg)	=	$P_{c(mm\ Hg)}$
Correction Factor (°F)	=	$C$
Theoretical O <sub>2</sub> Saturation (mg/L)	=	$S_T$
Adjusted Theoretical Saturation (mg/L)	=	$O_T$
Measured Dissolved Oxygen (mg/L)	=	$O_M$
Actual Dissolved Oxygen (mg/L)	=	$O_A$

### 2.0 DETERMINE AIR TEMPERATURE CORRECTION FACTOR

- 2.1 Use the following equation, which was calculated with data from Standard Methods, to determine the correction factor for the ambient air temperature (or laboratory temperature for the saturated standard):

$$C = -0.0434 T_a + 23.171$$

### 3.0 CORRECTED PRESSURE

- 3.1 Subtract the correction factor determined in the previous step, from the observed pressure to determine the corrected pressure given as:

$$P_c (mbar) = P_o - C$$

- 3.2 Use this result and convert the units of pressure with the following expression:

$$P_c (mm\ Hg) = 0.750062 P_c (mbar)$$

### 4.0 THEORETICAL O<sub>2</sub> SATURATION

- 4.0 Determine the theoretical O<sub>2</sub> saturation at  $T_w$  and 760 mm Hg with the following equation, which was calculated using temperature data from Standard Methods:

$$S_T = 0.0049 T_w - 0.3661 T_w + 14.534$$

- 4.1 Adjust this result using the corrected pressure in mm Hg,  $P_c (mm\ Hg)$ :

$$O_T = \frac{S_T \times P_C \text{ (mm Hg)}}{760}$$

## 5.0 ACTUAL DISSOLVED O<sub>2</sub> CONCENTRATION

5.1 Determine the actual dissolved O<sub>2</sub> concentration:

$$O_A = \frac{(60 O_M)}{(V_{bottle} - 0.8)}$$

## 6.0 STATISTICS

6.1 The relative percent deviation is given as:

$$RPD = \frac{|O_T - O_A|}{\left(\frac{O_T + O_A}{2}\right)} \times 100$$

6.2 While:

$$\% \text{ Saturation} = \left(\frac{O_A}{O_T}\right) \times 100$$