Standard Operating Procedure for the Determination of Dissolved O₂ Correction Factors

LG303

Revision 01, April 2000

TABLE OF CONTENTS

Section Number	Subject	<u>Page</u>
1.0	DEFINITIONS	1
2.0	DETERMINE AIR TEMPERATURE CORRECTION FACTOR	1
3.0	CORRECTED PRESSURE	1
4.0	THEORETICAL O ₂ SATURATION	1
5.0	ACTUAL DISSOLVED O ₂ CONCENTRATION	2
6.0	STATISTICS	2

Standard Operating Procedure for the Determination of Dissolved O₂ Correction Factors

1.0 DEFINITIONS

 $\begin{array}{llll} \text{H}_2\text{O Temperature (°C)} & = & T_w \\ \text{Air Temperature (°F)} & = & T_a \\ \text{Atmospheric Pressure (mbar)} & = & P_o \\ \text{Corrected Pressure (mbar)} & = & P_{c(mbar)} \\ \text{Corrected Pressure (mm Hg)} & = & P_{c(mm Hg)} \\ \text{Correction Factor (°F)} & = & C \\ \text{Theoretical O}_2 \text{ Saturation (mg/L)} & = & S_T \\ \text{Adjusted Theoretical Saturation (mg/L)} & = & O_T \\ \text{Measured Dissolved Oxygen (mg/L)} & = & O_M \\ \text{Actual Dissolved Oxygen (mg/L)} & = & O_4 \\ \end{array}$

2.0 DETERMINE AIR TEMPERATURE CORRECTION FACTOR

2.1 Use the following equation, which was calculated with data from Standard Methods, to determine the correction factor for the ambient air temperature (or laboratory temperature for the saturated standard):

$$C = -0.0434 T_a + 23.171$$

3.0 CORRECTED PRESSURE

3.1 Subtract the correction factor determined in the previous step, from the observed pressure to determine the corrected pressure, given as:

$$P_c$$
 (mbar) = $P_o - C$

3.2 Use this result and convert the units of pressure with the following expression:

$$P_c$$
 (mm Hg) = 0.750062 P_c (mbar)

4.0 THEORETICAL O₂ SATURATION

Determine the theoretical O_2 saturation at T_w and 760 mm Hg with the following equation, which was calculated using temperature data from Standard Methods:

$$S_T = 0.0049 \ T_w = -0.3661 \ T_w + 14.534$$

4.2 Adjust this result using the corrected pressure in mm Hg, P_c (mm Hg):

$$O_T = \frac{S_T \times P_c \ (mm \ Hg)}{760}$$

5.0 ACTUAL DISSOLVED O₂ CONCENTRATION

5.1 Determine the actual dissolved O_2 concentration:

$$O_A = \frac{(60 \ O_M)}{(V_{bottle} - 0.8)}$$

6.0 STATISTICS

6.1 The relative percent deviation is given as:

$$RPD = \frac{|O_T - O_A|}{\left(\frac{O_T + O_A}{2}\right)} \times 100$$

6.2 while:

%Saturation =
$$\left(\frac{O_A}{O_T}\right) \times 100$$