

## RELEASE NOTES – FORMS II Lite™ Version 5.1

This is a FORMS II Lite software upgrade for FORMS II Lite users who are currently using Version 5.1. Upon installation of this upgrade, immediately discontinue the use of all previous versions of FORMS II Lite.

### ***FORMS II Lite Technical Support***

DynCorp's FORMS II Lite Help Desk can be reached at 703-818-4200 (Mon. – Fri., 9AM to 5PM ET). Users may also request help on the Web at <http://dyncsdao1.dyncorp.com/itg/forms2lite/>. Finally, email can be sent to [f2lite@dyncorp.com](mailto:f2lite@dyncorp.com).

### ***The FORMS II Lite v5.1.47 Upgrade Includes:***

- The ability to view and print only the lab and Regional copies of the TR when the Show Archived TRs check box is selected.
- The ability to import and export the sample concentration and Case Complete status in the .F2L file.