

# Sensitivity enhancement for trace detection of carbonyls by submillimeter wave spectroscopy

Christopher D. Ball, Christopher N. Neese, Frank C. De Lucia, The Ohio State University

Ian C. MacGregor, Battelle

David M. Shelow, EPA/OAR/OAQPS

# Introduction

## Submillimeter wave (SMMW) spectroscopic gas sensor

- Developed by OSU and Battelle for DARPA
- Exceptional sensitivity, selectivity, and speed
- Adaptable to ambient air monitoring applications
  - Direct detection of formaldehyde, acrolein, NO<sub>x</sub>, SO<sub>x</sub>, etc.
  - Simultaneous detection of multiple air pollutants
- Reduce reliance on lab-based sample analysis



# Presentation outline


- Background
- Previous work
- Objective of current study
- Optical multipass design and assembly
- System characterization
- Diagnostic spectroscopy measurements
- Next steps

# Background

- EPA requires carbonyl detection in ambient air at sub-ppb concentrations
  - Formaldehyde, acrolein, acetaldehyde
- Current field monitoring methods need improvement
- SMMW spectroscopy offers near-real time detection of carbonyls
- Additional research needed to enhance SMMW sensitivity without relying on sorbent-based preconcentration



# Spectroscopic phenomenology

- High resolution rotational spectroscopy
  - Resolves individual spectral lines
  - Optimal sensitivity near Doppler limit ( $\sim 10$  mTorr)



# Previous work

- DARPA MACS program
  - Sensitivity: ~ppt (with preconcentration)
  - Selectivity: simultaneous detection of 30+ gases
  - False alarm rate:  $< 10^{-10}$
  - Speed: 10 min
  - Size: 1 cubic foot
- EPA/OAQPS feasibility study
  - Investigate sorbent-free MDL of formaldehyde, acrolein, acetaldehyde
  - Within factor of 200 for formaldehyde MDL



Neese, et al., *IEEE Sensors Journal* vol. 12, pp. 2565-2574, 2012

# Detection sensitivity

- Inherent molecular absorptivity
- Pressure/number density
  - Preconcentration
- SMMW power
- Signal averaging
- Optical path length




# Objective

- Develop optical multipass configuration to enhance sensitivity of SMMW spectrometer
  - Focus on carbonyl detection
    - Formaldehyde, acrolein, acetaldehyde
  - Demonstrate gains toward carbonyl minimum detection levels (MDLs)
  - Enable potential sorbent-free operation in a future fielded SMMW system








# Multipass design considerations


- Multipass vs. White cell
- Polarization
- Vacuum chamber
- Mirror selection
- Proof-of-concept design




# Multipass assembly



# Multipass assembly



Polarization-rotating retroreflector constructed in machine shop





Receiver


Wire grid beamsplitter

Transmitter

# Multipass characterization



Baseline changes with increasing mirrors  
Decreased baseline intensity  
Increased baseline fringe frequency  
Decreased baseline fringe intensity







**Average reflectivity = 99.3%**

# Test measurements: Carbonyl sulfide

- Simple linear rotor
- Conduct system diagnostic measurements
- Approx. 5 mTorr neat OCS in multipass cell
- Compare single pass vs. multipass
  - Estimate total gain from multipass
  - Estimate sensitivity, SNR
- Access to multiple isotopomers
  - OCS, OC<sup>34</sup>S (4%), O<sup>13</sup>CS (1%), <sup>18</sup>OCS (0.2%)

# OCS (neat, 5.5 mTorr)



~2,000 ppm



~10,000 ppm

# Test measurements: Acetonitrile

- Diagnostic tests at intermediate dilution
- 100 ppm  $\text{CH}_3\text{CN}$  in  $\text{N}_2$
- Demonstrate detection of  $^{13}\text{C}$  isotopomer (1% of normal species,  $\sim 1$  ppm)
  - 80 pass configuration only
- Results
  - Detected  $^{13}\text{C}$  isotopomer with decent SNR
  - Suggested that 1 ppm  $\text{H}_2\text{CO}$  is feasible

# CH<sub>3</sub>CN (100 ppm in N<sub>2</sub>)

Red = gas present  
Black = zero






# Test measurements: Formaldehyde

- Preliminary sensitivity tests
- Detect 1 ppm H<sub>2</sub>CO in N<sub>2</sub>
- Results
  - Detected strongest line at 20:1 SNR
  - Equivalent detection of approx. 50 ppb
  - No additional efforts taken to enhance signal or reduce noise

# Formaldehyde (1 ppm in N<sub>2</sub>)




- Baseline fringes x20 intensity of H<sub>2</sub>CO line
- Residual baseline fringes on subtracted signature
- SNR ~20
- Additional optimization needed

# Next steps

- Sensitivity optimization (ppb-level detection)
  - Compensation for baseline fringes
    - Algorithm approach
    - Leverage astrophysical procedures
  - Increased signal averaging
  - Decreased scan rate
  - Assess SMMW power and saturation effects
- Cold finger preconcentration
  - No consumables

# Technology development

- Development of ambient air monitoring prototypes
- Field testing at EPA monitoring stations
- Commercialization



# Conclusions

- SMMW spectroscopic sensor enables highly sensitive detection of carbonyls
- Optical multipass shows promise for approaching formaldehyde MDL (0.065 ppb)
- Additional optimization of sensitivity required
- Path to commercialization defined