



Analytical Method for PFAS in Environmental Media: CWA-1633

Marc A. Mills

Center for Environmental Solutions and Emergency Response
Office of Research and Development

Executive Meeting | Board of Scientific Counselors
September 29-30, 2021

The views expressed in this presentation are those of the author(s) and do not necessarily represent the views or policies of the U.S. Environmental Protection Agency.



PFAS Analytical Methods

- Drinking Water Methods
 - EPA 537.1
 - EPA 533
- Non-Drinking Water Methods
 - SW-846 Method 8327 – Direct Injection
 - Draft CWA-1633—Isotope Dilution
- PFAS Analysis in Marine Waters
- PFAS Analysis in Fish Tissue
- Total Organofluorine Analysis using Combustion Ion Chromatography (TOF)
- Total Oxidizable Precursors (TOP)
- Summary of EPA PFAS Methods as of April 2021



Types of Standard Methods

Three broad categories of EPA Standard Methods:

- [Safe Drinking Water Act Methods](#)
- [Clean Water Act Methods](#)
- [SW846 Methods](#)

These are generally targeted methods for solids and water.



Non-Drinking Water Sample Methods: CWA-1633/SW-846 Method—*Isotope Dilution*

Build in flexibility

- Columns
- Elution schemes


Single laboratory validated and released as draft CW-1633 in Sept 2021

- Collaborative effort among DoD, EPA Office of Water, EPA Office of Land and Emergency Management, and EPA ORD
- Multi-laboratory validation will start in 2021
- Method being developed in accordance SW-846 protocols for method development

Target Quantitation Limits: 1-10 nanogram/L

Matrices include:

- Wastewater (influent and effluent)
- Groundwater
- Surface water
- Landfill leachate
- Soil
- Sediment
- Biosolids
- Fish tissue





Non-Drinking Water Sample Methods: CWA-1633/SW-846 Method—*Isotope Dilution*

More complex method relative to direct injection; however, will:

- be more robust for complex matrices (e.g., wastewater influents, biosolids). Account for matrix effects (e.g., sorption) through isotopically-labelled standard recoveries;
- afford options to meet DoD requirements; and
- allow users to perform lower-level analyses based on screening results (e.g., 8327, TOF).

40 PFAS analytes - includes all analytes listed in 537.1, 533, and SW-846 8327

Non-drinking water samples

- Surface water, groundwater, wastewater
- Landfill leachates
- Solids (soils, sediments, biosolids, tissues)





Summary: EPA PFAS Methods, April 2021

EPA has validated Standard Methods complete or in development for PFAS in water

- Final SDWA Methods 533 and 537.1 for available for drinking water (29 PFAS)
- SW-846 Method 8327 validated for non-potable water (24 PFAS)
- Method in CWA-1633 completed single lab and undergoing multi lab validation for non-potable water/solids (40 PFAS). SW-846 determination to follow.

EPA has or is developing additional methods for partner use

- **Fish Tissue** – Isotope dilution method for 13 PFAS has been used in national surveys
- **Serum** – Isotope dilution method (targeted and non-targeted) used in biomonitoring
- **Ambient air and emissions** – Sampling and analysis methods undergoing development and testing
- **Total Organic Precursors (TOP)** – Identify total PFAS load which may degrade to most persistent PFAS
- **Total Organic Fluorine (TOF)** – Potential rapid screening tool to identify total PFAS presence/absence
- **Nontargeted analysis** – Continued development and application of HRMS methods for discovery of novel PFAS, suspect screening analysis, and identification of transformation and end products.



Contact

Marc A. Mills, PhD

US EPA Office of Research and Development

mills.marc@epa.gov

Acknowledgement:

Christopher Impellitteri, Andy Gillespie, and Susan Burden


US EPA ORD

The views expressed in this presentation are those of the individual author and do not necessarily reflect the views and policies of the US EPA.



EPA 823R18004 | February 2019 | www.epa.gov/pfas

EPA's Per- and Polyfluoroalkyl Substances (PFAS) Action Plan



U.S. Environmental Protection Agency