

Manufacturer: Neopac
Model: Twist 'n' Use CR
ASTM Type: VD

Description:

This is a unit dose non-reclosable rigid package. The package consists of a tube and a closure that is permanently attached to the tube neck.



Figure 1

The package is opened by pushing down on the closure and simultaneously turning it counterclockwise. (see arrows on tube closure on right in Figure 2)



Figure 2

A part of the tube neck is separated from the tube, but stays in a recessed area (see Figure 3). This permits the closure tip to be opened and have product dispensed through a hole in the closure (see drop at tip of tube on left in Figure 2 indicated with the arrow). The fill volume that this package can hold is from 3ml to 10ml. Once the unit dose package is opened it is non-reclosable. There are no manufacturer's instructions on the package.